**OHST**   
medical technology

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**IMPLANT****ZEN Knee System**

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party placing this product on the market accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for knee arthroplasty. Familiarity with the surgical technique recommended for this system and careful application thereof is essential in order to achieve the best possible result.

**1. Product description and implant materials**

Knee prostheses are implants that remain in the body permanently and either partially or fully replace the damaged knee. The knee joint in its entirety is replaced by three components: femoral component, tibial component and patellar component.

Irrespective of the shape and fixation manner of the individual joint components, they always take over the same functions (ISO 7207-1, ISO 21536):

- |                     |  |
|---------------------|--|
| Femoral component:  | Component of a knee replacement system that is attached to the femur to replace its gliding surfaces.<br>These implants can be comprised of a single component or of a set of components that have to be assembled by the user.                    |
| Tibial component:   | Component of a knee replacement system that is attached to the tibia to replace its gliding surfaces.<br>These implants can be comprised of a single component or of a set of components that have to be assembled by the user.                    |
| Patellar component: | Components of the total or partial knee replacement system that are used to replace the gliding surfaces of the patella.<br>These implants can be comprised of a single component or of a set of components that have to be assembled by the user. |

The product, the package contents and the materials used are specified on the product label. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. The explanations in the respective surgical instructions must be followed.

## 1.1. Overview of implants

### 1.1.1. Femoral components

#### Femur STD:

Size	Version	Reference number		
		Cemented		Cementless
		ISO 5832-4 CoCrMo	ISO 5832-4 CoCrMo TiNbn	ISO 5832-4 CoCrMo TPS, BONIT®
1	right	482-002	482-122	482-014
2	right	482-004	482-124	482-016
3	right	482-006	482-126	482-018
4	right	482-008	482-128	482-020
5	right	482-010	482-130	482-022
6	right	482-012	482-132	482-024
1	left	482-001	482-121	482-013
2	left	482-003	482-123	482-015
3	left	482-005	482-125	482-017
4	left	482-007	482-127	482-019
5	left	482-009	482-129	482-021
6	left	482-011	482-131	482-023

### 1.1.2. Tibial components

#### Tibia STD, mobile bearing:

Size	Version	Reference number		
		Cemented		Cementless
		ISO 5832-4 CoCrMo	ISO 5832-4 CoCrMo TiNbn	ISO 5832-4 CoCrMo TPS, BONIT®
1	right	482-098	482-134	482-026
2	right	482-100	482-136	482-028
3	right	482-102	482-138	482-030
4	right	482-104	482-140	482-032
5	right	482-106	482-142	482-034
6	right	482-108	482-144	482-036
1	left	482-097	482-133	482-025
2	left	482-099	482-135	482-027
3	left	482-101	482-137	482-029
4	left	482-103	482-139	482-031
5	left	482-105	482-141	482-033
6	left	482-107	482-143	482-035

**Tibia STD, fixed bearing:**

Size	Version	Reference number		
		Cemented		Cementless
		ISO 5832-4 CoCrMo	ISO 5832-4 CoCrMo TiNbN	ISO 5832-4 CoCrMo TPS, BONIT®
1	right	522-902	522-926	522-914
2	right	522-904	522-928	522-916
3	right	522-906	522-930	522-918
4	right	522-908	522-932	522-920
5	right	522-910	522-934	522-922
6	right	522-912	522-936	522-924
1	left	522-901	522-925	522-913
2	left	522-903	522-927	522-915
3	left	522-905	522-929	522-917
4	left	522-907	522-931	522-919
5	left	522-909	522-933	522-921
6	left	522-911	522-935	522-923

**Tibia Modular, with locking screw, mobile-bearing:**

Size	Reference number	
	Cemented	
	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V TiNbN
2	482-170	482-202
3	482-171	482-203
4	482-172	482-204
5	482-173	482-205
6	482-174	482-206

**Tibia Modular, with locking screw, fixed -bearing:**

Size	Reference number	
	Cemented	
	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V TiNbN
2	522-937	522-942
3	522-938	522-943
4	522-939	522-944
5	522-940	522-945
6	522-941	522-946

**Augment:**

Size	Height (mm)	Version	Reference number
			Cemented
			ISO 5832-3 Ti6Al4V
2	5	left lateral / right medial	482-182
3	5	left lateral / right medial	482-183
4	5	left lateral / right medial	482-184
5	5	left lateral / right medial	482-185
6	5	left lateral / right medial	482-186
2	10	left lateral / right medial	482-187
3	10	left lateral / right medial	482-188
4	10	left lateral / right medial	482-189
5	10	left lateral / right medial	482-190
6	10	left lateral / right medial	482-191
2	5	left medial / right lateral	482-192
3	5	left medial / right lateral	482-193
4	5	left medial / right lateral	482-194
5	5	left medial / right lateral	482-195
6	5	left medial / right lateral	482-196
2	10	left medial / right lateral	482-197
3	10	left medial / right lateral	482-198
4	10	left medial / right lateral	482-199
5	10	left medial / right lateral	482-200
6	10	left medial / right lateral	482-201

**Stem extension:**

Length	Version	Reference number
		Cemented
		ISO 5832-3 Ti6Al4V
40	straight	482-175
75	straight	482-176
110	straight	482-177
110	with offset	482-178

**Centraliser:**

Size	Reference number
	Cemented
	ISO 5834-2 UHMWPE
Ø10	482-179
Ø12	482-180
Ø14	482-181

**Insert STD, mobile-bearing:**

<b>Size</b>	<b>Height (mm)</b>	<b>Reference number ISO 5834-2 UHMWPE ISO 5832-11.4441</b>	<b>Size</b>	<b>Height (mm)</b>	<b>Reference number ISO 5834-2 UHMWPE ISO 5832-11.4441</b>
1	10	482-054	4	10	482-078
1	12	482-055	4	12	482-079
1	14	482-056	4	14	482-080
2	10	482-062	5	10	482-086
2	12	482-063	5	12	482-087
2	14	482-064	5	14	482-088
3	10	482-070	6	10	482-094
3	12	482-071	6	12	482-095
3	14	482-072	6	14	482-096

**Insert UC, mobile-bearing:**

<b>Size</b>	<b>Height (mm)</b>	<b>Reference number ISO 5834-2 UHMWPE ISO 5832-11.4441</b>	<b>Size</b>	<b>Height (mm)</b>	<b>Reference number ISO 5834-2 UHMWPE ISO 5832-11.4441</b>
1	10	482-049	4	10	482-073
1	12	482-050	4	12	482-074
1	14	482-051	4	14	482-075
1	17	482-052	4	17	482-076
1	20	482-053	4	20	482-077
2	10	482-057	5	10	482-081
2	12	482-058	5	12	482-082
2	14	482-059	5	14	482-083
2	17	482-060	5	17	482-084
2	20	482-061	5	20	482-085
3	10	482-065	6	10	482-089
3	12	482-066	6	12	482-090
3	14	482-067	6	14	482-091
3	17	482-068	6	17	482-092
3	20	482-069	6	20	482-093

**Insert STD, fixed bearing:**

Size	Height (mm)	Reference number		Size	Height (mm)	Reference number	
		ISO 5834-2 UHMWPE	Vitamin-E XL-PE			ISO 5834-2 UHMWPE	Vitamin-E XL-PE
1	10	522-947	522-801	4	10	522-962	522-816
1	12	522-948	522-802	4	12	522-963	522-817
1	14	522-949	522-803	4	14	522-964	522-818
1	17	522-950	522-804	4	17	522-965	522-819
1	20	522-951	522-805	4	20	522-966	522-820
2	10	522-952	522-806	5	10	522-967	522-821
2	12	522-953	522-807	5	12	522-968	522-822
2	14	522-954	522-808	5	14	522-969	522-823
2	17	522-955	522-809	5	17	522-970	522-824
2	20	522-956	522-810	5	20	522-971	522-825
3	10	522-957	522-811	6	10	522-972	522-826
3	12	522-958	522-812	6	12	522-973	522-827
3	14	522-959	522-813	6	14	522-974	522-828
3	17	522-960	522-814	6	17	522-975	522-829
3	20	522-961	522-815	6	20	522-976	522-830

**Insert UC, fixed bearing:**

Size	Height (mm)	Reference number	Size	Height (mm)	Reference number
		Vitamin-E XL-PE			Vitamin-E XL-PE
1	10	522-831	4	10	522-846
1	12	522-832	4	12	522-847
1	14	522-833	4	14	522-848
1	17	522-834	4	17	522-849
1	20	522-835	4	20	522-850
2	10	522-836	5	10	522-851
2	12	522-837	5	12	522-852
2	14	522-838	5	14	522-853
2	17	522-839	5	17	522-854
2	20	522-840	5	20	522-855
3	10	522-841	6	10	522-856
3	12	522-842	6	12	522-857
3	14	522-843	6	14	522-858
3	17	522-844	6	17	522-859
3	20	522-845	6	20	522-860

### 1.1.3. Patellar components

#### Patella STD, mobile-bearing:

Size	Reference number		
	Cemented		Cementless
	ISO 5832-12 CoCrMo ISO 5834-2 UHMWPE ISO 5832-1 1.4441	ISO 5832-12 CoCrMo TiNbN ISO 5834-2 UHMWPE ISO 5832-1 1.4441	ISO 5832-12 CoCrMo TPS, BONIT® ISO 5834-2 UHMWPE ISO 5832-1 1.4441
1	482-109	482-145	482-115
2	482-110	482-146	482-116
3	482-111	482-147	482-117
4	482-112	482-148	482-118
5	482-113	482-149	482-119
6	482-114	482-150	482-120

#### Patella FIX, fixed-bearing:

Size	Reference number
	Cemented
	ISO 5834-2 UHMWPE ISO 5832-1 1.4441
1	482-163
2	482-164
3	482-165
4	482-166
5	482-167
6	482-168

## 1.2. Overview of instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
ZEN size 1 & 6 instrument set	482-541
ZEN Fixed Bearing size 1 & 6 instrument set	522-120
ZEN Femur instrument set	482-542
ZEN General instrument set	482-543
ZEN Fixed Bearing General instrument set	522-121
ZEN Trial Components instrument set	482-544
ZEN Fixed Bearing Trial Components instrument set	522-122
ZEN Patella instrument set	482-545
ZEN Tibia instrument set	482-546
ZEN Tibia Modular instrument set	482-741
ZEN Anterior Cut First instrument set	482-742

<b>Designation</b>	<b>Reference number</b>
ZEN Optional instrument set	482-775

### 1.3. Miscellaneous accessories

<b>Designation</b>	<b>Reference number</b>
ZEN Femur + Insert X-ray templates	482-507
ZEN Femur X-ray templates	482-559
ZEN Tibia X-ray templates	482-508
ZEN Tibia Fixed Bearing with Insert X-ray templates	522-501
ZEN Tibia Modular X-ray templates	482-738
ZEN Tibia Modular Fixed Bearing with Insert and Stem Extension Straight X-ray templates	522-502
ZEN Tibia Modular Offset Stem X-ray templates	482-739
ZEN Tibia Modular Fixed Bearing with Insert and Stem Extension with offset X-ray templates	522-503
ZEN Knee System Mobile Bearing surgical technique	50000434
ZEN Tibia Modular Mobile Bearing surgical technique	50000435
ZEN Knee System Fixed Bearing surgical technique	50000656
ZEN Tibia Modular Fixed Bearing surgical technique	50000657
Implant passport	50000572

## 2. Handling

### 2.1. General information

This implant is part of a system and must only be used with the corresponding original system components. For implantation, only the above-mentioned instruments of the system are to be used. Before the instruments are used, the corresponding instructions for use must be followed.

**Caution:** Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When the implant is unpacked, a check must be performed to ensure that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the appropriate hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Machining or bending an implant cannot only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

**Caution:** The implants are intended for single use only! The individual loads on the functional surfaces of a patient affect the articulation surfaces to such an extent that re-use is excluded. The stress marks on the functional surfaces cannot be reliably seen with

visual methods alone. Therefore, after explanation it has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding side is marked on the implants with "L" for the left knee and "R" for the right knee. The orientation of the implants must absolutely correspond to the side of the body of the knee to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right knee.

## **2.2. Authorised component combinations**

Together with its instruments, the ZEN Knee System is designed as a complete system and does not allow for the substitution of components from other systems of manufacturers. Only the above-mentioned system components may be used for implantation.

### **Compatibility instructions ZEN Mobile Bearing:**

Each of the femoral components must be combined with inserts of the same size.

The inserts can be combined with tibial components of the same size as well as with all larger tibial components. This applies to both the STD inserts and the UC inserts.

### **Compatibility instructions ZEN Fixed Bearing:**

Each of the tibial components must be combined with inserts of the same size.

The inserts can be combined with femoral components of the same size as well as with the next larger and the next smaller size. This applies to both the STD inserts and the UC inserts.

### **General compatibility instructions:**

The femoral components, stem extensions and augments can be used for both the Fixed Bearing variant and for the Mobile Bearing variant. Each of the patellar components are to be used with the femoral components of the same size.

We only guarantee compatibility of our products in combination with our own CE-marked products and with the products we have approved for combined use, for which the competent authority has issued its approval. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST. Due to reasons relating to the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

## **2.3. Application instructions**

Before the instruments are used, the enclosed product information and instructions for use must be observed.

The surfaces with a porous coating (TPS, BONIT®) and the roughened implant surfaces must not come in contact with clothing or other fibre-shedding materials.

**Caution:** Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free latex gloves.

For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

Prior to introducing the cement (for cemented fixation) or before inserting the implant (for cementless fixation) the implant bed has to be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed. In the case of cemented anchoring, the implants must be inserted centrally and straight into the cement bed. Following the cementing step, all excess or loose cement particles must be removed from the wound area.

**Caution:** When using high-frequency surgical instruments (e. g. cauter knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e. g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

## **2.4. Surgical technique**

### **2.4.1. ZEN Knee**

#### Approach

- Parapatellar incision with medial arthrotomy or another approach, e.g. for a revision procedure or marked valgus malposition

#### Flexion gap - tibia processing step

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Mounting and positioning of the tibial goniometer
- Adjustment of rotation in the frontal plane
- Adjustment of the slope
- Adjustment of the cutting height
- Attachment of the tibial cutting block
- Performance of the tibial cut
- Verification of the tibial cut

#### Flexion gap - femur processing step

- Determination of the femur size
- Intramedullary opening
- Align and fix femoral cutting block in a.-p. position
- Anterior and posterior femoral resection

**Caution:** This step is a temporary step which has to be finalised in a later work step.

- Review of the flexion gap

#### Extension gap - distal femur processing step

- Distal positioning of the cutting block
- Review of the angle and distal cut
- Review of the extension gap

### Finalisation of the femur

- Attachment of the femoral cutting block finisher
- Execution of the cuts

**Caution:** It is important to carry out all cuts precisely, as the previous A/P incision was only temporary.

- Verification of the cuts
- Preparation of the trochlear groove

### Tibia processing step

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Selection of the plateau
- Tibia finalisation

### Patella processing step

For reduced bone thickness, the thinner full-poly patella should be employed, as it guarantees greater preservation of the bone. Prosthetic patellar replacement is contraindicated where the patella is too thin or in the case of a flat patella.

- Removal of the osteophytes
- Align and fix patellar resection guide perpendicular to the patellar tendon
- Measurement of the patellar thickness (at least 13 mm of bone must remain after deducting the thickness of the patellar component)

Components	Height of patella STD	Height of patella FIX
Size 1	11.2 mm	9.7 mm
Size 2	11.2 mm	9.7 mm
Size 3	11.2 mm	9.7 mm
Size 4	11.2 mm	9.7 mm
Size 5	11.7 mm	10.2 mm
Size 6	11.7 mm	10.2 mm

- Adjustment of the height of the drilling gauge and insertion of the pins
- Removal of the patellar resection guide
- Positioning of the patellar cutting block
- Performance of the cut
- Alignment of the base plate of the trial patella
- Clamping of the patellar drilling guide into the patella forceps
- Positioning of the patella base plate and compression of the patellar forceps
- Drilling of the pin holes
- Removal of the patellar forceps and drilling guide
- Insertion of the trial patella insert

### Attachment of the trial components

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Components are inserted in the following order: Trial tibia, trial insert, trial femur and, if applicable, the trial patella

**Caution:** The implant components must not be used together in any random combination. Compatibility instructions are provided in section 2.2.

**Caution:** Should the height of the insert have to be changed or the trial components removed after successfully checking the overall mobility of the joint, the components are removed in the reverse order of the insertion. In other words: first remove the trial femur, followed by the trial insert.

#### Implantation

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Components are inserted in the following order: Tibial components, insert, femoral component and, if applicable, the patellar component

**Caution:** The UC (ultra-congruent) inserts are to be used only after the resection of both cruciate ligaments. They do not provide any room for a posterior cruciate ligament that is still intact.

### **2.4.2. ZEN Tibia Modular**

#### Flexion gap - tibia processing step

Extramedullary alignment:

- Mounting and positioning of the tibial goniometer
- Adjustment of rotation in the frontal plane
- Adjustment of the slope
- Adjustment of the cutting height
- Attachment of the Tibia Modular cutting block
- Performance of the tibial cut
- Verification of the tibial cut

Intramedullary alignment:

- Positioning of the intramedullary tibial alignment gauge
- Medial-lateral centred opening of the tibia, aligned 1/3 anterior and 2/3 posterior
- Use of the medullary drill in ascending order, starting with the smallest diameter
- Mounting of the intramedullary tibial alignment gauge on the in-situ medullary drill
- Adjustment of rotation in the frontal plane
- Setting and fixing of the slope to 3°

**Caution:** Should a different slope setting than 3° be chosen, the further processing steps of the tibia are carried out after performing the cut according to the extramedullary alignment. Only straight stem extensions may be used in this case.

- Adjustment of the cutting height

- Attachment of the Tibia Modular cutting block
- Performance of the tibial cut
- Verification of the tibial cut

### Tibia processing step

Extramedullary alignment:

- Selection of the plateau
- Positioning of the Tibia Modular cutting block (optional, only when using augments)
- Performance of the cuts (optional, only when using augments)
- Tibia finalisation

Intramedullary alignment:

**Caution:** The following work steps for intramedullary alignment can only be carried out if the resection was previously carried out with a slope setting of 3° using the intramedullary tibial alignment gauge.

If a different slope setting was used, the further processing of the tibial is done as per the extramedullary alignment. Only straight stem extensions may be used in this case.

- Selection of the plateau
- Positioning of the tibial size gauge without offset or with offset
- Transfer of the offset angle to the trial cone (optional, only when using the stem extension with offset)
- Positioning of the tibia cutting block for augments (optional, only when using augments)
- Performance of the cuts (optional, only when using augments)
- Tibia finalisation

### Attachment of the trial components

- Components are inserted in the following order: Trial tibia with trial stem extension and optional trial augments and trial centralisers, trial insert, trial femur and, if applicable, the trial patella

**Caution:** The implant components must not be used together in any random combination. Compatibility instructions are provided in section 2.2.

**Caution:** Should the height of the insert have to be changed or the trial components removed after successfully checking the overall mobility of the joint, the components are removed in the reverse order of the insertion. In other words: first remove the trial femur, followed by the trial insert.

### Implantation

- Components are inserted in the following order: Tibial component with stem extension and optional augments and centralisers, insert, femoral component and, if applicable, the patellar component
- The stem extension is placed onto the tibia and fixed in place with the provided screw using the torque screwdriver
- The alignment of the stem extension with offset on the tibia is done using the alignment instrument

- Optional augments are attached with the torque screwdriver

**Caution:** The UC (ultra-congruent) inserts are to be used only after the resection of both cruciate ligaments. They do not provide any room for a posterior cruciate ligament that is still intact.

It is important to ensure that the surfaces are not damaged. Appropriate impactors are available for impacting the implants. Position the impactor at the centre and align it parallel to the mechanical axis. The patella is pushed into the bone using the patellar forceps and the pressing aid.

In order to prevent excessive heating of the bone, all resections must be performed while cooling with physiological saline solution. The thickness of the saw blades must be adjusted to the instruments used!

The smaller height of the fixed-bearing patellar component may provide greater bone preservation.

For further information see the surgical technique of the knee system (see 1.3. Miscellaneous accessories).

### **3. Packaging and sterility**

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use. The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

**Caution:** The implants must not be re-sterilised. Reprocessing of components which have not been implanted and for which the packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

### **4. Preoperative planning and postoperative care**

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the probable combination of components, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

Both conventional and digital X-ray templates (see 1.3. Miscellaneous accessories) are available for planning the surgery. To check the correct fit, trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognized procedures must be used for postoperative care.

## **5. Indication**

Total prosthesis is indicated in patients who have lost their mobility for any of the following reasons:

- Degenerative, post-traumatic or rheumatoid arthritis
- Avascular necrosis of the femoral condyle
- Post-traumatic loss of the joint configuration, in particular in the case of patellofemoral erosion, dysfunction or prior patellectomy
- Moderate varus, valgus or flexion malpositions

This system may also be indicated for the revision of earlier surgical attempts which have been unsuccessful.

Naturally, a total knee replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would, however, in many cases, it can restore mobility whilst simultaneously relieving pain. A total knee prosthesis should be the course of action only where all other treatment approaches, which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

## **6. Contraindication**

- Acute or chronic, local or systemic infection
- Severe muscle, nerve or vascular diseases endangering the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Hypersensitivity to the materials used
- Loss of the ligamenture
- Overweight
- Alcoholism or other addictions
- High-level/intense physical activity (e.g. competitive sports, heavy physical labour)
- Any underlying condition which could compromise the function of the implant

## **7. Risks and conditions that may affect the success of the operation**

Potential risks associated with the operation are:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Overweight
- Alcoholism or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Growth in children and adolescents
- Expected extreme stresses e. g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities which complicate anchoring of the implant
- Weakening of the supporting structures due to a tumour
- High-dose administration of cortisone or cytotoxic drugs
- History or risk of infectious diseases with potential joint manifestation

- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

**Caution:** If any of the conditions stated above apply to the patient, the latter must be informed of the likely impact of these factors on the success of the operation and the lifetime of the implant.

## **8. Undesirable effects**

The following adverse reactions may be associated with total knee arthroplasty, amongst others:

- Tibial or femoral fracture
- Sub-dislocation or dislocation of the patella
- Instability, movement or loosening of components
- Separation of components
- Insufficient freedom of movement
- Metal hypersensitivity reactions
- Patellofemoral complications
- Neurovascular complications
- Paralysis of nervus peroneus
- Fat embolism
- Arterial insufficiency / injury
- Arthrofibrosis / adhesions
- Implant breakage
- Infection
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Inflammation
- Oedem/fluids
- Metallosis
- Elevated metal ions in blood
- Pseudotumours
- Cardiac arrhythmia, increased pulmonary vascular resistance, cardiac arrest (caused by BCIS = bone cement implantation syndrome)

## **9. Patient information, documentation**

The serial numbers of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, although some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks. Patients receiving a knee replacement must be informed that the life of the implant will depend on the patient's weight and how active the patient is in everyday life. The patient must be informed of which activities he can undertake in order to reduce the effects of these aggravating circumstances.

All information provided to the patient must be documented in writing by the surgeon performing the operation. Following surgery, the patient must be given an implant passport containing all necessary

information about the implant. Undesirable effects that can harm the patient can occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant. Before use, the information provided by the equipment manufacturer should be studied. As part of an individual risk assessment, the suitability of comparable implants should be checked in the respective MRI machine. The patient is to be informed of the risks.

#### **10. Explanation of the label symbols**

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 90).

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Puede ver los símbolos utilizados por OHST Medizintechnik AG en el anexo (S. 90).