**OHST**   
medical technology

OHST Medizintechnik AG  
Grünauer Fenn 3  
14712 Rathenow  
Deutschland / Germany / Allemagne /  
Germania / Alemania  
Tel.: +49 (0)3385 - 54 20-0  
Fax: +49 (0)3385 - 54 20-99  
[www.ohst.de](http://www.ohst.de)

**CE** 0482

---

EFK Kniesystem	
Deutsch	3
EFK Knee System	
English	16
EFK Système de genou	
Français	28
EFK Sistema per ginocchio	
Italiano	40
Sistema de rodilla EFK	
Español	53

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## IMPLANT

### EFK Knee System

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party placing this product on the market accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for knee arthroplasty. Familiarity with the surgical technique recommended for this system and careful application thereof is essential in order to achieve the best possible result.

#### 1. Product description and implant materials

Knee prostheses are implants that remain in the body permanently and either partially or fully replace the damaged knee. The knee joint in its entirety is replaced by three components: femoral component, tibial component and patellar component.

Irrespective of the shape and fixation manner of the individual joint components, they always take over the same functions (ISO 7207-1, ISO 21536):

**Femoral component:** Component of a knee replacement system that is attached to the femur to replace its gliding surfaces.

These implants can be comprised of a single component or of a set of components that have to be assembled by the user.

**Tibial component:** Component of a knee replacement system that is attached to the tibia to replace its gliding surfaces.

These implants can be comprised of a single component or of a set of components that have to be assembled by the user.

**Patellar component:** Components of the total or partial knee replacement system that are used to replace the gliding surfaces of the patella.

These implants can be comprised of a single component or of a set of components that have to be assembled by the user.

The product, the package contents and the materials used are specified on the product label. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. The explanations in the respective surgical instructions must be followed.

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### 1.2. Overview of instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
Instrument Set EFK Tibia	10006776
Instrument Set EFK Femur Sz. 4-10	10006778
Instrument Set EFK Trials Sz. 4-10	10006779
Instrument Set EFK Additional Sz. 2+12	10006780
Instrument Set EFK Trial Inserts UC	10006781
Instrument Set Patella	10006821
Instrument Set EFK Optional	10007371
Instrument Set EFK Tibia Modular	10007372

### 1.3. Miscellaneous accessories

Designation	Reference number
X-ray templates EFK Femur	100000510296
X-ray templates EFK Tibia+Insert	100000510297
X-ray templates EFK Tibia Modular	100000510319
X-ray templates EFK S Femur	100000511064
Surgical technique EFK Knee System	50000673
Surgical technique EFK Tibia Modular	50000674
Implant passport	50000572

## 2. Handling

### 2.1. General information

This implant is part of a system and must only be used with the corresponding original system components. For implantation, only the above-mentioned instruments of the system are to be used. Before the instruments are used, the corresponding instructions for use (50000354) must be followed.

**Caution:** Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When the implant is unpacked, a check must be performed to ensure that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the appropriate hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Machining or bending an implant cannot only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

**Caution:** The implants are intended for single use only! The individual loads on the functional surfaces of a patient affect the articulation surfaces to such an extent that re-use is excluded. The stress marks on the functional surfaces cannot be reliably seen with visual methods alone. Therefore, after explantation it has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding side is marked on the implants with "L" for the left knee and "R" for the right knee. The orientation of the implants must absolutely correspond to the side of the body of the knee to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right knee.

### 2.2. Authorised component combinations

The femoral components are to be combined with inserts of the same size. The combination of inserts one size above or under nominal size is also allowed. Inserts of one size are only to be combined with

tibial components of the same size. Given standard for combination of standard inserts shall also apply for inserts UC. The patella component can be combined with all sizes of the femoral component.

We only guarantee compatibility of our products in combination with our own CE-marked products and with the products we have approved for combined use, for which the competent authority has issued its approval. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST. Due to reasons relating the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

### 2.3. Application instructions

Together with its instruments, the EFK Knee System is designed as a complete system and does not allow for the substitution of components from other systems of manufacturers. Only the above-mentioned system components may be used for implantation.

Before the instruments are used, the enclosed product information and instructions for use must be observed.

The surfaces with a porous coating (TPS, BONIT®) and the roughened implant surfaces must not come in contact with clothing or other fibre-shedding materials.

**Caution:** Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free latex gloves.

For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

Prior to introducing the cement (for cemented fixation) or before inserting the implant (for cementless fixation) the implant bed has to be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed. In the case of cemented anchoring, the implants must be inserted centrally and straight into the cement bed. Following the cementing step, all excess or loose cement particles must be removed from the wound area.

**Caution:** Inserts UC shall only be used when there is no intact posterior cruciate ligament, because there is no room for it.

**Caution:** When using high-frequency surgical instruments (e. g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e. g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

### 2.4. Surgical technique

Sequence of resections:

1. Distal femoral resection
2. A/P-femoral resections
3. Facet resections (including ablation of dorsal osteophytes)
4. Tibial resection (can be performed alternatively as resection 1)
5. Patellar resection (optionally)

It is important that the flexion gap and extension gap are approximately the same size.

Approach:

For information on accessing the knee joint, we refer to the relevant surgical procedures.

Preparation of the femur:

1. Distal femoral opening

2. Positioning of the distal femoral resection gauge

- At its distal end, the femoral implant is 9 mm thick.

- Bone pins must always be pre-drilled using a 3.2 mm drill.

- Before covering the patient, it is advisable to apply an electrode to the patient at the height of the femoral head centre using a C-arm X-ray unit, as the position depending on the soft tissue situation and the covering can otherwise not be precisely determined.

3. Distal femoral resection

4. Determination of the femoral size and rotation

- If you choose the next larger femoral size, this will increase the posterior offset and reduce the flexion gap. If you select the smaller size, this will increase the flexion gap and reduce the posterior offset. A reduced posterior offset can diminish the flexion ability.

5. A/P resection, facet resection

- If you perform the upper facet cut before the ventral cut, this will reduce the risk of the saw blade "wandering off", in particular when working on sclerotic bones.

6. Preparing the trochlear groove

- A lateral position of the femoral component is beneficial for the patella tracking.

- In case the position of the trial component is not stable enough to work on the trochlea, leave the stop drill in one of the drill holes.

Preparation of the tibia:

1. Installation of the extramedullary alignment gauge

2. Adjustment of the slope and varus/valgus alignment

3. Adjustment of the rotation

4. Adjustment of the resection height

Alternatively, intramedullary alignment: opening of the tibia, adjustment of rotation, resection height and slope

5. Axis check

6. Resection of the tibia

**Caution:** When displacing the revision cutting gauge to the distal or proximal direction (for thicker or thinner bone resections), this is done either to the right (+/- 2 mm) or parallel (+/- 4 mm). In the case of a left-hand shift (approx. pin width), make sure that the saw blade is applied to the outside of the right pin.

7. Determination of the tibial component size

8. Trial run

- After adjusting the rotation, use a cautery knife or a sterile pen to mark the ventral tibia just below the two anterior grooves of the tibial size gauge. This makes it easier to retrieve the correct rotation position for the tibial size gauge after removing the trial femur and the trial insert.

9. Preparation of the medullary cavity of the tibia for the stem

**Caution:** In combination with original implants, a cleaning of the trial insert is urgently necessary in order to avoid damages and scratches on the implant surfaces by third bodies.

**Patella preparation:**

1. Resection

- During resecting, make sure that the saw blade does not migrate, e. g. due to sclerotic bone segments.

2. Patella size determination

3. Milling

4. Drill anchoring holes

5. Manipulating

**Implantation:**

1. Insert the tibial component and impact by means of the impactor

2. Introduce the insert and click it ventrally

3. Fix the insert by tightening the screw using the screwdriver

4. Insert the femoral component and impact by means of the impactor

**Caution:** The sequence of implantation is generally the responsibility of the surgeon.

**Caution:** The tibia fixation screw must be tightened with the supplied screwdriver as firmly as possible without any additional tools.

**Caution:** Care must be taken not to damage the surfaces. Use the femur/tibia impactor for punching in, put it on centrally and align it parallel to the mechanical axis.

**Caution:** It should be noted that an insert removed from the tibial component or jammed when snapped into place is to be visually inspected for damage before re-use. In the event of damage, the parts may no longer be used.

In order to prevent excessive heating of the bone, all resections must be performed while cooling with physiological saline solution. The thickness of the saw blades must be adjusted to the instruments used!

For further information see the surgical technique of the knee system (see 1.3. Miscellaneous accessories).

### 3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

**Caution:** The implants must not be re-sterilised. Reprocessing of components which have not been implanted and for which the packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

#### 4. Preoperative planning and postoperative care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the probable combination of components, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

Both conventional and digital X-ray templates (see 1.3. Miscellaneous accessories) are available for planning the surgery. To check the correct fit, trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognized procedures must be used for postoperative care.

#### 5. Indication

Total prosthesis is indicated in patients who have lost their mobility for any of the following reasons:

- Degenerative, post-traumatic or rheumatoid arthritis
- Avascular necrosis of the femoral condyle
- Post-traumatic loss of the joint configuration, in particular in the case of patellofemoral erosion, dysfunction or prior patellectomy
- Moderate varus, valgus or flexion malpositions

This system may also be indicated for the revision of earlier surgical attempts which have been unsuccessful.

Naturally, a total knee replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would, however, in many cases, it can restore mobility whilst simultaneously relieving pain. A total knee prosthesis should be the course of action only where all other treatment approaches, which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

#### 6. Contraindication

- Acute or chronic, local or systemic infection
- Severe muscle, nerve or vascular diseases endangering the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Hypersensitivity to the materials used
- Loss of the ligamenture
- Overweight
- Alcoholism or other addictions
- High-level/intense physical activity (e.g. competitive sports, heavy physical labour)
- Any underlying condition which could compromise the function of the implant

## 7. Risks and conditions that may affect the success of the operation

Potential risks associated with the operation are:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Overweight
- Alcoholism or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Growth in children and adolescents
- Expected extreme stresses e. g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities which complicate anchoring of the implant
- Weakening of the supporting structures due to a tumour
- High-dose administration of cortisone or cytotoxic drugs
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Caution: If any of the conditions stated above apply to the patient, the latter must be informed of the likely impact of these factors on the success of the operation and the lifetime of the implant.

## 8. Undesirable effects

The following adverse reactions may be associated with total knee arthroplasty, amongst others:

- Tibial or femoral fracture
- Sub-dislocation or dislocation of the patella
- Instability, movement or loosening of components
- Separation of components
- Insufficient freedom of movement
- Metal hypersensitivity reactions
- Patellofemoral complications
- Neurovascular complications
- Paralysis of nervus peroneus
- Fat embolism
- Arterial insufficiency / injury
- Arthrofibrosis / adhesions
- Implant breakage
- Infection
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Inflammation
- Oedem/fluids
- Metallosis

- Elevated metal ions in blood
- Pseudotumours
- Cardiac arrhythmia, increased pulmonary vascular resistance, cardiac arrest (caused by BCIS = bone cement implantation syndrome)

#### 9. Patient information, documentation

The serial numbers of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, although some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks. Patients receiving a knee replacement must be informed that the life of the implant will depend on the patient's weight and how active the patient is in everyday life. The patient must be informed of which activities he can undertake in order to reduce the effects of these aggravating circumstances.

All information provided to the patient must be documented in writing by the surgeon performing the operation. Following surgery, the patient must be given an implant passport containing all necessary information about the implant. Undesirable effects that can harm the patient can occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant. Before use, the information provided by the equipment manufacturer should be studied. As part of an individual risk assessment, the suitability of comparable implants should be checked in the respective MRI machine. The patient is to be informed of the risks.

#### 10. Explanation of the label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 66).

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Puede ver los símbolos utilizados por OHST Medizintechnik AG en el anexo (S. 66).