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Grünauer Fenn 3

14712 Rathenow

Deutschland / Germany / Allemagne /

Germania / Alemania / Duitsland /

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IMPLANT

Burch-Schneider reinforcement cage

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance. These implants may be used only by physicians with appropriate experience and practice in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

1. Product description and implant materials

The primary function of the implants is to ensure stable bridging of bony defects of the acetabulum and to restore the centre of rotation. The aim is to restore and preserve the hip joint as well as its mobility and load-bearing capacity. The reinforcement cage is designed for the treatment of acetabulums that are already so badly damaged by bony defects that treatment with primary implants such as press-fit cups or other devices is no longer indicated. The reinforcement cage is fixed in place without cement and with bone screws. The reinforcement cage is made of titanium sheet metal (ISO 5832-2 Grade 1) and comes in four sizes, each with a left and right variant. The bone screws are made of Ti6Al4V alloy (ISO 5832-3). The outer surfaces are rough blasted, the inner surfaces are smooth blasted. This means that, on one hand, a high degree of secondary stability is achieved by growth of the bone on the outside of the cup, while the blasted surface on the inside of the implant provides a good adhesive surface for the cement fixation of a polyethylene cup.

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Reference number
Burch-Schneider reinforcement cage Sz. 44R (Outer-Ø50)	ISO 5832-2 unalloyed Titanium	367-1418
Burch-Schneider reinforcement cage Sz. 44L (Outer-Ø50)	ISO 5832-2 unalloyed Titanium	367-1419
Burch-Schneider reinforcement cage Sz. 50R (Outer-Ø56)	ISO 5832-2 unalloyed Titanium	367-1420
Burch-Schneider reinforcement cage Sz. 50L (Outer-Ø56)	ISO 5832-2 unalloyed Titanium	367-1421
Burch-Schneider reinforcement cage Sz. 56R (Outer-Ø62)	ISO 5832-2 unalloyed Titanium	367-1422

Name	Material	Reference number
Burch-Schneider reinforcement cage Sz. 56L (Outer-Ø62)	ISO 5832-2 unalloyed Titanium	367-1423
Burch-Schneider reinforcement cage Sz. 62R (Outer-Ø68)	ISO 5832-2 unalloyed Titanium	367-1424
Burch-Schneider reinforcement cage Sz. 62L (Outer-Ø68)	ISO 5832-2 unalloyed Titanium	367-1425
Flat-head screw Ø 6.5 x 15 mm	ISO 5832-3 Ti6Al4V	000-290-15
Flat-head screw Ø 6.5 x 20 mm	ISO 5832-3 Ti6Al4V	000-290-20
Flat-head screw Ø 6.5 x 25 mm	ISO 5832-3 Ti6Al4V	000-290-25
Flat-head screw Ø 6.5 x 30 mm	ISO 5832-3 Ti6Al4V	000-290-30
Flat-head screw Ø 6.5 x 35 mm	ISO 5832-3 Ti6Al4V	000-290-35
Flat-head screw Ø 6.5 x 40 mm	ISO 5832-3 Ti6Al4V	000-290-40
Flat-head screw Ø 6.5 x 45 mm	ISO 5832-3 Ti6Al4V	000-290-45
Flat-head screw Ø 6.5 x 50 mm	ISO 5832-3 Ti6Al4V	000-290-50
Flat-head screw Ø 6.5 x 55 mm	ISO 5832-3 Ti6Al4V	000-290-55
Flat-head screw Ø 6.5 x 60 mm	ISO 5832-3 Ti6Al4V	000-290-60

1.2 Instrument overview

The instruments of the OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Reference number
Instruments Burch-Schneider Reinforcement Cage	367-1482
Instrument Set Trial Burch-Schneider Reinforcement Cage	367-1483
Instrument Set Acetabulum Reamer	367-147

1.3 Accessories

Name	Reference number
Surgical technique Burch-Schneider Reinforcement Cage	50000295
X-ray template Burch-Schneider Reinforcement Cage	367-1485
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance is required with the appropriate hygiene regulations during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must not therefore come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. Therefore, damage after explantation must be assumed which excludes any reuse.

2.2 Authorised component combinations

We only guarantee compatibility of our products in combination with our own CE-marked products and with the products we have approved for combined use, for which the competent authority has issued its approval. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

2.3 Information for use

Before inserting the implant, the implant bed must be irrigated. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed. The porous coating (TPS, Bonit®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

Prepare the implant bed by removing any necrotic tissue and osteophytes near the edge of the acetabulum and then gradually ream the acetabulum using acetabulum reamers in ascending order.

Make sure to always let the reamer move freely and never apply force. The last reamer used has to be 6 mm larger than the size of the reinforcement cage to be implanted.

Using the trial reinforcement cage, determine the size, alignment and bending requirement of the reinforcement cage and mark the caudal tab in the ischium.

After that, open up the passage for the caudal tab in the ischium and fill up any defects in the acetabulum with bony mass from the reaming of the acetabulum. Do not fill up defects of the acetabulum with bone cement.

Drive the caudal tab of the reinforcement cage into the prepared ischium using the insertion ram until the implant sits on the base of the cup. Screw the cranial tab and the cage of the reinforcement cage into the ilium and fill up any defects in the base of the cup with cancellous bone. Alternatively, the caudal tab can also be screwed into the ischium. If necessary, the tabs can be adjusted with the aid of the bending instruments. In doing so, make sure to only perform the bending in one direction and that a maximum total bending angle of 30° is permissible. Once the adaptation is completed, the reinforcement cage has to be optically inspected for notches in the region of the bend.

Continue the surgical procedure by inserting the polyethylene cup. Select the cup according to the size description of the reinforcement cage and cement it with an inclination of 40° and anteversion of 10-15°. A cement coating of 1 mm thickness is thus created. If a thicker cement coating is desired, the next smaller polyethylene cup can be used. Continue the surgical procedure according to the surgical technique of the components being used.

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar is indispensable and provides important information about suitable implants, placing and possible component combinations and enables the size of the implant to be used to be preselected. Surgery may only be performed once it has been established that the patient is able to tolerate the implant material. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1,15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available should another size be required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

The Burch-Schneider reinforcement cage is indicated in cases of total hip replacement if the remaining acetabular bone bed is insufficient and the product can provide support, while protecting the bone graft during the bone resorption phase. This generally applies to:

- contained acetabular defects affecting the entire acetabulum where contact with 50% of the host bone is not possible, even if further reaming is undertaken (with bone graft).
- uncontained acetabular defects affecting > 50% of the acetabulum (with bone graft).
- if acetabular roof reinforcement cages do not come into contact with the inferomedial aspect of the host acetabulum.
- acetabular defects classified as AAOS types III and IV, Paprosky types IIIA and IIIB, or Gross type IV.

Underlying diagnoses/causes may be:

- Avascular or tumour-related necrosis
- Degenerative, post-traumatic or rheumatoid arthritis
- Degenerative, post-traumatic osteoarthritis
- Osteoporotic/low-energy trauma
- Acetabular fractures
- Acetabular dysplasia
- Osteoporosis
- Necrosis of the femur
- Malignancy

6. Contraindications

- Large segmental defects that do not provide sufficient support for the support cup.
- Hypersensitivity to the materials used
- Insufficient bone substance or quality that does not allow for stable fixation of the implant
- Severe malformations
- Local bone tumours
- Acute or chronic, local or systemic infection
- Systemic diseases and metabolic disorders
- Expected overload of the implant e.g. due to overweight, heavy physical activities

7. Conditions that can impair the success of the operation

Potential risks in connection with the operation are:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Growth in children and adolescents
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

8. Possible negative effects

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Muscle spasm
- Stiffness,
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Inflammation
- Oedem/fluids
- Metallosis
- Elevated metal ions in blood
- Coxa vara
- Osteolysis
- Heterotopic ossification
- Pseudotumours
- Cardiac arrhythmia, increased pulmonary vascular resistance, cardiac arrest (caused by BCIS = bone cement implantation syndrome)

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, although the contraindications described above partially apply to the patient, it is particularly important to point out to the patient the effects of these circumstances on the success. It must be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. The patient has to be informed about activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adhesive

labels are enclosed for documenting the implant used. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 133).

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