OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

DE

Tel.: +49 (0)3385 - 54 20-0

Fax: +49 (0)3385 - 54 20-99

www.ohst.de

CE 0482

ELEC®plus Inlay

Deutsch

3

ELEC®plus Inlay

English

10

Inlay ELEC®plus

Français

17

Inserto ELEC®plus

Italiano

24

Inserto ELEC®plus

Español

31

ELEC®plus inlay

Nederlands

38

Interior ELEC®plus

Português

45

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IMPLANT

ELEC®plus Inlay

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance. These implants may be used only by physicians with appropriate experience and practice in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

1. Product description and implant materials

The ELEC®plus Inlay is a ceramic inlay for clamping in metallic acetabular cups and is made from high-purity implant ceramic (ISO 6474-2 Al₂O₃ + ZrO₂). The ceramic is bioinert, biocompatible, biostable, mechanically stable and corrosion-resistant. It combines excellent resistance to fatigue, high shock resistance and tensile strength, excellent fracture toughness, an extremely high degree of hardness and there are no ions to penetrate the body.

The ELEC®plus Inlays are offered in the different articulation diameters 28 mm, 32 mm and 36 mm.

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Article number
ELEC®plus Inlay Ø28/37-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-012
ELEC®plus Inlay Ø28/39-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-013
ELEC®plus Inlay Ø28/44-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-014
ELEC®plus Inlay Ø28/48-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-015
ELEC®plus Inlay Ø28/52-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-016
ELEC®plus Inlay Ø32/39-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-017
ELEC®plus Inlay Ø32/44-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-018
ELEC®plus Inlay Ø32/48-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-019
ELEC®plus Inlay Ø32/52-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-020
ELEC®plus Inlay Ø36/44-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-021
ELEC®plus Inlay Ø36/48-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-022
ELEC®plus Inlay Ø36/52-18	ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - ceramic	013-023

1.2 Instrument overview

The instruments of the OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Article number
Insertion instrument, curved	510-000
Insertion aid Ø37 XLW 18	510-004
Insertion aid Ø39 XLW 18	510-006
Insertion aid Ø44 XLW 18 / XLW 18 zero	510-010
Insertion aid Ø48 XLW 18 / XLW 18 zero	510-012
Insertion aid Ø52 XLW 18 / XLW 18 zero	510-014
Connector Ø28	510-100
Connector Ø32	510-200
Connector Ø36	510-300

1.3 Accessories

Name	Article number
Instructions for use of insertion instrument	200815300886
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance is required with the appropriate hygiene regulations during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must not therefore come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. Therefore, damage after explantation must be assumed which excludes any reuse.

2.2 Authorised component combinations

We guarantee compatibility of our products only in combination with our own CE-marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

Caution: The ELEC®plus Inlays may only be used in combination with ELEC® or ELEC®plus Femoral Heads or with ELEC®plus Revision Femoral Heads. The combination with a different femoral head or with a different ceramic ball of other manufacturers is prohibited.

2.3 Information for use

The clamping cone connection allows for a secure fixation of the inlay inside the implanted hip cup. It ensures a uniform transmission of power into the cup. The ELEC®plus Inlays must fit the implanted cup perfectly. Especially when using press-fit fixation systems, minor deformation of the cup following implantation into the osseous bed can occur, which can in turn make it more difficult to achieve the correct fixation of an ELEC®plus Inlay.

Caution: The ball diameter of the prosthetic femoral head must be absolutely identical to the spherical nominal diameter of the ELEC®plus Inlay.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Before inserting the implant, the implant bed must be irrigated. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

Attaching the inlay in the implanted hip cup:

To ensure reliable function, utmost care is absolutely indispensable. Note the following points before attaching the ELEC®plus Inlay:

1. Ensure an inclination (40°–45°) and anteversion (10°–20°) of the hip cup in line with the indication
2. Thoroughly rinse the clamping cone of the implanted hip cup with water to remove any tissue fragments, bone splinters or residual cement
3. Dry the clamping cone of the implanted hip cup and check the position and function of the joint by inserting a trial implant; make sure that any screws that may have been used have been screwed in sufficiently deep
4. Remove the trial implant and rinse and dry it again
5. Grab one part of the rim of the ELEC®plus Inlay using your thumb and index finger and carefully allow it to slide into the clamping cone of the hip cup; do not apply any additional pressure in this state
6. Check the correct flush position of the ELEC®plus Inlay by running a finger along the rim of the cup. The edge of the face of the cup and of the face of the ELEC®plus Inlay must be flush; if necessary, correct any uneven areas by carefully applying targeted pressure
7. Once this is done, firmly press the ELEC®plus Inlay centrally into the dome using your thumb and wedge it into the clamping cone. Again check that it is positioned flush by running your finger along the rim and - where applicable - reposition it using the instrument authorised by the manufacturer of the endoprosthesis system

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar is indispensable and provides important information about suitable implants, placing and possible component combinations and enables the size of the implant to be used to be preselected. Surgery may only be performed once it has been established that the patient is able to tolerate the implant material. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1,15:1. In addition, X-ray

templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available should another size be required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

- Total hip replacement, primary surgery in combination with prosthetic stems and cups
- Revision (= removal of the endoprosthesis) with new prosthetic stems and cups

6. Contraindications

- Revision involving leaving the stem and/or hip cup in situ (risk of component fracture due to damaged fixation cones)
- Inclination outside 40°–45°, anteversion outside 10°–20° (risk of dislocation, subluxation, impingement)
- Retroversion

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Growth in children and adolescents
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

8. Possible negative effects

The negative effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedem/fluids

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Muscle spasm
- Stiffness
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Metallosis
- Elevated metal ions in blood
- Coxa vara
- Osteolysis
- Heterotopic ossification
- Pseudotumours

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, although the contraindications described above partially apply to the patient, it is particularly important to point out to the patient the effects of these circumstances on the success. It must be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. The patient has to be informed about activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adhesive labels are enclosed for documenting the implant used. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 112).

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