OHST

medical technology

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IMPLANT

Expersus Femoral Stem

Before using the product, the user is under obligation to carefully study and follow the recommendations and information below, together with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance.

These implants may be used only by physicians with appropriate experience, practice and skills in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website www.ohst.de/ifu/.

1. Product description and implant materials

The Expersus femoral stem serves to restore the anatomical conditions of the hip joint in primary hip replacement to relieve pain and restore functionality. The Expersus femoral stem is a prosthesis for cementless and cemented anchorage in the femur.

The cementless variant is manufactured from titanium alloy Ti6Al4V (ISO 5832-3) and the cemented variant from implant steel (ISO 5832-9). The Expersus femoral stem comes with a 12/14 cone. The stem can be combined with both metal and ceramic femoral heads. Apart from the polished neck area, the cementless Expersus femoral stem has an all-round titanium plasma spray and hydroxyapatite coating. The body of the cemented femoral stem has a satin-finished, lightly blasted surface.

To account for the various femoral anatomies of the patients, the femoral stem comes in 9 sizes with a standard 125° and a standard 135° variant for cementless anchorage as well as in 9 sizes with a standard 125° and in 10 sizes with a standard 135° variant for cemented anchorage.

The Expersus femoral stem can generally be used in patients in whom total hip replacement surgery is indicated. The Expersus femoral stem may only be used to treat skeletally mature patients.

All implant variants of the Expersus femoral stem can be inserted with the same instrument set.

The product, the package contents and the materials used are specified on the labels. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Reference number
Expersus femoral stem cemented		
Expersus femoral stem size 1 STD, 135°, cemented (*)	ISO 5832-9 Implant steel	367-1556
Expersus femoral stem size 2 STD, 135°, cemented (*)	ISO 5832-9 Implant steel	367-1557
Expersus femoral stem size 3 STD, 135°, cemented (**)	ISO 5832-9 Implant steel	367-1558
Expersus femoral stem size 4 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1559
Expersus femoral stem size 5 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1560
Expersus femoral stem size 6 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1561
Expersus femoral stem size 7 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1562
Expersus femoral stem size 8 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1563
Expersus femoral stem size 9 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1564
Expersus femoral stem size 10 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1565
Expersus femoral stem size 2 STD, 125°, cemented (*)	ISO 5832-9 Implant steel	367-1547
Expersus femoral stem size 3 STD, 125°, cemented (**)	ISO 5832-9 Implant steel	367-1548
Expersus femoral stem size 4 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1549
Expersus femoral stem size 5 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1550
Expersus femoral stem size 6 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1551
Expersus femoral stem size 7 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1552
Expersus femoral stem size 8 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1553
Expersus femoral stem size 9 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1554
Expersus femoral stem size 10 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1555

(*) Weight limited to 65 kg (see Contraindications)

(**) Weight limited to 90 kg (see Contraindications)

Expersus femoral stem cementless		
Expersus femoral stem size 2 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1400
Expersus femoral stem size 3 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1401
Expersus femoral stem size 4 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1402
Expersus femoral stem size 5 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1403
Expersus femoral stem size 6 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1404
Expersus femoral stem size 7 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1405
Expersus femoral stem size 8 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1406
Expersus femoral stem size 9 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1407
Expersus femoral stem size 10 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1408
Expersus femoral stem size 2 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1409
Expersus femoral stem size 3 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1410
Expersus femoral stem size 4 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1411
Expersus femoral stem size 5 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1412
Expersus femoral stem size 6 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1413
Expersus femoral stem size 7 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1414
Expersus femoral stem size 8 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1415
Expersus femoral stem size 9 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1416
Expersus femoral stem size 10 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1417

Cement restrictor		
Cement restrictor size 1, external Ø24mm	ISO 5834-2 UHMWPE	506-101
Cement restrictor size 2, external Ø27mm	ISO 5834-2 UHMWPE	506-102

1.2 Instrument overview

The instruments of OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Reference number
Expersus femoral stem instrument set	367-1493
Expersus femoral stem rasp instrument set	367-1492
Femoral Stem Extraction Instrument Set	206-010

1.3 Accessories

Name	Reference number
Surgical technique Expersus femoral stem cementless	50000351
Surgical technique Expersus femoral stem cemented	50000352
X-ray template Expersus femoral stem cementless STD 135° KD28	367-246
X-ray template Expersus femoral stem cementless STD 125° KD28	367-247
X-ray template Expersus femoral stem cementless STD 135° KD32	367-074
X-ray template Expersus femoral stem cementless STD 125° KD32	367-075
X-ray template Expersus femoral stem cemented STD 135° KD28	367-244

Name	Reference number
X-ray template Expersus femoral stem cemented STD 125° KD28	367-245
X-ray template Expersus femoral stem cemented STD 135° KD32	367-072
X-ray template Expersus femoral stem cemented STD 125° KD32	367-073
Instructions for Use Extraction Instrument Set	50000428
Instructions for use stem reducer	50000427
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website www.ohst.de/ifu/, must be observed.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance with the appropriate hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must therefore not come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on the functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. After explantation, previous damages which exclude any reuse must therefore be assumed.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint. The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

2.2 Authorised component combinations

The following femoral heads must be used for assembly with the implant:

Nominal Ø	Name	Neck length	Reference number
Ø22 mm	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L	-4 / 0 / +4 mm	030-2200 to 030-2202
Ø28 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL	-3,5 / 0 / +3,5 / +7 mm	020-2800 to 020-2803
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL	-3,5 / 0 / +3,5 / +7 mm	030-2800 to 030-2803
	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-907 to 367-909
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-1140 to 367-1142
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	384-001 to 384-003
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	013-001 to 013-003
Ø32 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	020-3200 to 020-3203
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	030-3200 to 030-3203
	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	367-920 to 367-912
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	367-1143 to 367-1145, 367-1149
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	384-004 to 384-006
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	013-004 to 013-007

Nominal Ø	Name	Neck length	Reference number
Ø36 mm	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	367-930 to 367-932
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	367-1146 to 367-1148, 367-1150
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	384-007 to 384-009
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	013-008 to 013-011
Ø40 mm to Ø60 mm	Fracture head size S implant steel ISO 5832-9, 12/14	-4 mm	155-140 to 155-160
	Fracture head size M implant steel ISO 5832-9, 12/14	0 mm	155-040 to 155-060
	Fracture head size L implant steel ISO 5832-9, 12/14	+4 mm	155-240 to 155-260

We guarantee compatibility of our products only in combination with our own CE-marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

2.3 Information for use

The application of the implant is done with or without cement. A cement restrictor must be used for cemented applications. Before using the cement restrictor, the most up-to-date version of the respective instructions for use, which is available for download on the website www.ohst.de/ifu/, must be observed. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads, the most up-to-date version of the respective instructions for use, which is available for download on the website www.ohst.de/ifu/, must also be observed. After reduction, the definitive stability, mobility and muscle tension must be checked.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Caution: The patient weight limitation stated in the contraindications must be observed for the Expersus femoral stem sizes listed under 1.1.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e.g., bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

The porous coating (TPS, Bonit®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

After opening the joint capsule and dislocating the femoral head from the cup, the femoral head should be resected as determined during the preoperative planning and completely removed.

A femoral neck resection guide is optionally available for marking the resection plane. Then open the medullary cavity using the cavity chisel. Apply the chisel to the far lateral and dorsal side in order to facilitate the subsequent driving in of the rasps in the direction of the femoral axis. The cavity chisel should be positioned relative to the desired antetorsion. Avoid fracturing the greater trochanter. The medullary cavity opening can be enlarged with the surgical awl.

Once the medullary cavity has been opened, rasp the femoral stem. The first rasp determines the alignment of the subsequent sizes. For this reason, take care to ensure that the antetorsion is already correct during the first rasping step. The antetorsion can be checked using the guide rod and rasp handle and is usually 10°-15°. Start the procedure with the smallest rasp size, which is connected to the rasp handle. After that, use the subsequent rasps in ascending order until the preoperatively determined size is reached. The sizes of the rasps match the sizes of the implants. Verify the correct position of the rasp inside the femur with an image converter. If during the rasping procedure the stem size does not correspond to the preoperatively determined size and there is a difference of two or more sizes, the alignment of the axis may be incorrect or there may

be an osseous obstruction present. In this case, the selected stem might be too small and would therefore not provide the required stability. On request, special rasp handles are available for the various access routes. The rasps are designed in such a way that they can also be used for trial reductions. The leg length, the range of motion and the tension of the ligaments can be checked by attaching the various trial cones and trial heads. Remove the handle from the rasp, which remains inside the femur, and attach the trial cone to the rasp. Trial cones with two offsets are offered for trial reductions of the various stem versions. The trial cones snap into place as soon as they are in their correct position. Then attach the trial head to the trial cone by hand. Trial heads with various diameters are available for trial reductions in the neck lengths S to XL. After reduction, the definitive stability, mobility and muscle tension should be checked.

With the cemented Expersus femoral stem, the cement restrictor is implanted first and the cement inserted afterwards; this surgical step is not applicable to the cementless variant of the Expersus femoral stem.

The cement restrictor is used depending on the preoperatively determined diaphyseal diameter of the medullary cavity. A special insertion instrument is available for controlled insertion to the optimal depth. The size 1 cement restrictor is recommended for medullary cavity diameters between 14 mm and 19 mm, or size 2 for diameters between 18 mm and 22 mm. Screw the cement restrictor onto the threaded rod of the insertion instrument and insert it into the medullary cavity. The insertion depth can be read off the scale on the threaded rod and should exceed the stem length of the prosthetic femoral stem to be implanted by at least 5 mm, but never more than 20 mm. It is also possible to gauge the position of the cement restrictor inside the femoral stem by means of X-ray control based on the position of the radiopaque ring. Thereafter, remove the insertion instrument by rotating it anti-clockwise. After inserting the cement restrictor, rinse and dry the implant bed. The cement restrictor should be implanted before the bone cement is mixed. The cement is to be inserted using a modern cementing technique; in doing so, the instructions for use of the cement manufacturer must be observed.

Insert the prosthetic device into the implant bed up to the resection plane using the stem impactor. When using the cemented variant, complete this step immediately after injecting the bone cement and gently apply pressure while waiting for the bone cement to harden. Then thoroughly clean and dry the prosthetic cone and attach the femoral head with the previously determined neck length (S to XL).

Caution: Do not use femoral heads with a neck length greater than XL!

Then reduce the stem with the femoral head into the cup and check the range of motion, leg length and tension of the ligaments. The operation is routinely completed with layer-by-layer wound closure.

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website www.ohst.de/ifu/. The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar data is indispensable and provides important information about suitable implants, placement, possible component combinations and enables the size of the implant to be used to be preselected. The operation may be performed only if the patient's material compatibility has been determined. Use the X-ray templates for planning the operation. In addition, the patient weight limitation stated in the contraindications must be observed for the Expersus femoral stem sizes named under 1.1. These are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available in case another size is required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

- Advanced degeneration of the hip joint due to degenerative, post-traumatic rheumatoid arthritis or congenital hip dysplasia
- Fracture (e.g. femoral neck fracture) or avascular necrosis of the femoral head
- Sequelae of earlier surgical procedures, e.g. internal fixation, articular reconstruction, arthrodesis
- Hemiarthroplasty or total hip prosthesis
- Certain cases of ankylosis

Total hip replacements are not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, they can restore mobility whilst simultaneously relieving pain. It should be the course of action only where all other treatment approaches, which could preserve the joint, whether surgically or conservatively, and which are deemed to be timely and appropriate in medical terms, have failed.

6. Contraindications

- Acute or chronic local or systemic infection
- Severe muscle, nerve or vascular diseases that endanger the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Any underlying condition that might compromise the function of the implant
- Hypersensitivity to the materials used
- Patients weighing more than 65 kg for the cemented variant of the Expersus femoral stem size 1 STD 135°, size 2 STD 135° and size 2 STD 125° (REF 367-1556, 367-1557 and 367-1547)
- Patients weighing more than 90 kg for the cemented variant of the Expersus femoral stem size 3 STD 135° and size 3 STD 125° (REF 367-1558 and 367-1548)

If femoral heads with a neck (XL and XXL) are used, the range of motion is reduced by about 30° and achieves values between 80° and 100° in flexion and extension.

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Growth in children and adolescents
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

8. Possible negative effects

The negative effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedem/fluids

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Stiffness
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Metallosis
- Elevated metal ions in blood
- Osteolysis
- Heterotopic ossification
- Pseudotumours
- Cardiac arrhythmia, increased pulmonary vascular resistance, cardiac arrest (caused by BCIS = bone cement implantation syndrome)
- Deep infections

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient information, documentation

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, even though the contraindications described above partially apply to the patient, the patients should be informed about the possible effects of these circumstances and the anticipated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. Based on the current state of technology, the expected service life is approx. 10 to 15 years. The patient must be informed about activities and, if applicable, the weight limitation with which the effects of these aggravating circumstances can be minimised.

It is also recommended that the patient is informed of activities with which he can reduce the effects of these aggravating circumstances. All the information given to the patient must be documented in writing by the surgeon. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

A short report about the safety and clinical performance is available from the Eudamed database. Until the database is launched, the short report can be provided on request.

10. Implant passport

Following surgery, the patient must be given an implant passport containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant passport must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website. The patient's details (patient name or patient ID) must be documented in the implant passport, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the passport.

The user must inform the patient that additional or updated information to ensure the safe use of the product by the patient is available on the indicated website.

11. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 191).

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