OHST

medical technology

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IMPLANT

BIOLOX[®]delta Inlay

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance. These implants may be used only by physicians with appropriate experience and practice in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result. Care must always be taken to use the most up-to-date version of the instructions for use provided with the system components. The date of revision is shown on the cover sheet in the footer.

1. Product description and implant materials

The BIOLOX[®]delta inlay is a ceramic inlay for clamping in metallic acetabular cups. It is made from high-purity aluminium oxide ceramic stabilised with zirconia (ISO 6474-2 Al₂O₃ + ZrO₂).

The BIOLOX[®]delta inlays are offered in the various articulation diameters 28 mm, 32 mm and 36 mm.

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Reference number
BIOLOX [®] delta inlay XLW 18 28/37G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-894
BIOLOX [®] delta inlay XLW 18 28/39G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-1307
BIOLOX [®] delta inlay XLW 18 28/44G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-1308
BIOLOX [®] delta inlay XLW 18 28/48G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-1309
BIOLOX [®] delta inlay XLW 18 28/52G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-1310
BIOLOX [®] delta inlay XLW 18 32/39G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-895
BIOLOX [®] delta inlay XLW 18 32/44G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-896
BIOLOX [®] delta inlay XLW 18 32/48G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-897
BIOLOX [®] delta inlay XLW 18 32/52G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-1311
BIOLOX [®] delta inlay XLW 18 36/44G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-898
BIOLOX [®] delta inlay XLW 18 36/48G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-899
BIOLOX [®] delta inlay XLW 18 36/52G	ISO 6474-2 Al ₂ O ₃ + ZrO ₂	367-1312

1.2 Instrument overview

The instruments of the OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Reference number
Liner Handle (curved)	510-000

Name	Reference number
Mounting Device Ø37 XLW 18	510-004
Mounting Device Ø39 XLW 18	510-006
Mounting Device Ø44 XLW 18 / XLW 18 zero	510-010
Mounting Device Ø48 XLW 18 / XLW 18 zero	510-012
Mounting Device Ø52 XLW 18 / XLW 18 zero	510-014
Connector Ø28	510-100
Connector Ø32	510-200
Connector Ø36	510-300

1.3 Accessories

Name	Reference number
Instructions for Use of insertion instrument	200815300886
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance is required with the appropriate hygiene regulations during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must not therefore come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. Therefore, damage after explantation must be assumed which excludes any reuse.

2.2 Authorised component combinations

We guarantee compatibility of our products only in combination with our own CE-marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

Caution: The use of prosthetic components that have not been approved by the endoprosthesis manufacturer in combination with a BIOLOX®delta inlay may lead to fracture of the BIOLOX®delta inlay.

Caution: The BIOLOX®delta inlays may only be used in combination with BIOLOX®delta –or BIOLOX®forte femoral heads. The combination with a different femoral head or with a different ceramic ball of other manufacturers is prohibited.

2.3 Information for use

Care must be taken to use only an unused and intact BIOLOX®delta inlay which was removed from its original packaging immediately prior to insertion. A BIOLOX®delta inlay that has been used before must not be reused. Due to the need for an exact fit of the BIOLOX®delta inlay with the cup housing, only new and intact prosthetic components may be combined. This also means, for example, that a BIOLOX®delta inlay which was already inserted into a cup housing before and was then removed again must not be reinserted into a cup housing. By the same token, a BIOLOX®delta inlay that is in any way damaged must not be reused, but should be discarded. This also applies, for example, to a BIOLOX®delta inlay that has been dropped.

A BIOLOX®delta inlay must only be used with a brand-new, unused cup housing that is free from damage.

Caution: The ball diameter of the femoral head prosthesis must be absolutely identical to the spherical nominal diameter of the BIOLOX®delta inlay.

Caution: The positioning of the implant components has a direct influence on the range of motion, and therefore potentially on the risk of impingement. An angle of inclination that is either too steep or too flat increases the load on the edge of the cup. This can lead to increased wear. In addition, movements may be restricted, which in turn may lead to subluxation and dislocation of the head from the BIOLOX®delta inlay.

Caution: The anteversion of the cup housing should not display any negative values. This can result in increased wear. In addition, movements may be restricted, which in turn may lead to subluxation and dislocation of the head from the BIOLOX®delta inlay.

Caution: Do not use a BIOLOX®delta inlay if the cup housing is in a position which is not recommended in the two paragraphs above (exceptions may be special systems intended for this purpose with elements to prevent dislocation, e.g. asymmetrical or tilted inserts).

Caution: In the case of retroverted cup housings, no BioloX®delta inlay should be used. Otherwise, failure of the mechanical function of the BIOLOX®delta inlay may occur. Excessive wear of the ceramic may result in intolerances or loosening of the prosthesis.

Caution: Revision using a BIOLOX®delta inlay with an endoprosthesis stem that remains in situ is contraindicated. In this case, as a rule, an inlay made of plastic can generally be used.

- Caution:** In very rare cases, the BIOLOX®delta inlay may break. To minimise the risk as much as possible, every individual component of the BIOLOX®delta inlay has been inspected prior to delivery. Fractures may be caused, among other things, by incorrect fixing the BIOLOX®delta inlay inside the metallic cup housing or by an incorrect fit (or no fit at all) between the BIOLOX®delta inlay and the metallic hip cup. The use of prosthetic components that have not been approved by the manufacturer of the endoprosthesis system in combination with a BIOLOX®delta inlay may also lead to a fracture of the BIOLOX®delta inlay. The same applies if the recommended positioning of the cup housing is disregarded. In the very rare event of a fracture of a ceramic component, synovectomy may have to be carried out.
- Caution:** In the very rare event that a ceramic insert breaks, there is a risk of injury from the sharp edges of ceramic fragments during revision surgery
- Caution:** Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing. Furthermore, a pairing of metal (femoral head) to plastic (inlay) as well as metal to metal is contraindicated during revision surgery.
- Caution:** If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Any noise during movement after implanting a BIOLOX®delta inlay is not a sufficient indication of a malfunction or change in the performance of the endoprosthesis system. It is nevertheless recommended to check the intactness of the endoprosthesis system.

Before inserting the implant, the implant bed must be irrigated. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

- Caution:** When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

As a vital prerequisite for safe operation, maximum care must be taken when inserting the BIOLOX®delta inlay into the pan housing.

In order to ensure central mounting of the BIOLOX®delta inlay in the metal acetabular cup, the following procedure should be adopted during installation:

The metal acetabular cup should be implanted into the acetabulum in accordance with the surgical technique. Trial implants should then be used to verify the correct fit of the implant and the movements of the joint. It must be ensured that any screws used are completely recessed in the metal cup.

Afterwards, the trial implant must be removed and the metal cup rinsed and dried thoroughly. There must be no bone or tissue residues remaining in the clamping area of the metal cup.

The BIOLOX®delta inlay is inserted by hand or with the insertion instrument. Whenever the insertion instrument is used, the respective Instructions for Use (see 1.3 Accessories) must be complied with.

When inserted by hand, the inlay is held with two fingers and inserted into the cup housing. The positioning of the BIOLOX®delta inlay is checked by feeling the edge of the cup and correcting it if necessary. The connection between the edge of the metal and edge of the ceramic must be flush!

Once the BIOLOX®delta inlay is positioned correctly, press down with the thumb on the base of the insert. To fixate the implant in its final position, tightly clamp the BIOLOX®delta inlay into the cup housing with one or more moderate hammer blows on the impactor in the direction of the axis of symmetry using an impactor that is suitable for BIOLOX®delta inlays and recommended by the endoprosthesis manufacturer.

Caution: Never use a metal hammer to hit the BIOLOX®delta inlay!

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar is indispensable and provides important information about suitable implants, placing and possible component combinations and enables the size of the implant to be used to be preselected. Surgery may only be performed once it has been established that the patient is able to tolerate the implant material. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1,15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available should another size be required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

As a general rule, all indications and recommendations that are known for total hip replacement with other tribological pairings and that are accepted as the generally recognised standard apply when using a BIOLOX®delta inlay.

In addition, the indications and recommendations for the prosthetic system being used, of which BIOLOX®delta inlays are integral components, are to be taken into account.

6. Contraindications

As a general rule, all contraindications and recommendations that are known for total hip replacement with other tribological pairings and that are accepted as the generally recognised standard apply when using a BIOLOX®delta inlay.

- Growth in children and adolescents

In addition, the contraindications and recommendations for the prosthetic system being used, of which BIOLOX®delta inlays are integral components, are to be taken into account.

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

8. Possible negative effects

The negative effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedem/fluids

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Stiffness
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Metallosis
- Elevated metal ions in blood
- Osteolysis
- Heterotopic ossification
- Pseudotumours

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, although the contraindications described above partially apply to the patient, it is particularly important to point out to the patient the effects of these circumstances on the success. It must be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. The patient has to be informed about activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adhesive labels are enclosed for documenting the implant used. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 131).

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