**OHST**   
medical technology

OHST Medizintechnik AG  
Grünauer Fenn 3  
14712 Rathenow  
Germany  
Tel.: +49 (0)3385 - 54 20-0  
Fax: +49 (0)3385 - 54 20-99  
[www.ohst.de](http://www.ohst.de)

**CE** 0482

---

EFK Kniesystem Deutsch (de)	3
EFK Knee System English (en)	17
EFK Système de genou Français (fr)	31
EFK Sistema per ginocchio Italiano (it)	45
Sistema de rodilla EFK Español (es)	59

## IMPLANTAT

### EFK Kniesystem

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## IMPLANT

### EFK Knee System

Before using the product, the user is under obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for knee arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website [www.ohst.de/fu/](http://www.ohst.de/fu/).

#### 1. Product Description and Implant Materials

Knee prostheses are implants that remain in the body permanently and either partially or fully replace the damaged knee. The knee joint in its entirety is replaced by three components: the femoral component, the tibial component and the patellar component.

Irrespective of the shape and fixation manner of the individual joint components, they always perform the same functions (ISO 7207-1, ISO 21536):

**Femoral component:** the component of a total knee joint replacement intended to be secured to the femur to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.

**Tibial component:** the component of a total knee joint replacement intended to be secured to the tibia to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.

**Patellar component:** the component of a total or partial knee joint replacement which is used to replace the articulating surface of the patella. These implants can be manufactured as one component or a set of components to be assembled by the user.

The EFK Knee System is a tricompartmental surface replacement for cemented, cementless or hybrid use within the scope of primary or revision knee replacement surgery. The implant system enables the posterior cruciate ligament to be preserved or requires the anterior and posterior cruciate ligament to be removed.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

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## 1.2 Overview of Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
Instrument Set EFK Tibia	10006776
Instrument Set EFK Femur Sz. 4-10	10006778
Instrument Set EFK Trials Sz. 4-10	10006779
Instrument Set EFK Additional Sz. 2+12	10006780
Instrument Set EFK Trial Inserts UC	10006781
Instrument Set Patella	10006821
Instrument Set EFK Optional	10007371
Instrument Set EFK Tibia Modular	10007372

## 1.3 Miscellaneous Accessories

Designation	Reference number
X-ray templates EFK Femur	100000510296
X-ray templates EFK Tibia+Insert	100000510297
X-ray templates EFK Tibia Modular	100000510319
X-ray templates EFK S Femur	100000511064
Surgical technique EFK Knee System	50000673
Surgical technique EFK Tibia Modular	50000674
Implant passport	50000572

## 2. Handling

### 2.1 General information

This implant is part of a system and must only be used with the corresponding original system components.

Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website [www.ohst.de/ifu/](http://www.ohst.de/ifu/), must be observed.

**Caution:** implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise.

Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

**Caution:** the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explanation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding side is marked on the implants with "L" for the left knee and "R" for the right knee. The orientation of the implants must always correspond to the side of the body of the knee to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right knee.

The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

## **2.2** Authorised Component Combinations

Each of the tibial components must be combined with inserts of the same size.

The inserts can be combined with femoral components of the same size as well as with the next larger and the next smaller size. This applies to both the inserts and the UC inserts. The patella component can be combined with all sizes of the femoral component.

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

## **2.3** Instructions for Use

Before the instruments are used, the enclosed product information and instructions for use must be observed.

The surfaces with a porous coating (TPS, BONIT®) and the roughened implant surfaces must not come in contact with clothing or other fibre-shedding materials.

**Caution:** Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free latex gloves.

Follow the manufacturer's instructions on using the bone cement to ensure correct performance of the cementing step. The use of bone cement mixed under vacuum is recommended.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed. In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

**Caution:** Inserts UC shall only be used when there is no intact posterior cruciate ligament, because there is no room for it.

**Caution:** When using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be

replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

## 2.4 Surgical Technique

Sequence of resections:

1. Distal femoral resection
2. A/P-femoral resections
3. Facet resections (including ablation of dorsal osteophytes)
4. Tibial resection (can be performed alternatively as resection 1)
5. Patellar resection (optionally)

It is important that the flexion gap and extension gap are approximately the same size.

Approach:

For information on accessing the knee joint, we refer to the relevant surgical procedures.

Preparation of the femur:

1. Distal femoral opening
2. Positioning of the distal femoral resection gauge
  - At its distal end, the femoral implant is 9 mm thick.
  - Bone pins must always be pre-drilled using a 3.2 mm drill.
  - Before covering the patient, it is advisable to apply an electrode to the patient at the height of the femoral head centre using a C-arm X-ray unit, as the position depending on the soft tissue situation and the covering can otherwise not be precisely determined.

**Caution:** Clinical experience has shown that restoration of a physiological leg axis and a correct alignment of the implants have an influence on the service life of the knee replacement. Poor alignment can lead to non-physiological stress on the implants and can therefore shorten the service life.

3. Distal femoral resection
4. Determination of the femoral size and rotation

**Caution:** Clinical experience has shown that the correct size selection influences the service life of the knee replacement. The poorer the coverage of the prosthesis parts on the bone, the shorter the service life. A lack of coverage or an excessive implant overhang of  $\geq 2$  mm must therefore be avoided.

- If you choose the next larger femoral size, this will increase the posterior offset and reduce the flexion gap. If you select the smaller size, this will increase the flexion gap and reduce the posterior offset. A reduced posterior offset can diminish the flexion ability.

5. A/P resection, facet resection

- If you perform the upper facet cut before the ventral cut, this will reduce the risk of the saw blade "wandering off", in particular when working on sclerotic bones.

6. Preparing the trochlear groove

- A lateral position of the femoral component is beneficial for the patella tracking.

- In case the position of the trial component is not stable enough to work on the trochlea, leave the stop drill in one of the drill holes.

#### Preparation of the tibia:

1. Installation of the extramedullary alignment gauge
2. Adjustment of the slope and varus/valgus alignment
3. Adjustment of the rotation
4. Adjustment of the resection height

Alternatively, intramedullary alignment: opening of the tibia, adjustment of rotation, resection height and slope

5. Axis check

**Caution:** Clinical experience has shown that restoration of a physiological leg axis and a correct alignment of the implants have an influence on the service life of the knee replacement. Poor alignment can lead to non-physiological stress on the implants and can therefore shorten the service life.

#### 6. Resection of the tibia

**Caution:** When displacing the revision cutting gauge to the distal or proximal direction (for thicker or thinner bone resections), this is done either to the right (+/- 2 mm) or parallel (+/- 4 mm). In the case of a left-hand shift (approx. pin width), make sure that the saw blade is applied to the outside of the right pin.

#### 7. Determination of the tibial component size

**Caution:** Clinical experience has shown that the correct size selection influences the service life of the knee replacement. The poorer the coverage of the prosthesis parts on the bone, the shorter the service life. A lack of coverage or an excessive implant overhang of  $\geq 2$  mm must therefore be avoided.

#### 8. Trial run

- After adjusting the rotation, use a cautery knife or a sterile pen to mark the ventral tibia just below the two anterior grooves of the tibial size gauge. This makes it easier to retrieve the correct rotation position for the tibial size gauge after removing the trial femur and the trial insert.

9. Preparation of the medullary cavity of the tibia for the stem

#### Patella preparation:

1. Resection

- During resecting, make sure that the saw blade does not migrate, e. g. due to sclerotic bone segments.

2. Patella size determination

3. Milling

4. Drill anchoring holes

5. Manipulating

#### Implantation:

1. Insert the tibial component and impact by means of the impactor

2. Introduce the insert and click it ventrally

3. Fix the insert by tightening the screw using the screwdriver

4. Insert the femoral component and impact by means of the impactor

- Caution: In combination with original implants, a cleaning of the trial insert is urgently necessary in order to avoid damages and scratches on the implant surfaces by third bodies.
- Caution: The sequence of implantation is generally the responsibility of the surgeon.
- Caution: The tibia fixation screw must be tightened with the supplied screwdriver as firmly as possible without any additional tools.
- Caution: Care must be taken not to damage the surfaces. Use the femur/tibia impactor for punching in, put it on centrally and align it parallel to the mechanical axis.
- Caution: It should be noted that an insert removed from the tibial component or jammed when snapped into place is to be visually inspected for damage before re-use. In the event of damage, the parts may no longer be used.

In order to prevent excessive heating of the bone, all resections must be performed while cooling with a physiological saline solution. The thickness of the saw blades must be adjusted to the instruments used! The cement is inserted using a modern cementing technique; in doing so, the instructions for use of the cement manufacturer must be observed.

- Caution: Clinical experience has shown that cementation influences the service life of the knee replacement. Cementation of < 1 mm was assessed as poor, 1-2 mm as suboptimal and 3-4 mm as optimal.

### **3. Packaging and Sterility**

Depending on the sterilisation method used, implants are packaged in a triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website [www.ohst.de/ifu/](http://www.ohst.de/ifu/). The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

- Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

#### **4. Preoperative Planning and Postoperative care**

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognised procedures must be used for postoperative care.

#### **5. Indications**

Total prosthesis is indicated in patients who have lost their mobility for any of the following reasons:

- Degenerative, post-traumatic or rheumatoid arthritis
- Avascular necrosis of the femoral condyle
- Post-traumatic loss of the joint configuration, in particular in the case of patellofemoral erosion, dysfunction or prior patellectomy
- Moderate varus, valgus or flexion malpositions

This system may also be indicated for the revision of earlier surgical attempts which have been unsuccessful.

Total knee replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total knee replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

#### **6. Contraindications**

- Acute or chronic, local or systemic infection
- Severe muscle, nerve or vascular diseases endangering the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Hypersensitivity to the materials used
- Loss of the ligamenture
- Overweight
- Alcoholism or other addictions
- High-level/intense physical activity (e.g. competitive sports, heavy physical labour)
- Any underlying condition which could compromise the function of the implant

#### **7. Risks and conditions that may affect the success of the surgery**

**Caution:** Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of knee arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcoholism or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Intake of high-dose cortisone or cytotoxicics
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to knee arthroplasty:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Growth in children and adolescents
- Expected extreme stresses, e.g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities that complicate fixation of the implant
- Weakening of the supporting structures due to a tumour

## **8.** Undesirable Effects

The undesirable effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematoma
- Paraesthesia
- Numbness
- Swelling
- Nerve damage
- Oedema

The undesirable effects listed below are amongst the most typical and commonly occurring consequences of total knee replacement surgery:

- Tibial or femoral fracture
- Sub-dislocation or dislocation of the patella
- Instability, movement or loosening of components
- Separation of components
- Insufficient freedom of movement
- Metal hypersensitivity reactions
- Patellofemoral complications
- Neurovascular complications
- Paralysis of nervus peroneus
- Fat embolism
- Arterial insufficiency / injury

- Arthrofibrosis / adhesions
- Implant breakage
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Inflammation
- Metallosis
- Elevated metal ions in blood
- Pseudotumors

Caution: Revision surgery may be required if specific undesirable effects occur.

## **9. Patient Information and Documentation**

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks. A service life of at least five years can be expected, even with difficult conditions (e.g. if there are multiple risk factors and comorbidities). Patients who receive a knee replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant. Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

A short report about the safety and clinical performance is available from the Eudamed database. Until the database is launched, the short report can be provided on request.

## **10. Implant passport**

Following surgery, the patient must be given an implant passport containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant passport must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website. The patient's details (patient name or patient ID) must be documented in the implant passport, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the passport.

The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

## **11.** Explanation of the Label Symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 73).

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Puede ver los símbolos utilizados por OHST Medizintechnik AG en el anexo (página 73).