OHST
medical technology 

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IMPLANT

NANOS femoral neck prosthesis

Before using the product, the user is under obligation to carefully study and comply with the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance.

These implants may be used only by physicians with appropriate experience, practice and skills in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

1. Product description and implant materials

The NANOS femoral neck prosthesis for cementless fixation is manufactured from Ti6Al4V forge alloy (ISO 5832-3) and has a proximal osteoconductive titanium plasma coating (TPS) and an additional calcium phosphate coating (BONIT®). The NANOS femoral neck prosthesis enables metaphyseal fixation and load distribution. The implant requires sparing and minimal bone resection. A total of 11 perfectly matched sizes are available.

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Art. no.	SAP no.
NANOS femoral neck prosthesis size 0	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425000	75008154
NANOS femoral neck prosthesis size 1	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425001	75008155
NANOS femoral neck prosthesis size 2	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425002	75008156
NANOS femoral neck prosthesis size 3	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425003	75008157
NANOS femoral neck prosthesis size 4	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425004	75008158
NANOS femoral neck prosthesis size 5	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425005	75008159
NANOS femoral neck prosthesis size 6	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425006	75008160
NANOS femoral neck prosthesis size 7	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425007	75008161

Name	Material	Art. no.	SAP no.
NANOS femoral neck prosthesis size 8	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425008	75008162
NANOS femoral neck prosthesis size 9	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425009	75008163
NANOS femoral neck prosthesis size 10 (*)	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	75103109	75103109

(*) Special size

1.2 Instrument overview

The instruments of OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Art. no.	SAP no.
Instrument Set NANOS Rasp Tray	-	75210420
Instrument Set NANOS Tray	-	75200170

1.3 Accessories

Name	Art. no.	SAP no.
Surgical technique femoral neck prosthesis	50000698	-
X-ray template femoral neck prosthesis D22	12330	-
X-ray template femoral neck prosthesis D28	09777	-
X-ray template femoral neck prosthesis D32	09778	-
X-ray template femoral neck prosthesis D36	09779	-
Instructions for use stem reducer	50000427	-
Implant passport	50000572	-
Rasp handle double offset right hook coupling (*)	367-155	75103113
Rasp handle double offset left hook coupling (*)	367-156	75103114
Rasp adapter IMT double offset right hook coupling (*)	367-158	75103115
Rasp adapter IMT double offset left hook coupling (*)	367-159	75103116
Rasp adapter IMT spring hook NANOS (*)	506-1210	75094887
Rasp handle spring hook trigger mechanism medial NANOS (*)	506-1245	75030476
Rasp handle spring hook straight NANOS (*)	506-1204	75094886
Awl extended with handle, 90° rotation (*)	506-1220	75030323
Trial head Ø22 S/+0 (*)	75103040	75103040
Trial head Ø22 M/+4 (*)	75103041	75103041
Trial head Ø22 L/+8 (*)	75103042	75103042
Trial head Ø22 XL/+12 (*)	75103043	75103043
Trial head Ø28 XS/-3 (*)	75103044	75103044
Trial head Ø32 XS/-3 (*)	75103050	75103050
Trial head Ø36 XS/-3 (*)	75103056	75103056
NANOS compactor size 10 (*)	75103110	75103110
NANOS rasp size 10 (*)	75103111	75103111

(*) optional instruments

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance with the appropriate hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must therefore not come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on the functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. After explantation, previous damages which exclude any reuse must therefore be assumed.

2.2 Authorised component combinations

Only femoral heads that are compatible with the system may be used for the assembly with the implant.

The approved combinations with Smith & Nephew femoral heads are listed in a compatibility matrix. The matrix is available at: www.smith-nephew.com/compatibilitymatrix (Lit.No. 04758).

In general, the guidelines for use and handling as specified in the manufacturer's instructions for use of the endoprotheses must be followed and implemented.

Due to reasons relating the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

2.3 Information for use

The application of the implant is done without cement. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads the corresponding instruction for use is also to be considered. After reduction, the definitive stability, mobility and muscle tension must be checked.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available. Before inserting the implant, the implant bed must be irrigated sufficiently. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

The porous coating (TPS, Bonit®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

Caution: Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free latex gloves.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

Either the anterior-lateral or transgluteal approach can be used. Other approaches that are considered appropriate by the surgeon can also be used. The prosthesis type is also suitable for minimally invasive approaches. Resect the bone as sparingly as possible according to the planning. A subcapital resection of 0.5 to 1.0 cm with a straight oblique osteotomy of the femoral neck is usually sufficient. Prepare the prosthetic bed with the curved universal rasp that is used to prepare the path for the shaping rasps. Insert the rasp into the inner cortical bone using gentle rotations until you reach the lateral cortical bone near the lower edge of the lesser trochanter. Connect the rasp size 0 with the curved rasp handle for the minimally invasive preparation and advance it along the curvature of the calcar bone. This rasp can later be used to finely prepare the lateral cortical ring. After that, finish preparing the bone using the cancellous bone compactors. The aim of preparing the bone with compactors is to compress the cancellous bone. Prepare the bone step by step until you reach the planned size or until the compactor comes into contact with the cortical bone in the load-bearing areas. Mechanical stability and a change in the sound during impaction signals contact with the cortical bone. Guide the rasp and the compactors using small circular motions to prevent creating free spaces that cannot be filled by the prosthesis. The compactor should end outside or be flush with the resection line. Use the next size if necessary. Once the suitable compactor is in the correct position, remove the rasp handle and attach the trial head. Carry out a trial reduction. Clinically ensure that there is unobstructed freedom of movement and check the telescoping and the leg length. Trial heads in the diameters 22*, 28, 32 and 36 mm with the neck lengths XS*, S, M, L and XL (*optional) are available. After reduction, it is recommended to use a C-arm X-ray unit to

assess the fit, centre of rotation and offset of the implants. Remove the trial head and the compactor. If necessary, continue to prepare the prosthetic bed using the rasps and a larger compactor. On the lateral side, where a bony conflict can develop between the compactor and the cortical femoral neck ring, widen the bone using the awl until the compactor can be inserted as desired.

Manually insert the original NANOS prosthesis that corresponds to the size of the last compactor used. After that, drive in the prosthesis to the depth prepared by the compactor using the insertion instrument. Carefully clean the cone before attaching and driving in the original prosthetic head in the correct size. Check the reduction, stability and range of motion. The wound is closed in the usual way routinely used by the surgeon. During stem implantation it is advisable to leave the protective cap on the cone. This will protect the cone against damage. Should it become necessary during surgery to remove an implanted original NANOS prosthesis, an extraction instrument is available. The extraction instrument is attached to the cone and connected to the rasp handle for the removal.

The wound is closed in the usual way routinely used by the surgeon.
For further information see the NANOS surgical technique (see 1.3. Accessories).

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar data is indispensable and provides important information about suitable implants, placement, possible component combinations and enables the size of the implant to be used to be preselected. The operation may be performed only if the patient's material compatibility has been determined. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available in case another size is required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

- Primary and secondary osteoarthritis of the hip
- Dysplasia osteoarthritis
- Avascular and post-traumatic femoral head necrosis
- Patients with good bone density and intact femoral neck

6. Contraindications

- Marked osteoporosis
- Dysmorphia and loss of bone substance, as well as prior surgical procedures that do not guarantee the intended support
- Marked coxa valga
- Marked coxa vara with a femoral neck implant bed $< 125^\circ$
- Intolerance to the materials

If femoral heads with a neck XL are used, the range of motion is reduced by about 30° .

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Obesity (Body-Mass-Index BMI > 30)
- Anticipated extreme stresses (e. g. through sports or work), especially in patients with a body weight over 100 kg
- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Growth in children and adolescents
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

8. Possible negative effects

The negative effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedem/fluids

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Muscle spasm
- Stiffness
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Metallosis
- Elevated metal ions in blood
- Coxa vara
- Osteolysis
- Heterotopic ossification
- Pseudotumours

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, even though the contraindications described above partially apply to the patient, the patients should be informed about the possible effects of these circumstances and the anticipated risks. It should be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. It is also recommended that the patient is informed of activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked

for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 125).

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OHST Medizintechnik AG izmantotie simboli ir atrodami pielikumā (125. lpp.).