OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 (0)3385 - 54 20-0

Fax: +49 (0)3385 - 54 20-99

www.ohst.de

CE 0482

ELEC® und ELEC®plus Hüftkopf Deutsch (de)	3
ELEC® and ELEC®plus Femoral Head English (en)	10
Tête fémorale ELEC® et ELEC®plus Français (fr)	17
Testa femorale ELEC® e ELEC®plus Italiano (it)	24
Cabeza femoral ELEC® y ELEC®plus Español (es)	31
ELEC® en ELEC®plus dijbeenkop Nederlands (nl)	38
Cabeça femoral ELEC® e ELEC®plus Português (pt)	45
Κεφαλή ισχίου ELEC® και ELEC®plus Ελληνικά (el)	52

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IMPLANT

ELEC® and ELEC®plus Femoral Head

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance. These implants may be used only by physicians with appropriate experience and practice in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result. Care must always be taken to use the most up-to-date version of the instructions for use that is provided with the system components. The date of revision is shown on the cover sheet in the footer.

1. Product description and implant materials

The ELEC® and ELEC®plus Femoral Heads are designed for fixation onto endoprosthetic hip stems with a 12/14 cone authorised by the manufacturer. ELEC® (ISO 6474-1) and ELEC®plus (ISO 6474-2) are high-purity implant ceramics. The materials are bioinert, biocompatible, biostable, mechanically stable, corrosion-resistant and there are no ions to penetrate the body. The product properties also combine high shock resistance and tensile strength, excellent fracture toughness, an extremely high degree of hardness and an excellent resistance to fatigue. Compared to metal heads, there is less polyethylene wear, as is clinically proven when the ELEC® or ELEC®plus Femoral Head is used in combination with a PE component in the joint.

Caution: The ELEC® and ELEC®plus Femoral Heads are to be used exclusively with a PE cup, a PE inlay, an ELEC®plus inlay or the bipolar head.

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Reference number
ELEC® Femoral Heads		
ELEC® Femoral Head Ø 28-12/14 S	ISO 6474-1 Al ₂ O ₃ - ceramic	384-001
ELEC® Femoral Head Ø 28-12/14 M	ISO 6474-1 Al ₂ O ₃ - ceramic	384-002
ELEC® Femoral Head Ø 28-12/14 L	ISO 6474-1 Al ₂ O ₃ - ceramic	384-003
ELEC® Femoral Head Ø 32-12/14 S	ISO 6474-1 Al ₂ O ₃ - ceramic	384-004
ELEC® Femoral Head Ø 32-12/14 M	ISO 6474-1 Al ₂ O ₃ - ceramic	384-005
ELEC® Femoral Head Ø 32-12/14 L	ISO 6474-1 Al ₂ O ₃ - ceramic	384-006
ELEC® Femoral Head Ø 36-12/14 S	ISO 6474-1 Al ₂ O ₃ - ceramic	384-007
ELEC® Femoral Head Ø 36-12/14 M	ISO 6474-1 Al ₂ O ₃ - ceramic	384-008

Name	Material	Reference number
ELEC® Femoral Head Ø 36-12/14 L	ISO 6474-1 Al2O3 - ceramic	384-009

ELEC®plus Femoral Heads		
ELEC®plus Femoral Head Ø 28 S	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-001
ELEC®plus Femoral Head Ø 28 M	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-002
ELEC®plus Femoral Head Ø 28 L	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-003
ELEC®plus Femoral Head Ø 32 S	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-004
ELEC®plus Femoral Head Ø 32 M	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-005
ELEC®plus Femoral Head Ø 32 L	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-006
ELEC®plus Femoral Head Ø 32 XL	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-007
ELEC®plus Femoral Head Ø 36 S	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-008
ELEC®plus Femoral Head Ø 36 M	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-009
ELEC®plus Femoral Head Ø 36 L	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-010
ELEC®plus Femoral Head Ø 36 XL	ISO 6474-2 Al2O3 / ZrO2 - ceramic	013-011

1.2 Instrument overview

The instruments of the OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Reference number
Head End Ø28/32mm for Femoral Head Impactor Handle	506-062
Head End Ø36mm for Femoral Head Impactor Handle	506-1212
Femoral Head Impactor Handle with Silicone Handle grey, L=192mm	506-060

1.3 Accessories

Name	Reference number
X-ray template ELEC® & ELEC®plus Femoral Head Ø 28	013-024
X-ray template ELEC® & ELEC®plus Femoral Head Ø 32	013-025
X-ray template ELEC® & ELEC®plus Femoral Head Ø 36	013-026
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance is required with the appropriate hygiene regulations during removal

of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must not therefore come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. Therefore, damage after explantation must be assumed which excludes any reuse.

2.2 Authorised component combinations

We guarantee compatibility of our products only in combination with our own CE-marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

2.3 Information for use

The cone plug connection provides for reliable and torsion-proof clamping of the modular femoral head. It ensures a uniform transmission of power on the stem. The ELEC® and ELEC®plus femoral heads have to exhibit a perfect fit with the cone of the prosthetic stem. The size of the cone is marked on the label and also on the implant if possible.

Only implants that are brand-new may be used. Always compare the identification of the cone dimensions of the femoral stem with those of the ceramic femoral heads.

Caution: The ball diameter of the femoral head prosthesis must agree perfectly with the spherical nominal diameter of the part of the acetabular cup implant which articulates as the sliding partner.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Before inserting the implant, the implant bed must be irrigated. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

Caution: When using high-frequency surgical instruments (e.g. cauterly knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

Plug the femoral head onto the prosthetic stem. The goal is to have it fit perfectly on the stem taper. To ensure reliable function, utmost care is an essential prerequisite. Do not remove the plastic cap that protects the stem taper from damage until the trial femoral head is ready to be attached.

Note the following points before attaching the ELEC® or ELEC®plus Femoral Head:

1. Thoroughly rinse the stem taper with water to remove any tissue fragments, bone splinters or residual cement.
2. Dry the stem taper with a clean longuette and check the position and function of the joint after attaching the trial head.
3. Remove the trial head and rinse and dry it again.
4. Attach the femoral head using light axial pressure while simultaneously screwing it onto the stem taper until it sits tightly in place without moving.
5. Next, place the plastic head impactor onto the head pin and tightly clamp the femoral head onto the stem with a gentle hammer blow in the axial direction. The hammer blow causes a plastic deformation of the surface structure of the stem taper, thus ensuring an optimal distribution of pressure. This results in a torsion-proof fixation.

Caution: Never use a metal hammer to hit a ceramic femoral head! Only use the plastic head impactor intended for this purpose.

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch

is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar is indispensable and provides important information about suitable implants, placing and possible component combinations and enables the size of the implant to be used to be preselected. Surgery may only be performed once it has been established that the patient is able to tolerate the implant material. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1,15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available should another size be required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

- Total hip replacement, primary surgery in combination with prosthetic stems
- Revision (= removal of the endoprosthesis), repeated surgery with new prosthetic stems and cup
- Hemiarthroplasty (partial prosthesis)

6. Contraindications

- Revision involving leaving the stem and/or hip cup in situ (risk of component fracture due to damaged fixation cones)

In case of breakage of a ceramic femoral head, the pairing of metal (femoral head) to polyethylene (acetabular component) as well as metal to metal is contraindicated in revision surgery.

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures

- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Growth in children and adolescents
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

8. Possible negative effects

The negative effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedem/fluids

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Muscle spasm
- Stiffness
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Metallosis
- Elevated metal ions in blood
- Coxa vara
- Osteolysis
- Heterotopic ossification
- Pseudotumours

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, although the contraindications described above partially apply to the patient, it is particularly important to point out to the patient the effects of these circumstances on the success. It must be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. The patient has to be informed about activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adhesive labels are enclosed for documenting the implant used. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 120).

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