**OHST**  
medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 (0)3385 - 54 20-0

Fax: +49 (0)3385 - 54 20-99

www.ohst.de

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## Reusable Instruments and Trays

Before using the product, the user is under obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance. Detailed information on compatibility with other medical devices, product-specific risks, indications or contraindications is shown in the instructions for use or surgical technique (of the implant systems) specific to the respective system.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website [www.ohst.de/ifu/](http://www.ohst.de/ifu/).

These instructions for use apply to reusable instruments and trays (hereafter referred to as "instrument") of OHST Medizintechnik AG, which are used for the implantation of endoprostheses and for surgical interventions. These instruments may be used only by medical specialists with appropriate experience and practice in the relevant specialist field. In case of questions, please contact the manufacturer directly.

### 1. Product Description and Materials

#### 1.1. General Information

Instruments manufactured by OHST Medizintechnik AG are always part of a system. They may only be used with the original parts that belong to the respective system and only with the original implants that form part of these systems. Use of the instruments for other purposes is not permitted. Any manipulation of instruments is prohibited. If not used as intended, instruments can wear out faster, break or otherwise lose their function. If a product of OHST Medizintechnik AG is passed on to others, anyone doing so must ensure that the product can be traced at any time (LOT tracking) and that these instructions for use are known.

#### Custom-made devices

If an instrument is a custom-made device, this is noted on the packaging. These instruments do not carry a CE mark.

A custom-made device may only be used by the doctor if there is no commercially available instrument that can be used instead. The doctor must inform the patient in good time that the instrument is a custom-made device with the corresponding inherent risks. In addition, the doctor is responsible for obtaining the written consent of the patient in good time.

### Instruments connected to an active drive

In the case of instruments that are intended for use with an active driving component, it must be ensured before use that the connections of the instrument and the drive component match. The use of drills and milling cutters is to be carried out in conjunction with irrigation with Ringer's solution in order to prevent overheating and avoid damage to the tissue. The information provided by the equipment manufacturer must be adhered to in order to ensure proper use.

#### 1.2. Materials

The following materials are used:

- Metallic alloys and coatings
- Plastics such as PPSU, POM, PP, PEEK, silicone rubber

More information about the chemical and mechanical properties of the materials used is contained in DIN EN ISO 16061 and the additional material standards referred to therein, and is also obtainable from the manufacturer.

#### 1.3. General Risk Factors

During the use of instruments, allergic reactions to the material used are possible, as is loosening, wear, corrosion, ageing and fracturing of the instrument or instrument parts.

During the use of milling cutters, drills and other cutting instruments frictional heat may arise, which can lead to cell damage. Instruments may have sharp edges. When using synthetic gloves there is a danger that the gloves will be destroyed – note the risk of infection! Rough, sharp or cutting surfaces of the instruments must not come into contact with clothing or other materials that release fibres. Instruments with a large lever arm may transmit considerable forces and break if not handled properly. In this case it must be ensured that no fragments are left in the wound.

Impact plates may be deformed by blows with the hammer. There is a risk that parts may be split off. Any fragments that have fallen into the wound must be removed.

## 2. Preparation of Instruments

**Caution:** Only instruments that have been put through a validated preparation process may be used or reused!

Brand-new instruments and instruments returned after being repaired must be removed from the transport packaging before being put into storage and/or introduced into the instrument cycle. Protective caps and films must be removed and disposed of properly.

The instruments are supplied unsterile. Before initial use, brand-new instruments and instruments returned after being repaired must undergo full preparation in the same manner as the instruments already in use.

All instruments must be cleaned, disinfected and sterilised after every use. Effective cleaning and disinfection is an essential prerequisite for effective sterilisation.

The operator is responsible for the cleaning/disinfection and sterilisation process, exclusively by qualified personnel, and for the regular maintenance and care of the currently qualified cleaning and sterilisation equipment (e.g., according to EN ISO 15883, EN ISO 17665). This includes in particular ensuring compliance with validated parameters and processes.

As part of the responsibility for the sterility of the instruments during use, the operator must ensure in all cases that the prescribed procedure – appropriately validated for the specific device and product – for cleaning/disinfection and sterilisation is properly applied.

The user must also comply with the applicable laws in his country, as well as the hygiene regulations of the doctor's surgery or hospital. This applies in particular to the various requirements with respect to effective prion inactivation.

It should also be noted that some products require additional aspects, which are listed in separate product-specific instructions for use.

The recommendations listed under section 2 are for information purposes only. No liability is accepted with regard to the sterility of the instruments cleaned, disinfected or sterilised by the purchaser or user, nor with respect to re-sterilised instruments.

### 2.1. Cleaning and Disinfection

#### Basics

If possible, a validated automated procedure in the cleaning disinfection device (CDD) should be used for the cleaning and disinfection of the instruments. Owing to the significantly lower level of efficacy and reproducibility, an exclusively manual method – also involving use of an ultrasonic bath – should be followed only if an automated procedure is not available.

The use of ultrasound (applying the parameters specified by the manufacturer of the ultrasonic device) for pre-cleaning is permissible if the instruments have no joints / movable parts / silicone handles or have been dismantled. If ultrasound is used, individual instruments or instrument parts must not touch one another, in order to avoid secondary damage.

Pre-treatment of the instruments at the place of use and pre-cleaning can be carried out with regular tap water. For all other cleaning steps, deionised water (purified water) must be used.

For the manual removal of contamination on instruments and trays, only soft brushes or clean, soft and lint-free cloths may be used which are used for this purpose only, but never metal brushes or steel wool.

If manual cleaning steps in Table 1 or Table 3 for instruments with multiple and varying construction features contradict one another, as a general rule the construction feature of the more critical group and the corresponding cleaning workflow applies.

### Pre-treatment of the instruments at the place of use

Immediately after surgery, coarse visible impurities must be removed from the products by means of water or a cloth. The instruments and trays must be brought to the cleaning department as soon as the internal workflows allow. Drying up of impurities is to be prevented. Metal instruments should never be stored in physiological saline solution (NaCl solution) as prolonged contact leads to pitting and stress corrosion. It must be noted that storing the instruments in moist towels can encourage microbial growth.

### Pre-cleaning

To achieve an optimal final result after manual or automated cleaning, effective pre-cleaning is indispensable, as this helps to keep the microbial and protein load low.

The following steps must be carried out for this:

1. Demountable instruments must be taken apart for the cleaning (note product-specific directions for use, assembly/disassembly instructions, if available).
2. The instruments must be placed in water at room temperature (16 °C to 30 °C) for at least 5 minutes.

Please ensure the following:

- All surfaces are wetted with water (if necessary, use a syringe to reach poorly accessible places).
- Drill-holes/internal contours are filled with water.
- Seals must be opened.
- Move movable parts while soaking.
- Instruments must not touch one another.

It should be noted that any disinfectants used are for personal safety only and cannot replace the subsequent disinfection step to be performed upon completion of the cleaning.

Pre-cleaning is followed by the automated or manual cleaning/disinfection procedure.

### Automated cleaning/disinfection procedure

In the selection of the CDD it must be ensured that:

- the CDD generally has a qualified state (e.g., according to EN ISO 15883);
- a validated program for thermal disinfection ( $A_0$ -value  $\geq 3000$ ) is used (chemical disinfection harbours the risk of disinfectant residues on the instruments);
- the program used is suitable for the instruments and includes sufficient rinse cycles;
- only sterile water or water with a low bacterial count (max. 10 bacteria/ml) and a low endotoxin count (max. 0.25 endotoxin units/ml) is used (e.g., purified water);
- the air used for drying is filtered accordingly;
- the CDD is maintained on a regular basis and kept in the qualified state.

In the selection of the process chemicals used, it must be ensured that:

- these are fundamentally suitable for the cleaning of the instruments;
- if no thermal disinfection is used, a suitable disinfectant with proven efficacy (e.g., DGHM or FDA approval or CE mark) is also used and that this is compatible with the process chemicals used;
- the process chemicals used are compatible with the instruments (see section 2.3);
- all surfaces of the instruments are accessible to the process chemicals.

Prior to automated cleaning, the instruments must be pre-treated according to the applicable construction features (see Table 1).

1. Clean the instruments in a water bath at room temperature (16 °C to 30 °C) with a soft brush or with a clean, soft and lint-free cloth until all visible impurities have been removed. All places that can be reached with a brush must be brushed; the use of a bottle brush (see Table 1) may be expedient.
2. Rinse the instruments with a water pistol (static pressure at least 3 bar) with water at room temperature (16 °C to 30 °C) for at least 15 seconds per instrument. The use of a bottle brush (see Table 1) may be expedient. Pay particular attention to poorly accessible areas (e.g., drill-holes, covered surfaces).
3. According to Table 1 the instruments are to be placed in an alkaline cleaning solution heated to between 37 °C and 40 °C (pH range: 10.4 – 10.8) with or without ultrasound for at least 5 minutes. The instruments must be fully immersed. The use, dosing and exposure time of the cleaning solution must be in accordance with the manufacturer's instructions.
4. After that, rinse the instruments with deionised (purified) water (see Table 1).
5. If impurities are still visible after the cleaning procedure, repeat the entire cleaning procedure.
6. Only when no more impurities are visible with the naked eye are the instruments sent for automated cleaning and disinfection (Table 2).

Table 1: Pre-treatment for automated cleaning and disinfection

| Design features<br>Construction features   | Brushing                 | Rinsing<br>(Water pistol) | Place in<br>cleaner (with<br>ultrasound) | Place in<br>cleaner<br>(without<br>ultrasound) | Further rinsing |
|--|--------------------------|---------------------------|--|--|-----------------|
| <b>Group 1:</b><br>- Joints<br>- Movable parts that cannot be disassembled<br>- Instruments with silicone handle   | Step 1                   | Step 2                    | Not permissible                          | Step 3   | Step 4          |
| <b>Group 2:</b><br>- Through-holes with inner workings (parts that protrude into the drill-hole)<br>- Blind holes with and without a thread  | Not permissible          | Step 1                    | Step 2                                   | Not permissible                                | Step 3          |
| <b>Group 3:</b><br>- Through-holes without inner workings<br>- Immovable parts that cannot be disassembled<br>- Through-holes with a thread<br>- Permeable opening slumens<br>- Structured, rough surfaces | Step 1<br>(Bottle brush) | Step 2                    | Step 3                                   | Not permissible                                | Step 4          |
| <b>Group 4:</b><br>- Simple smooth surface (without internal contours)<br>- All areas visible  | Step 1                   | Step 2                    | Step 3                                   | Not permissible                                | Step 4          |

Legend:

Step = Sequence of the cleaning procedure which must be followed

Not permissible = Cleaning step which must not be used for this construction feature

Note:

The criticality decreases from Group 1 to Group 4.

The concentrations stated by the manufacturer of the process chemicals must be strictly adhered to. When immersing the instruments it must be ensured that they do not touch one another. If applicable, the instruments must be connected to flexible rinsing tubes.

Position drill-holes and lumens in the CDD such that flushing is possible. An example of a validated workflow for automated cleaning and disinfection is given in Table 2. A qualified CDD was used for this.

Table 2: Example of a validated workflow for automated cleaning and disinfection

| Program block        | Parameter  |
|----------------------|--|
| Pre-cleaning         | <ul style="list-style-type: none"> <li>• Cold water influx</li> <li>• At least 2 min. exposure</li> <li>• Emptying</li> </ul>  |
| Cleaning             | <ul style="list-style-type: none"> <li>• Cold-warm water influx</li> <li>• Dosing of alkaline cleaning solution (pH range: 11.6 – 11.8) at min. 40 °C according to manufacturer's instructions</li> <li>• Heating to 55 °C and at least 5 min. exposure at min. 55 °C</li> <li>• Emptying</li> </ul> |
| Neutralisation       | <ul style="list-style-type: none"> <li>• Cold water influx</li> <li>• At least 3 min. exposure</li> <li>• Emptying</li> </ul>  |
| Rinsing              | <ul style="list-style-type: none"> <li>• Cold water influx</li> <li>• At least 2 min. exposure</li> <li>• Emptying</li> </ul>  |
| Thermal disinfection | <ul style="list-style-type: none"> <li>• Purified water influx</li> <li>• Heating to 90 °C and at least 5 min. exposure at min. 90 °C</li> <li>• Emptying</li> </ul>   |
| Drying               | <ul style="list-style-type: none"> <li>• Drying for 30 min. at min. 110 °C</li> </ul>  |

The drying time can vary depending on how the CDD is loaded and is to be determined individually by the operator performing the preparation.

After the end of the program, the instruments are disconnected from the rinsing hoses and removed from the CDD under low-bacteria conditions.

Proof of the general suitability of the instruments for effective automated cleaning and disinfection has been provided by independent, accredited test laboratories using a CDD. In this context, the method described above was applied.

Upon conclusion of the automated cleaning / disinfection procedure, the instruments are inspected in accordance with the "Checks" section.

### Manual Cleaning/Disinfection Procedure

During manual cleaning, the following steps must be carried out in consideration of the construction features (see Table 1):

1. Clean the instruments in a water bath at room temperature (16 °C to 30 °C) with a soft brush or with a clean, soft and lint-free cloth until all visible impurities have been removed. All parts that can be reached with a brush must also be scrubbed. The use of a bottle brush may be expedient (see Table 3).
2. Rinse the instruments with a water pistol (static pressure at least 3 bar) with water at room temperature (16 °C to 30 °C) for at least 15 seconds per instrument. The use of a bottle brush (see Table 3) may be expedient. Pay particular attention to poorly accessible areas (e.g., drill-holes, covered surfaces).
3. According to Table 3 the instruments are to be placed in an alkaline cleaning solution heated to between 37 °C and 40 °C (pH range: 10.4 – 10.8) with or without ultrasound for at least 5 minutes. The instruments must be fully immersed. The use, dosing and exposure time of the cleaning solution must be in accordance with the manufacturer's instructions.
4. After that, rinse the instruments with deionised (purified) water (see Table 3).
5. If impurities are still visible after the cleaning procedure, repeat the entire cleaning procedure.
6. Only when no more impurities are visible with the naked eye are the instruments placed fully in an immersion bath with a disinfectant solution (pH value approx. 8.5) (see Table 3). The use, dosing and exposure time of the disinfectant solution must be in accordance with the manufacturer's instructions.

Please ensure the following:

- All surfaces are wetted with the disinfectant solution.
  - Drill-holes/internal contours are filled with water.
  - Seals must be opened.
  - Instruments must not touch one another.
7. After the exposure time, the instruments are rinsed thoroughly with deionised (purified) water and then dried (see Table 3).

Upon conclusion of the manual cleaning, the instruments are inspected according to the "Checks" section.

Table 3: Manual cleaning and disinfection instructions for instruments

| Design features<br>Construction features  | Brushing                    | Rinsing<br>(water<br>pistol) | Place in<br>cleaner<br>(with<br>ultrasound) | Place in cleaner (without<br>ultrasound) | Further<br>rinsing | Disinfecting | Further<br>rinsing |
|---|-----------------------------|------------------------------|---|--|--------------------|--------------|--------------------|
| <b>Group 1:</b> <ul style="list-style-type: none"> <li>- Joints</li> <li>- Movable parts that cannot be disassembled</li> <li>- Instruments with silicone handle</li> </ul>   | Step 1                      | Step 2                       | Not permissible                             | Step 3                                   | Step 4             | Step 5       | Step 6             |
| <b>Group 2:</b> <ul style="list-style-type: none"> <li>- Through-holes with inner workings (parts that protrude into the drill-hole)</li> <li>- Blind holes with and without a thread</li> </ul>  | Not permissible             | Step 1                       | Step 2                                      | Not permissible                          | Step 3             | Step 4       | Step 5             |
| <b>Group 3:</b> <ul style="list-style-type: none"> <li>- Through-holes without inner workings</li> <li>- Immovable parts that cannot be disassembled</li> <li>- Through-holes with a thread</li> <li>- Permeable openings/lumens</li> <li>- Structured, rough surfaces</li> </ul> | Step 1<br>(Bottle<br>brush) | Step 2                       | Step 3                                      | Not permissible                          | Step 4             | Step 5       | Step 6             |
| <b>Group 4:</b> <ul style="list-style-type: none"> <li>- Simple smooth surface (without internal contours)</li> <li>- All areas visible</li> </ul>  | Step 1                      | Step 2                       | Step 3                                      | Not permissible                          | Step 4             | Step 5       | Step 6             |

Legend:

Step = Sequence of the cleaning procedure which must be followed

Not permissible = Cleaning step which must not be used for this construction feature

Note:

The criticality decreases from Group 1 to Group 4.

## Checks

**Caution:** The instruments are subject to natural wear through use and preparation. Therefore, instruments must be checked for damage and function before every use!

After cleaning, all instruments must be visually inspected with respect to corrosion, damaged surfaces, chips, impurities and function. Damaged instruments are to be rejected and replaced. Instruments which are still contaminated must be reprocessed.

## Maintenance

Disassembled instruments are to be reassembled prior to sterilisation and their function checked (see product-specific instructions for use).

If the use of an oil for maintenance or care of the instruments is required, it must be ensured that only instrument oils (white oil) are used which – taking into account the maximum applied sterilisation temperature – are approved for steam sterilisation and have proven biocompatibility. Excess oil must be wiped off. After oiling, the instrument must be sterilised.

## Packaging

After cleaning and before sterilisation, the instruments must be placed in the tray and packaged together in disposable sterilisation packaging (double packaging) and sterilisation containers that meet the following requirements:

- According to EN 868/EN ISO 11607
- Suitable for steam sterilisation (temperature stability up to at least 137 °C, sufficient steam permeability)
- Adequate protection of the instruments or sterilisation packaging against mechanical damage
- Regularly serviced according to the manufacturer's instructions (sterilisation containers)

## 2.2. Sterilisation

**Caution:** The instruments must not be sterilised in the protective packaging provided. New instruments are also to be cleaned prior to sterilisation! Storage of the instruments for sterilisation takes place in the intended mounting brackets of the corresponding tray. If no respective tray is available, make sure that the instruments do not touch one another and position the instruments in such a way that residual moisture after drying is avoided.

Only steam sterilisation is permissible for the sterilisation. The points listed below must be observed for this. Other sterilisation methods are not permitted.

## Steam sterilisation

- Fractionated vacuum procedure at least three times (with sufficient product drying)
- Steam steriliser in accordance with EN 13060 or EN 285
- Validated according to EN 554/EN ISO 17665 (valid commissioning and product-specific performance assessment)
- Sterilisation temperature of 134 °C (273 °F; plus tolerance according to EN 554/EN ISO 17665). Apply a hold time of at least 3 min. Please note the nationally applicable instructions.
- Minimum drying time: 20 min
- The drying time can vary depending on the load and is to be determined individually by the operator performing the preparation.

All surfaces must be accessible to the water vapour. The instruments may be sterilised only loose or in their original retainers and not one on top of the other in a sterilisation container that satisfies the above requirements (see section 2.1).

Proof of general suitability of the instruments for effective sterilisation has been provided by independent, accredited test laboratories. Here, the sterilisation procedures described above were included and a sufficient sterilisation safety level (SAL  $>10^{-6}$ ) was achieved.

### 2.3. Material Compatibility

Instruments must not come into contact with agents containing chlorine or fluorine. The process chemicals used must not contain the following ingredients:

- Mineral acid, with the exception of phosphoric acid
- Oxidizing acids
- Strong lyes (pH  $>12.5$ )
- Aromatic hydrocarbons, benzines
- Strong oxidizing agents
- Trichlorethylene/perchloroethylene

If in doubt, please contact the manufacturer of the process chemicals.

All instruments and trays must not be exposed to temperatures higher than 137 °C (279 °F)!

### 2.4. Reusability

The instruments can, with appropriate care and provided that they are free from damage and contamination, be reused. However, they must be inspected for proper function and damage before each use!

Any further use or the use of damaged or contaminated instruments is the responsibility of the user. If this is ignored, all liability is excluded.

## 3. Storage and Handling

Instruments may only be stored in dust-proof dry areas/cabinets at room temperature. Under no circumstances may instruments be stored in the immediate vicinity of chemicals which because of their ingredients can give off corrosive vapours (e.g., active chlorine).

Instruments are sensitive to damage. They must therefore be treated carefully. Dents, scratches or other mechanical damage to the surfaces cause excessive wear and can cause corrosion and improper use. Instruments are to be transported in their respective trays, if available, or in other appropriate containers.

Instruments must be examined before use to ensure they are functioning properly. Any instruments that are no longer functioning correctly must be withdrawn from circulation.

Flexible shafts and drills may only be applied with a maximum torque of 0.2 Nm at a maximum inflexion of 45°. When using the instrument, the user must ensure that the bending radius of the flexible part is as uniform as possible. If the instrument is used with an uneven or S-shaped curve, the useful life may be shortened. The dimensional stability of the flexible coil is guaranteed up to a torque of 0.5 Nm.

When using flexible drills, the drilling gauge belonging to the system must be used.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

#### 4. Repair / Disposal

Damaged or improperly treated instruments or instruments altered without authorisation must no longer be used.

The packaging components and instruments are to be sent for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer these instruments may also be returned to the manufacturer for proper disposal at no cost. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

#### 5. Explanation of the Label Symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 187).

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