

## Knee Joint Replacement

### Patient Information



## Dear patient,

You have received a knee joint replacement from OHST Medizintechnik AG. In the following we would like to give you important information in relation to the product. Please read the information carefully and follow the instructions. This document is intended to fulfil the requirements arising from Article 18 and Section 23 of Annex I to Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

## Intended purpose

Knee prostheses are implants that remain in the body permanently and either partially or fully replace the damaged knee. The knee joint in its entirety is replaced by three components: the femoral component, the tibial component and the patellar component.

Irrespective of the shape and fixation manner of the individual joint components, they always perform the same functions (ISO 7207-1, ISO 21536):

**Femoral component:** the component of a total knee joint replacement intended to be secured to the femur to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.

**Tibial component:** the component of a total knee joint replacement intended to be secured to the tibia to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.

**Patella component:** the component of a total or partial knee joint replacement which is used to replace the articulating surface of the patella. These implants can be manufactured as one component or a set of components to be assembled by the user.

The aim is to relieve pain and improve the function of the knee joint in patients who have undergone surgery.

## Safety and performance

Register data show<sup>1</sup> that about 82 % of all total knee arthroplasties can last 25 years. Depending on the type of fixation in the bone, the following failure probabilities are possible in the course of the wearing time of the prosthesis<sup>2</sup> :

Type of fixation in the bone	Probability of failure within...							
	1 year	2 years	3 years	4 years	5 years	6 years	7 years	8 years
cemented	1,7 %	2,5 %	3,0 %	3,3 %	3,6 %	3,8 %	4,0 %	4,3 %
hybrid	1,8 %	2,7 %	3,1 %	3,4 %	3,6 %	3,9 %	4,1 %	4,2 %
cementless	1,9 %	3,0 %	3,8 %	4,0 %	4,1 %	4,3 %	4,3 %	4,3 %

The lifespan of the artificial knee joint replacement can be significantly influenced by factors such as age, gender, body mass index (BMI) and concomitant diseases. The earlier a prosthesis is inserted, the more likely it is that the implant will be changed. Male patients have a higher probability of failure than female patients. A high BMI as well as the number of concomitant diseases can reduce the lifespan of the prosthesis. Hospitals and doctors with greater experience due to higher treatment numbers tend to reduce the risk of prosthesis replacement.

Reasons for a subsequent intervention on the knee joint can be: Infection, loosening, osteolysis, fracture, instability, implant wear, implant malalignment or movement restrictions.

### Summary of safety and clinical performance

The EU is in the process of setting up a European Database on Medical Devices (EUDAMED). As soon as this is available, you will be able to view a so-called "summary of safety and clinical performance" for your specific implant, which will be updated regularly. Until EUDAMED is launched, this report can be requested from OHST Medizintechnik AG.

<sup>1</sup> Evans JT et al: How long does a knee replacement last? A systematic review and meta-analysis of case series and national registry reports with more than 15 years of follow-up. Lancet. 2019 Feb 16;393(10172):655-663

<sup>2</sup> Endoprosthesis Register Germany (EPRD) Annual Report 2023. DOI: 10.36186/reportepd092024

## Important information related to your implant

Clinical experience shows that the presence of one or more of the following concomitant circumstances (risk factors) can lead to shorter service life, more frequent complications or an overall poorer outcome of knee arthroplasty. This list is not exhaustive.

General risk factors and conditions:

- Overweight
- Alcoholism or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Intake of high-dose cortisone or cytostatics
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to knee arthroplasty:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle dysfunction of the affected joint
- Growth in children and adolescents
- Expected extreme stresses, e.g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

Undesirable effects

The negative effects listed below are among the most typical and common consequences of surgery:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedema

The negative effects listed below are among the most typical and common consequences of knee arthroplasty:

- Tibial or femoral fracture
- Sub-dislocation or dislocation of the patella
- Instability, migration or loosening of components
- Detachment of components
- Insufficient freedom of movement
- Metal hypersensitivity reactions
- Patellofemoral complications
- Neurovascular complications
- Paralysis/palsy of the peroneal nerve
- Fat embolism
- Arterial insufficiency and damage
- Arthrofibrosis / accretions
- Implant fractures
- Implant noises
- Reduced quality of life (pain, sleep disorders and a restricted range of motion, especially when lying down)
- Inflammations
- Metallosis
- Increase in the level of metal ions in the blood
- Pseudotumours

The occurrence of specific adverse effects may necessitate revision surgery.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the Competent Authority of the member state<sup>3</sup> in which the user and/or patient is established.

#### Activities of daily living

Depending on the profession, the prosthesis used and the course of recovery, it usually takes 2 to 3 months before it is possible to work again. Sometimes a gradual reintegration is possible, starting with a few hours per day and then gradually increasing the amount<sup>4</sup>.

Strong muscles, tendons and ligaments provide support for the prosthesis, which is why movement is also important with a joint replacement. However, the lifespan of a prosthesis also depends on how much it is used. Therefore, there are activities and sports that are more suitable for people with an artificial knee joint than others.

Activities that put less stress on the joint include:

- Walks or easy (not too mountainous) hikes
- Swimming
- Dance
- Light gymnastics

Activities and sports that put a lot of stress on the artificial joint are:

- Jogging
- Sports such as tennis, squash or football where the joints are exposed to impacts, twists or jerky movements
- Sports with physical contact and corresponding risk of falling, for example martial arts

For many sports, it depends on how well they were practised before the joint replacement<sup>5</sup>. Anyone starting a new sport is not yet familiar with it and therefore has an increased risk of unfavourable movements and accidents. Jobs that require a lot of kneeling, squatting or heavy lifting are not well suited with an artificial knee joint. The same applies to activities that require frequent standing on ladders, scaffolding, roofs or uneven surfaces.

It is best to consult your doctor about which activities are suitable.

#### Follow-up examinations

The implantation of the prosthesis is followed by clinical and radiological checks of the patient. The aim of these follow-up examinations is to detect any emerging complications at an early stage and thus to be able to treat them before more extensive damage to the joint or bone has occurred<sup>6</sup>.

There are no fixed intervals for a follow-up examination. Your treating doctor will therefore have established their own routine. Clinical examinations may be more frequent than radiological examinations in the first few months. The reason for this is the radiation exposure associated with each X-ray examination.

With regard to radiological follow-up examinations, we refer to the following recommendations of the German Society for Endoprosthetics (AE)<sup>7</sup>:

- Since aseptic loosening and wear of cemented total knee arthroplasties only occur in late failures, these prostheses only need to be checked at 2-3 year intervals from the 5<sup>th</sup> postoperative year onwards in patients with no or few complaints after an unremarkable postoperative check-up.
- Cementless total knee arthroplasties should also be checked one year after surgery. If no abnormalities are identified during this time, the check-up intervals can be the same as for cemented prostheses in patients with no or few complaints.
- The recommended check-ups should be communicated to the patients after the operation so that they can make use of them in a self-determined manner.

<sup>3</sup> [https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts\\_en](https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en) (accessed on 22.10.2024)

<sup>4</sup> Aresti N et al: Primary care management of postoperative shoulder, hip, and knee arthroplasty. *BMJ*. 2017 Oct 18;359:j4431

<sup>5</sup> Witjes S et al: Return to Sports and Physical Activity After Total and Unicondylar Knee Arthroplasty: A Systematic Review and Meta-Analysis. *Sports Med*. 2016 Feb;46(2):269-92

<sup>6</sup> Hass H: Stellungnahme zur Durchführung von Ganzbeinstandaufnahmen postoperativ April 2016.

[https://www.clarcert.com/\\_Resources/Persistent/0/e/1/5/0e15a2aaf9449835f76c642566839a75c8aaf7ae/\\_stellungnahme\\_ganzbeinstandaufnahme%20%28160401%29.pdf](https://www.clarcert.com/_Resources/Persistent/0/e/1/5/0e15a2aaf9449835f76c642566839a75c8aaf7ae/_stellungnahme_ganzbeinstandaufnahme%20%28160401%29.pdf) (accessed on 22.10.2024)

<sup>7</sup> Von Roth et al: Die radiologische Verlaufskontrolle von primären Hüft- und Knieendoprothesen – Empfehlung der Deutschen Gesellschaft für Endoprothetik (AE). *Z Orthop Unfall* 2020; 158(03): 276-279

Patients with leg axis deviations after implantation of the prosthesis may require shorter check intervals than patients without axis deviations. The reason for this is that asymmetrical loading of the artificial knee joint can lead to increased wear and thus to the premature occurrence of osteolysis or loosening of the prosthesis.
























In your own interest, you should therefore attend all agreed follow-up appointments, even if you have no complaints! If you have any complaints, contact your doctor immediately.

If you are to undergo magnetic resonance imaging (MRI), inform the examining doctor about all your implants. MRI examinations can have undesirable effects that can harm the patient. Possible effects include artefacts, heating of the implant, induction of electrical currents and loosening of the implant. These effects can vary greatly depending on the characteristics of the implant (material, size, position in the body, etc.) and the combination with other implants. Furthermore, there are different manufacturers and generations of MRI systems. OHST is therefore unable to make any statements about the safety of OHST implants with a specific MRI system. A patient-specific risk assessment by the examining physician prior to performing the MRI is therefore essential in order to decide whether and under what conditions the examination can be performed safely.

## Information about your prosthesis / endoprosthesis passport

When you are discharged from the clinic after the operation, you will receive an implant card (endoprosthesis passport). This contains all the important information about your implant. You should always carry this document with you.

The meaning of the symbols used, e.g. in the instructions for use, on the labels of OHST Medizintechnik AG or on the implant card, is listed below.

 ♦ Read instructions for use	 ♦ Batch code
 ♦ Caution	 ♦ Manufacturer
 ♦ Do not reuse	 ♦ Do not use if package is damaged
 ♦ Manufacturing date (year-month)	 ♦ Keep dry
 ♦ Use-by date (year-month)	 ♦ Keep away from sunlight
 ♦ Do not resterilise	 ♦ Unique Device Identifier
 ♦ Sterilised using irradiation	 ♦ Device name
 ♦ Sterilised using ethylene oxide	 ♦ Patient name or patient ID
 ♦ Double sterile barrier system with protective packaging outside	 ♦ Date of implantation
 ♦ Double sterile barrier system	 ♦ Name and address of implanting health care facility
 ♦ Reference number	 ♦ Website
 ♦ Serial number	

Knee joint prostheses from OHST Medizintechnik AG are made of the following materials. You can find detailed information on this on your implantation card.

<b>Stainless steel according to ISO 5832-1 (chemical composition)</b>	
<b>Element</b>	<b>Mass fraction in %</b>
Carbon	max. 0.030
Silicon	max. 1.0
Manganese	max. 2.0
Phosphorus	max. 0.025
Sulphur	max. 0.010
Nitrogen	max. 0.10
Chromium	17.0 to 19.0
Molybdenum	2.25 to 3.00
Nickel	13.0 to 15.0
Copper	max. 0.50
Iron	Balance

<b>Titanium 6-aluminium 4-vanadium wrought alloy according to ISO 5832-3 (chemical composition)</b>	
<b>Element</b>	<b>Mass fraction in %</b>
Aluminium	5.5 to 6.75
Vanadium	3.5 to 4.5
Iron	max. 0.3
Oxygen	max. 0.2
Carbon	max. 0.08
Nitrogen	max. 0.05
Hydrogen	max. 0.015
Titanium	Balance

<b>Cobalt-chromium-molybdenum casting alloy according to ISO 5832-4 (chemical composition)</b>	
<b>Element</b>	<b>Mass fraction in %</b>
Chromium	26.5 to 30.0
Molybdenum	4.5 to 7.0
Nickel	max. 1.0
Iron	max. 1.0
Carbon	max. 0.35
Manganese	max. 1.0
Silicon	max. 1.0
Cobalt	Balance

<b>Cobalt-chromium-molybdenum forging alloy according to ISO 5832-12 (chemical composition)</b>	
<b>Element</b>	<b>Mass fraction in %</b>
Chromium	26.0 to 30.0
Molybdenum	5.0 to 7.0
Iron	0.75 Maximum value
Manganese	1.0 Maximum value
Silicon	1.0 Maximum value
Carbon	0.35 Maximum value
Nickel	1.0 Maximum value
Nitrogen	0.25 Maximum value
Cobalt	Balance

<b>Ultra-high molecular weight polyethylene (UHMWPE) according to ISO5834-1/2 (ash and trace elements)</b>	
<b>Element</b>	<b>Maximum permitted quantity in mg/kg</b>
Ash	125
Titanium	40
Calcium	5
Chlorine	30
Aluminium	20

<b>Highly cross-linked ultra-high molecular weight polyethylene with the addition of vitamin E (alpha-tocopherol) (ash and trace elements)</b>	
<b>Element</b>	<b>Quantity</b>
Ash	max. 125 mg/kg
Titanium	max. 40 mg/kg
Calcium	max. 5 mg/kg
Chlorine	max. 30 mg/kg
Aluminium	max. 20 mg/kg
Vitamin E	1000 ±150 ppm

The following coatings can also be applied:

<b>Titanium Plasma Spray (TPS) according to ASTM F1580 (chemical composition)</b>	
<b>Element</b>	<b>Mass fraction in %</b>
Aluminium	max. 6.75
Vanadium	max. 4.5
Oxygen	max. 0.4
Iron	max. 0.5
Carbon	max. 0.08
Hydrogen	max. 0.05
Nitrogen	max. 0.05
Copper	max. 0.1
Tin	max. 0.1
Silicon	max. 0.04
Chlorine	max. 0.2
Yttrium	max. 0.005
Titanium	Balance

<b>BONIT®</b>	
Phase composition	≥ 70 % Brushite [CaH(PO <sub>4</sub> ) x 2 H <sub>2</sub> O] ≤ 30 % Hydroxyapatite [Ca <sub>5</sub> (PO <sub>4</sub> ) <sub>3</sub> OH]
Ca/P ratio	1,1 ± 0,1
Trace elements	Arsenic: max. 3 ppm Cadmium: max. 5 ppm Mercury: max. 30 ppm Lead: max. 50 ppm

<b>Titanium niobium nitride (TiNbN)</b>	
<b>Element</b>	<b>Mass fraction in %</b>
Titanium	max. 70
Niobium	max. 30
Iron	< 0,05
Oxygen	< 0,2
Carbon	< 0,2
Nitrogen	<0,05
Hydrogen	< 0,02

## Further information

AE - German Society for Endoprosthetics e.V.

The AE - German Society for Endoprosthetics was founded as a non-profit association. Its members are leading orthopaedic surgeons and trauma surgeons as well as scientists who deal with questions of endoprosthetics and alternative joint-preserving treatment methods.

It is a section of the German Society for Orthopaedics and Trauma Surgery (DGOU) and is thus responsible for all questions concerning arthroplasty.

Here you will find answers to frequently asked questions about diseases and treatments, including the provision of an artificial joint (endoprosthesis).

Internet address (German language only): <https://www.ae-germany.com>

Federal Institute for Drugs and Medical Devices (BfArM)

The Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM) is an independent federal higher authority within the portfolio of the Federal Ministry of Health.

The BfArM's tasks with regard to medical devices are mainly the central collection, analysis, and evaluation of risks arising from use or application of the devices and the co-ordination of the necessary measures to be taken. This work is based on incoming incident reports on medical devices.

Internet address: [https://www.bfarm.de/EN/Home/\\_node](https://www.bfarm.de/EN/Home/_node)

German Society for Orthopaedics and Orthopaedic Surgery

The Society sees itself as responsible for scientific issues vis-à-vis the public and the medical profession and thus as a contact for other societies and associations, institutions and authorities.

The purpose of the association is the promotion of orthopaedic science in research, teaching and practical application including the rehabilitation of physically disabled persons.

Here you will find information about orthopaedic issues.

Internet address (German language only): <https://dgooc.de>

EPRD German Endoprosthesis Register gGmbH

The German Arthroplasty Registry (EPRD) has been established to ensure and improve the quality of endoprosthetic surgery and thereby increase patient safety. The EPRD's extensive data pool will make it easier than ever to identify causes of failure in endoprosthetic interventions in future. By this, it will be easier to sort out whether the implants used, the surgical procedure and/or patient-specific features are responsible for a revision.

Since 2020, the EPRD has published its own patient information annually and in addition to its annual report. The publication contains key results of the data evaluation from the respective annual report - summarised briefly and in a way that is easy for patients to understand.

Internet address: <https://www.eprd.de/en>

EUDAMED - European Database on Medical Devices

EUDAMED will provide a living picture of the lifecycle of medical devices that are made available in the European Union (EU). It will integrate different electronic systems to collate and process information about medical devices and related companies (e.g. manufacturers). In doing so, EUDAMED aims to enhance overall transparency, including through better access to information for the public and healthcare professionals, and to enhance coordination between the different Member States in the EU.

Internet address: <https://ec.europa.eu/tools/eudamed>

### Institute for Quality and Efficiency in Health Care (IQWiG)

The Institute for Quality and Efficiency in Health Care (IQWiG) is an independent, scientific institution of the private and non-profit Foundation for Quality and Efficiency in Health Care. The foundation aims to support evidence-based decision making.

InformedHealth.org is the English-language version of the German website Gesundheitsinformation.de. By publishing this bilingual website, the Institute for Quality and Efficiency in Health Care (IQWiG, Germany) fulfills part of its legal mandate to inform the public in matters of health. The website addresses both patients and (healthy) consumers by offering a wide range of different topics.

IQWiG health information is written with the aim of helping people understand the advantages and disadvantages of the main treatment options and health care services.

Internet address: <https://www.informedhealth.org>

### Association of Statutory Health Insurance Funds Germany (vdek)

The data basis of the vdek Kliniklotsen are the quality reports provided by the individual hospitals, which were prepared in accordance with the valid regulations of the Federal Joint Committee. Thus, the Kliniklotse is able to provide you with comprehensive information on, for example, treatment, equipment and quality in hospitals.

Via the online portal [www.vdek-arztlotse.de](http://www.vdek-arztlotse.de), users can search specifically for established doctors, dentists, psychological psychotherapists and emergency outpatient clinics in Germany. In addition, the vdek doctor's guide provides comprehensive information on the respective practice, such as accessibility and office hours or the degree of accessibility, but also on the doctor's therapeutic focus.

Internet address (German language only): <https://www.vdek-kliniklotse.de>  
<https://www.vdek-arztlotse.de>

### Weisse Liste

The core component of the project is the portal [www.weisse-liste.de](http://www.weisse-liste.de) which has been online since 2008. The title stands for an orderly overview (Liste = "List") and indicates a strong affinity with the fields of medicine and health (Weisse = "White"). The portal provides highly practical help for its users in their search for the right provider in their area or as the case may be, nationwide. In this, it differs from other mostly commercial services: first of all, it carries absolutely no advertising and is free of charge. But more than this: there is a strong emphasis throughout on the quality of the information in all areas, underscored for example by the use of scientifically developed questionnaires or measures designed to prevent the manipulation of information.

The Weisse Liste aims to provide citizens with direction and enable them to make informed choices. Its purpose is to prepare and provide access to information on the services and quality offered by health care providers in a way that non-specialists can understand. It aims to create transparency and thereby to ensure that providers compete fairly to provide the best quality. In turn, this competition benefits the populace through improved quality across the board.

Internet address: <https://www.bertelsmann-stiftung.de/en/our-projects/weisse-liste>

**Manufacturer:**  
**OHST Medizintechnik AG**  
Grünauer Fenn 3  
14712 Rathenow  
Germany  
Phone+49 (0) 3385 5420 0  
Fax+49 (0) 3385 5420 99  
E-mail [info@ohst.de](mailto:info@ohst.de)  
Internet [www.ohst.de](http://www.ohst.de)

Disclaimer

This information is intended for medical laypersons.  
The explanations of the products contained in the information are of a general nature and do not constitute medical advice.

The information has been compiled by medical experts and professionally qualified employees of OHST AG to the best of their knowledge.

No liability or guarantee is assumed for the topicality, correctness and completeness of the information provided.  
Any liability for material or immaterial damage arising from the use of this information is excluded.

**CE** 0482