OHST

medical technology

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IMPLANT

NANOS femoral neck prosthesis

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for hip arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the websites <https://www.ohst.de/ifu-instructions/> and <https://ifu.smith-nephew.com/>.

1. Product Description and Implant Materials

The NANOS femoral neck prosthesis is an endoprosthesis (femoral component) for cementless fixation in total hip arthroplasty (THA). The NANOS femoral neck prosthesis is manufactured from a Ti6Al4V forge alloy (ISO 5832-3) and has a proximal osteoconductive titanium plasma coating (TPS) and an additional calcium phosphate coating (BONIT®). The distal tip is polished. The implant requires sparing and minimal bone resection. In all, 11 perfectly matched sizes are available.

The NANOS femoral neck prosthesis may only be used in patients who require total hip arthroplasty and who fulfil the relevant indications and contraindications. In addition, the NANOS femoral neck prosthesis may only be used to treat skeletally mature patients.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

1.1 Overview of Implants

Name	Material	Art. no.	SAP no.
NANOS femoral neck prosthesis size 0	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425000	75008154
NANOS femoral neck prosthesis size 1	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425001	75008155
NANOS femoral neck prosthesis size 2	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425002	75008156
NANOS femoral neck prosthesis size 3	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425003	75008157
NANOS femoral neck prosthesis size 4	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425004	75008158

Name	Material	Art. no.	SAP no.
NANOS femoral neck prosthesis size 5	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425005	75008159
NANOS femoral neck prosthesis size 6	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425006	75008160
NANOS femoral neck prosthesis size 7	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425007	75008161
NANOS femoral neck prosthesis size 8	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425008	75008162
NANOS femoral neck prosthesis size 9	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	425009	75008163
NANOS femoral neck prosthesis size 10 (*)	Ti6Al4V (ISO 5832-3) + TPS + BONIT®	75103109	75103109

(*) Special size

1.2 Overview of accessories

1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Name	Art. no.	SAP no.
Instrument Set NANOS Rasp Tray	-	75210420
Instrument Set NANOS Tray	-	75200170
Rasp handle double offset right hook coupling (*)	367-155	75103113
Rasp handle double offset left hook coupling (*)	367-156	75103114
Rasp adapter IMT double offset right hook coupling (*)	367-158	75103115
Rasp adapter IMT double offset left hook coupling (*)	367-159	75103116
Rasp adapter IMT spring hook NANOS (*)	506-1210	75094887
Rasp handle spring hook trigger mechanism medial NANOS (*)	506-1245	75030476
Rasp handle spring hook straight NANOS (*)	506-1204	75094886
Awl extended with handle, 90° rotation (*)	506-1220	75030323
Trial head Ø22 S/+0 (*)	75103040	75103040
Trial head Ø22 M/+4 (*)	75103041	75103041
Trial head Ø22 L/+8 (*)	75103042	75103042
Trial head Ø22 XL/+12 (*)	75103043	75103043
Trial head Ø28 XS/-3 (*)	75103044	75103044
Trial head Ø32 XS/-3 (*)	75103050	75103050
Trial head Ø36 XS/-3 (*)	75103056	75103056
NANOS compactor size 10 (*)	75103110	75103110
NANOS rasp size 10 (*)	75103111	75103111
Head End Ø22 / Ø26mm (*)	506-061	75029275
Head End Ø28 / Ø32mm (*)	506-062	75008166
Head End Ø36mm (*)	506-1212	75029345
Femoral Head Impactor Handle (*)	506-060	75008165

Name	Art. no.	SAP no.
Shaft extractor for taper 12/14 (*)	206-009	-
IMT adapter (*)	334-106	-
Extractor rod M8, L=430mm (*)	367-1250	-
Modular impact weight (*)	506-2107	75006921

(*) optional instruments

1.2.2 Accessories

Name	Art. no.	SAP no.
X-ray template femoral neck prosthesis D22	12330	-
X-ray template femoral neck prosthesis D28	09777	-
X-ray template femoral neck prosthesis D32	09778	-
X-ray template femoral neck prosthesis D36	09779	-

1.3 Accompanying reference documents

Name	Art. no.	SAP no.
Surgical technique femoral neck prosthesis	50000698	-
Instructions for use stem reducer	50000427	-
Implant Card	50000572	-
Patient information Symbol Meaning	50000837	-
Patient information Hip joint replacement	50000841	-

2. Handling

2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>, must be observed.

Caution: Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

Caution: the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint. The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

2.2 Authorised Component Combinations

Only femoral heads that are compatible with the system may be used for the assembly with the implant. The approved combinations with Smith & Nephew femoral heads are listed in a compatibility matrix. The matrix is available at: www.smith-nephew.com/compatibilitymatrix (Lit.No. 04758).

In general, the guidelines for use and handling as specified in the manufacturer's instructions for use of the endoprostheses must be followed and implemented.

Due to reasons relating the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

2.3 Instructions for Use

The application of the implant is done without cement. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads, the most up-to-date version of the respective instructions for use, which is available for download on the websites <https://www.ohst.de/ifu-instructions/> and <https://ifu.smith-nephew.com/>, must also be observed. After reduction, the definitive stability, mobility and muscle tension must be checked.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If an already implanted femoral head is to be replaced, for example in case of a revision, we specifically refer to the approved combinations with Smith&Nephew femoral heads. These are available at: www.smith-nephew.com/compatibilitymatrix (Lit.No. 04758). The

use of heads with a very long neck can result in increased mechanical stress and later stem fracture.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

The porous coating (TPS, BONIT®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

Caution: Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free gloves.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

If femoral heads with a neck (XL and XXL) are used, the range of motion is reduced by about 30° and achieves values between 80° and 100° in flexion and extension.

Caution: Component malposition or the use of a constrained hip or a total hip with a modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement, premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

2.4 Surgical Technique

Either the anterior-lateral or transgluteal approach can be used. Other approaches that are considered appropriate by the surgeon can also be used. The prosthesis type is also suitable for minimally invasive approaches. Resect the bone as sparingly as possible according to the planning. A subcapital resection of 0.5 to 1.0 cm with a straight oblique osteotomy of the femoral neck is usually sufficient. Prepare the prosthetic bed with the curved universal rasp that is used to prepare the path for the shaping rasps. Insert the rasp into the inner cortical bone using gentle rotations until you reach the lateral cortical bone near the lower edge of the lesser trochanter. Connect the rasp size 0 with the curved rasp handle for the minimally invasive preparation and advance it along the curvature of the calcar bone. This rasp can later be used to finely prepare the lateral cortical ring. After that, finish preparing the bone using the cancellous bone compactors. The aim of preparing the bone with compactors is to compress the cancellous bone. Prepare the bone step by step until you reach the planned size or until the compactor comes into contact with the cortical bone in the load-bearing areas. Mechanical stability and a change in the sound during

impaction signals contact with the cortical bone. Guide the rasp and the compactors using small circular motions to prevent creating free spaces that cannot be filled by the prosthesis. The compactor should end outside or be flush with the resection line. Use the next size if necessary. Once the suitable compactor is in the correct position, remove the rasp handle and attach the trial head. Carry out a trial reduction. Clinically ensure that there is unobstructed freedom of movement and check the telescoping and the leg length. Trials heads in the diameters 22*, 28, 32 and 36 mm with the neck lengths XS*, S, M, L and XL (*optional) are available. After reduction, it is recommended to use a C-arm X-ray unit to assess the fit, centre of rotation and offset of the implants. Remove the trial head and the compactor. If necessary, continue to prepare the prosthetic bed using the rasps and a larger compactor. On the lateral side, where a bony conflict can develop between the compactor and the cortical femoral neck ring, widen the bone using the awl until the compactor can be inserted as desired.

Manually insert the original NANOS prosthesis that corresponds to the size of the last compactor used. After that, drive in the prosthesis to the depth prepared by the compactor using the insertion instrument. Carefully clean the cone before attaching and driving in the original prosthetic head in the correct size. Check the reduction, stability and range of motion. The wound is closed in the usual way routinely used by the surgeon. During stem implantation it is advisable to leave the protective cap on the cone. This will protect the cone against damage. Should it become necessary during surgery to remove an implanted original NANOS prosthesis, an extraction instrument is available. The extraction instrument is attached to the cone and connected to the rasp handle for the removal.

The wound is closed in the usual way routinely used by the surgeon.

For further information see the NANOS surgical technique (see 1.3 Accompanying reference documents).

3. Packaging and Sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative Planning and Postoperative Care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognised procedures must be used for postoperative care.

5. Indications

- Primary and secondary osteoarthritis of the hip
- Dysplasia osteoarthritis
- Avascular and post-traumatic femoral head necrosis
- Patients with good bone density and intact femoral neck

Total hip replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total hip replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

6. Contraindications

- Marked osteoporosis
- Dysmorphia and loss of bone substance, as well as prior surgical procedures that do not guarantee the intended support
- Marked coxa valga
- Marked coxa vara with a femoral neck implant bed < 125°
- Intolerance to the materials

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Use of high-dose cortisone or cytotoxic drugs
- History or risk of infectious diseases with potential joint manifestation

- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Obesity (Body-Mass-Index BMI > 30)
- Anticipated extreme stresses (e. g. through sports or work), especially in patients with a body weight over 100 kg
- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Growth in children and adolescents
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour

8. Undesirable Effects

The undesirable effects and risks listed below can occur in the course of a total hip arthroplasty (THA) or hemiarthroplasty:

- Foreign body reactions, osteolysis, loosening
- ARMD, trunnionosis, metallosis
- Toxic reactions
- Luxation, dislocation, dissociation
- Implant failure
- Limited ROM
- Loosening
- Difference in leg length
- Bone fracture
- Instability
- Stem subsidence
- Tissue damage
- Heterotopic ossification
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack / stroke

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient Information and Documentation

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. Based on the current state of technology, the expected service life is approx. 15 years.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant.

Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

The patient can download the patient information from the website <https://www.ohst.de/patient-information/>.

A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

10. Implant Card

Following surgery, the patient must be given an implant card containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant card must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website. The patient's details (patient name or patient ID) must be documented in the implant card, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the card.

The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

11. Explanation of the Label Symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 149).

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