OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 (0)3385 - 54 20-0

Fax: +49 (0)3385 - 54 20-99

www.ohst.de

CE 0482

Primaro Pfanne	
Deutsch (de)	3
Primaro Cup	
English (en)	19
Cotyle Primaro	
Français (fr)	34
Cotile Primaro	
Italiano (it)	50
Cótilo Primaro	
Español (es)	66
Caixa Primaro	
Português (pt)	83

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IMPLANT

Primaro Cup

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for hip arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website <https://www.ohst.de/ifu-instructions/>.

1. Product Description and Implant Materials

The Primaro Cup is a modular acetabular hip cup for total hip arthroplasty (THA).

A Primaro Cup may only be used in patients who require total hip arthroplasty and who fulfil the relevant indications, taking into account the contraindications. In addition, the Primaro Cup may only be used to treat skeletally mature patients.

The Primaro Cup system consists of an outer shell made of Ti6Al4V alloy (ISO 5832-3) for cementless fixation in the acetabulum, an inlay made of highly cross-linked PE (XL-PE 75) or ceramic (ISO 6474-2 Al₂O₃+ZrO₂) and screw lids made of Ti6Al4V (ISO 5832 3). The outer shell is available in the variants rough-blasted surface, TPS coating and TPS+BONIT® coating. If necessary, self-tapping flat-head screws made of Ti6Al4V alloy (ISO 5832-3) with a cancellous bone screw thread of 6.5 mm diameter are available in varying lengths from 15 mm to 60 mm for additional fixation of the outer shell inside the acetabulum.

The outer shell comes in outside diameters of 44 mm to 68 mm. The internal contour of the outer shell is designed so that for every cup size, XL-PE inlays are used that are specially designed for this system. The XL-PE inlays are available in standard and dysplasia variants for femoral head diameters measuring 28 mm, 32 mm and 36 mm. In addition, ceramic inlays are available for femoral head diameters of 28 mm, 32 mm and 36 mm.

Ceramic inlays are inlays for clamping in the metallic hip cups and are made from ultrapure implant ceramic (ISO 6474-2 Al₂O₃ + ZrO₂). The ceramic is bioinert, biocompatible, biostable, mechanically stable and corrosion-resistant. It combines excellent resistance to fatigue, high shock resistance and tensile strength, excellent fracture toughness, an extremely high degree of hardness and no ions penetrate the body.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

1.1 Overview of Implants

Designation	Material	Reference number
Primaro Cup with screw fixation		
Primaro SF Cup size 44 (XX/37G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1216
Primaro SF Cup size 46 (XX/39G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1217
Primaro SF Cup size 48 (XX/39G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1218
Primaro SF Cup size 50 (XX/39G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1219
Primaro SF Cup size 52 (XX/44G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1220
Primaro SF Cup size 54 (XX/44G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1221
Primaro SF Cup size 56 (XX/48G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1222
Primaro SF Cup size 58 (XX/48G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1223
Primaro SF Cup size 60 (XX/52G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1224
Primaro SF Cup size 62 (XX/52G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1225
Primaro SF Cup size 64 (XX/52G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1226
Primaro SF Cup size 66 (XX/52G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1227
Primaro SF Cup size 68 (XX/52G) coating: TPS	ISO 5832-3 Ti6Al4V	367-1228
Primaro SF Cup size 44 (XX/37G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1229
Primaro SF Cup size 46 (XX/39G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1230
Primaro SF Cup size 48 (XX/39G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1231
Primaro SF Cup size 50 (XX/39G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1232
Primaro SF Cup size 52 (XX/44G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1233
Primaro SF Cup size 54 (XX/44G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1234
Primaro SF Cup size 56 (XX/48G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1235
Primaro SF Cup size 58 (XX/48G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1236
Primaro SF Cup size 60 (XX/52G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1237
Primaro SF Cup size 62 (XX/52G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1238
Primaro SF Cup size 64 (XX/52G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1239
Primaro SF Cup size 66 (XX/52G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1240
Primaro SF Cup size 68 (XX/52G) coating: TPS+BONIT®	ISO 5832-3 Ti6Al4V	367-1241
Primaro SF Cup size 44 (XX/37G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1534
Primaro SF Cup size 46 (XX/39G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1535
Primaro SF Cup size 48 (XX/39G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1536
Primaro SF Cup size 50 (XX/39G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1537
Primaro SF Cup size 52 (XX/44G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1538
Primaro SF Cup size 54 (XX/44G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1539
Primaro SF Cup size 56 (XX/48G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1540
Primaro SF Cup size 58 (XX/48G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1541
Primaro SF Cup size 60 (XX/52G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1542

Designation	Material	Reference number
Primaro SF Cup size 62 (XX/52G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1543
Primaro SF Cup size 64 (XX/52G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1544
Primaro SF Cup size 66 (XX/52G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1545
Primaro SF Cup size 68 (XX/52G), rough-blasted	ISO 5832-3 Ti6Al4V	367-1546
Primaro XL-PE Inlay		
XL-PE Inlay Primaro Cup, size 44 KD28 STD	XL-PE 75	367-1276
XL-PE Inlay Primaro Cup, size 46-50 KD28 STD	XL-PE 75	367-1277
XL-PE Inlay Primaro Cup, size 52-54 KD28 STD	XL-PE 75	367-1278
XL-PE Inlay Primaro Cup, size 56-58 KD28 STD	XL-PE 75	367-1279
XL-PE Inlay Primaro Cup, size 60-68 KD28 STD	XL-PE 75	367-1280
XL-PE Inlay Primaro Cup, size 46-50 KD32 STD	XL-PE 75	367-1281
XL-PE Inlay Primaro Cup, size 52-54 KD32 STD	XL-PE 75	367-1282
XL-PE Inlay Primaro Cup, size 56-58 KD32 STD	XL-PE 75	367-1283
XL-PE Inlay Primaro Cup, size 60-68 KD32 STD	XL-PE 75	367-1284
XL-PE Inlay Primaro Cup, size 52-54 KD36 STD	XL-PE 75	367-1285
XL-PE Inlay Primaro Cup, size 56-58 KD36 STD	XL-PE 75	367-1286
XL-PE Inlay Primaro Cup, size 60-68 KD36 STD	XL-PE 75	367-1287
XL-PE Inlay Primaro Cup, size 44 KD28 DYS	XL-PE 75	367-1288
XL-PE Inlay Primaro Cup, size 46-50 KD28 DYS	XL-PE 75	367-1289
XL-PE Inlay Primaro Cup, size 52-54 KD28 DYS	XL-PE 75	367-1290
XL-PE Inlay Primaro Cup, size 56-58 KD28 DYS	XL-PE 75	367-1291
XL-PE Inlay Primaro Cup, size 60-68 KD28 DYS	XL-PE 75	367-1292
XL-PE Inlay Primaro Cup, size 46-50 KD32 DYS	XL-PE 75	367-1293
XL-PE Inlay Primaro Cup, size 52-54 KD32 DYS	XL-PE 75	367-1294
XL-PE Inlay Primaro Cup, size 56-58 KD32 DYS	XL-PE 75	367-1295
XL-PE Inlay Primaro Cup, size 60-68 KD32 DYS	XL-PE 75	367-1296
XL-PE Inlay Primaro Cup, size 52-54 KD36 DYS	XL-PE 75	367-1297
XL-PE Inlay Primaro Cup, size 56-58 KD36 DYS	XL-PE 75	367-1298
XL-PE Inlay Primaro Cup, size 60-68 KD36 DYS	XL-PE 75	367-1299
Ceramic Inlay		
BIOLOX®delta Inlay XLW 18 28/37G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-894
BIOLOX®delta Inlay XLW 18 28/39G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-1307
BIOLOX®delta Inlay XLW 18 28/44G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-1308
BIOLOX®delta Inlay XLW 18 28/48G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-1309
BIOLOX®delta Inlay XLW 18 28/52G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-1310
BIOLOX®delta Inlay XLW 18 32/39G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-895
BIOLOX®delta Inlay XLW 18 32/44G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-896

Designation	Material	Reference number
BIOLOX®delta Inlay XLW 18 32/48G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-897
BIOLOX®delta Inlay XLW 18 32/52G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-1311
BIOLOX®delta Inlay XLW 18 36/44G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-898
BIOLOX®delta Inlay XLW 18 36/48G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-899
BIOLOX®delta Inlay XLW 18 36/52G	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	367-1312
ELEC®plus Inlay Ø28/37-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-012
ELEC®plus Inlay Ø28/39-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-013
ELEC®plus Inlay Ø28/44-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-014
ELEC®plus Inlay Ø28/48-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-015
ELEC®plus Inlay Ø28/52-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-016
ELEC®plus Inlay Ø32/39-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-017
ELEC®plus Inlay Ø32/44-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-018
ELEC®plus Inlay Ø32/48-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-019
ELEC®plus Inlay Ø32/52-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-020
ELEC®plus Inlay Ø36/44-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-021
ELEC®plus Inlay Ø36/48-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-022
ELEC®plus Inlay Ø36/52-18	ISO 6474-2 Al ₂ O ₃ +ZrO ₂	013-023

1.2 Overview of accessories

1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Name	Reference number
Acetabulum reamer instrument set	367-147
Equatorial size tester instrument set	367-1464
Cementless cup instrument set	367-149
Primaro Cup instrument set	367-1341
Optional:	
Primaro Cup insertion instrument set (not to be used in combination with ceramic inlays)	367-1477
CERA-Guide insertion instrument set, Primaro Cup (not to be used with ELEC®plus inlays)	367-1479
Cup Positioning Guide Ø8 to Ø24	506-010
Primaro trial cups instrument set	367-1342
Primaro Trial Insert DYS Ø 44 (28/37G) PPSU green, with pole screw	367-1494
Primaro Trial Insert DYS Ø 46-50 (28/39G) PPSU green, with pole screw	367-1495
Primaro Trial Insert DYS Ø 52-54 (28/44G) PPSU green, with pole screw	367-1496
Primaro Trial Insert DYS Ø 56-58 (28/48G) PPSU green, with pole screw	367-1497
Primaro Trial Insert DYS Ø 60-68 (28/52G) PPSU green, with pole screw	367-1498
Primaro Trial Insert DYS Ø 46-50 (32/39G) PPSU blue, with pole screw	367-1499

Name	Reference number
Primaro Trial Insert DYS Ø 52-54 (32/44G) PPSU blue, with pole screw	367-1500
Primaro Trial Insert DYS Ø 56-58 (32/48G) PPSU blue, with pole screw	367-1501
Primaro Trial Insert DYS Ø 60-68 (32/52G) PPSU blue, with pole screw	367-1502
Primaro Trial Insert DYS Ø 52-54 (36/44G) PPSU yellow, with pole screw	367-1503
Primaro Trial Insert DYS Ø 56-58 (36/48G) PPSU yellow, with pole screw	367-1504
Primaro Trial Insert DYS Ø 60-68 (36/52G) PPSU yellow, with pole screw	367-1505

1.2.2 Accessories

Designation	Material	Reference number
X-ray template Primaro Cup	-	367-1338

1.3 Accompanying reference documents

Designation	Reference number
Surgical technique Primaro Cup	50000349
Implant Card	50000572
Instructions for use of insertion instrument	200815300886
Instructions for use of flat-head screw	50000585
Patient information Meaning of symbols	50000837
Patient information Hip joint replacement	50000841

2. Handling

2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>, must be observed.

Caution: Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

Caution: the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint.

The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

2.2 Authorised Component Combinations

Only the following flat-head screws must be used for screw fixation:

Designation	Material	Reference number
Flat-head screw Ø 6.5 x 15 mm	ISO 5832-3 Ti6Al4V	000-290-15
Flat-head screw Ø 6.5 x 20 mm	ISO 5832-3 Ti6Al4V	000-290-20
Flat-head screw Ø 6.5 x 25 mm	ISO 5832-3 Ti6Al4V	000-290-25
Flat-head screw Ø 6.5 x 30 mm	ISO 5832-3 Ti6Al4V	000-290-30
Flat-head screw Ø 6.5 x 35 mm	ISO 5832-3 Ti6Al4V	000-290-35
Flat-head screw Ø 6.5 x 40 mm	ISO 5832-3 Ti6Al4V	000-290-40
Flat-head screw Ø 6.5 x 45 mm	ISO 5832-3 Ti6Al4V	000-290-45
Flat-head screw Ø 6.5 x 50 mm	ISO 5832-3 Ti6Al4V	000-290-50
Flat-head screw Ø 6.5 x 55 mm	ISO 5832-3 Ti6Al4V	000-290-55
Flat-head screw Ø 6.5 x 60 mm	ISO 5832-3 Ti6Al4V	000-290-60

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

Caution: The ELEC®plus inlays may only be used in combination with ELEC®plus femoral heads. The BIOLOX®delta inlays may only be combined with BIOLOX®delta femoral heads.

The combination with a different femoral head or with a different ceramic ball of other manufacturers is prohibited.

2.3 Instructions for Use

The implant is used without cement.

The clamping cone connection allows for a secure fixation of the inlay inside the implanted hip cup. It ensures a uniform transfer of load to the cup. Care must be taken to use only an unused and intact inlay which was removed from its original packaging immediately prior to insertion. An inlay that has been used before must not be reused. Since an exact fit between the ceramic inlay and the cup housing is essential, only new and intact prosthetic components may be combined. This also means, for example, that an inlay that was already inserted into a cup housing before and was then removed again must not be reinserted into a cup housing. By the same logic, an inlay that is in any way damaged must not be reused and instead must be discarded. This also applies, for example, to an inlay that has been dropped. An inlay must only be used with a brand-new, unused cup housing that is not damaged.

Caution: The ceramic inlays must fit the implanted cup perfectly. Especially when using press-fit fixation systems, minor deformation of the cup following implantation into the osseous bed can occur, which can in turn make it more difficult to achieve the correct fixation of a ceramic inlay.

Caution: The ball diameter of the femoral head prosthesis must agree perfectly with the spherical nominal diameter of the part of the acetabular cup implant which articulates as the sliding partner.

Caution: The positioning of the implant components has a direct influence on the range of motion, and therefore potentially on the risk of impingement. An angle of inclination that is either too steep or too flat increases the load on the edge of the cup. This can lead to increased wear. In addition, movements may be restricted, which in turn may lead to subluxation and dislocation of the head from the ceramic inlay.

Caution: The anteversion of the cup housing should not have any negative values. This can result in increased wear. In addition, movements may be restricted, which in turn may lead to subluxation and dislocation of the head from the ceramic inlay.

Caution: Do not use a ceramic inlay if the cup housing is in a position that is not recommended in the two paragraphs above (with the exception of special systems intended for this purpose, with elements to prevent dislocation, e.g. asymmetrical or tilted inserts).

Caution: In the case of retroverted cup housings, a ceramic inlay must not be used. Otherwise, mechanical malfunctioning of the inlay may occur. Excessive wear of the ceramic may result in intolerances or loosening of the prosthesis.

Caution: Revision using a ceramic inlay with a cup housing that remains in situ is contraindicated. In this case, in general, an inlay made of plastic can be used.

Caution: In very rare cases, the ceramic inlay may fracture. Fractures may be caused, among other things, by incorrect fixation of the ceramic inlay inside the metallic cup housing or by an incorrect or absent fit of the inlay with the metallic cup housing. Prosthetic components that have not been approved by the endoprosthesis manufacturer to be combined with a ceramic inlay may also lead to fracture of the inlay. The same applies

if the recommended positioning of the cup housing is disregarded. In the very rare event of a fracture of a ceramic component, synovectomy may have to be carried out.

Caution: In the very rare event that a ceramic insert fractures, there is a risk of injury from the sharp edges of ceramic fragments during revision surgery.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

The porous coating (TPS, BONIT®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

Caution: Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free gloves.

Caution: when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

When using Snap or Dysplasia cups or cup inserts, the movement range in flexion and extension is reduced by approximately 60° and in abduction and adduction by approximately 22° as compared to the standard application.

Caution: Component malposition or the use of a constrained hip or a total hip with a modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement, premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

2.4 Surgical Technique

The resection of the femoral neck is carried out in accordance with the preoperative planning of the hip stem.

The acetabulum is prepared using cup reamers; the last reamer size used must be the same as the preoperative nominal size of the implant. Never apply force, always allow the reamer to move freely! In particular, the final reaming must not be distorted as a result of lateral pressure. The bony mass from the last reaming can be reserved to fill in any potential gaps between the implant and the acetabulum.

Caution: When implanting the Primaro Cup only the acetabulum reamers from size 44 to 68 mm must be used.

Equatorial size testers are available to determine the size of the cup needed. The size of the implant must be identical to the size of the last equatorial size tester used that fitted accurately. Trial cups can optionally be used for this. The outside diameter of the outer shell of the Primaro Cup is oversized compared to the size indication of the acetabulum reamer to achieve the desired primary stability. This oversize is also taken into account in the trial cups. The size of the implant must absolutely correspond to the size of the last trial cup used.

To implant the cup, connect the insertion head that fits the previously determined implant size with the straight handle or the MIS handle. After that, retract the spreader with the hook and position the cup onto the insertion head. The cup is locked when the hook is removed. Use moderate hammer blows to insert the cup into the pre-reamed dome. The insertion heads are not intended to be used for removing the cup. The cup is optimally positioned inside the acetabulum at an inclination of 35-45° and an anteversion of 10°. This is merely a recommendation. The individual anatomical conditions have to be taken into account for the exact alignment. If necessary, defects in the bony area are to be filled with cancellous bone. Alignment rods for the MIS handle are available as a guide for positioning the cup.

Caution: If the Primaro Cup requires fixation with screws owing to the condition of the bone, the user must ensure that the screw holes point towards the cranio-lateral side. Markings on the edge of the cup indicate the position of the screw holes for better guidance.

Caution: When inserting the Primaro Cup in combination with ceramic inlays, only impactors that support the cone from the inside must be used. In this case, driving in the cup with an insertion instrument that is connected to the cup only via the pole drill-hole is not permitted!

In order to place the screws, the covers of the screw holes (screw lids) have to be removed from the cup with the Cardan screwdriver. When using screws, pre-drill with a drill bit of Ø3.2 mm. A drilling gauge is available for the alignment of the drill.

A Cardan screwdriver is provided for insertion of the screws.

Caution: Only the flat-head screws (see 2.2 Authorised Component Combinations) belonging to this system must be used for fixation of the screws.

Caution: The drill-holes and screws must be placed in such a way as to avoid damage of the pelvic vessels!

The femoral stem is prepared according to the surgical instructions of the hip stem system being used. Trial inlays corresponding to the inlay implant sizes that must be used together with the previously implanted cup are available for trial repositioning of the entire joint. First position the trial inlay into the cup manually and then tighten the pole screw with a screwdriver. In combination with a rasp, trial head and trial inlays, the function of the entire joint can be assessed. Afterwards, all trial components must be removed and the inlay inserted.

Before inserting the inlay, the cup must be thoroughly rinsed and then dried. The inlay is then driven in with the insertion instrument. Insertion of the PE inlay is completed with an audible snap into the cup. The snap mechanism of the PE inlay is deformed, ensuring a tight fit inside the cup.

Dysplasia inlays have a raised edge for positioning in the cranial portion of the acetabulum.

Caution: When correctly installed, the PE inlay protrudes approximately 1.3 mm beyond the cup edge. In the event that the PE inlay ever needs to be removed, it must never be driven

in a second time. This applies to both the standard and the dysplasia variants of the inlay.

Position and insert the ceramic inlays into the cup; when positioned correctly, the inlay is flush with the edge of the cup. Place the ceramic inlay with the impactor and corresponding impactor head according to the manufacturer's instructions. The CERA-Guide insertion instrument set is optionally available for ceramic inlays BIOLOX®delta.

Caution: When a ceramic inlay is used, only a ceramic femoral head may be used as an articulation partner.

The insertion of the cup is concluded and the latter can be covered by small sterile compresses to avoid damage during the further course of the surgical procedure. The surgical procedure is continued with the implantation of a hip stem as per the surgical instructions of the hip step system being used. After implantation of the hip stem, the femoral head has to be reduced into the acetabulum. If a compress was used for covering, it has to be removed again prior to the reduction. The operation is routinely completed with the layer-by-layer wound closure.

Instruction for using optional Primaro Cup insertion instruments via pole drill hole

The insertion instrument must be screwed as far as it will go into the drill hole on the pole of the cup. Use moderate hammer blows to insert the cup into the pre-reamed dome. If necessary, defects in the bony area are to be filled with cancellous bone. Refer to the surgical technique for the cup with respect to optimally positioning the cup in the acetabulum. However, for exact alignment the patient's anatomy must be taken into account.

The Primaro Cup MIS insertion instrument with screwdriver is designed for the anterior surgical approach.

Caution: The Primaro Cup insertion instruments via pole drill hole may be used only in the case of PE inlays. Use with ceramic inlays is not permitted.

Additional instructions for use to ensure the safe functioning of ELEC®plus inlays

Note the following points before using the ELEC®plus inlay:

1. Depending on the indication, ensure an inclination (40° – 45°) and an anteversion (10° – 20°) of the hip cup the
2. Thoroughly rinse the clamping cone of the implanted hip cup with water to remove any tissue fragments, bone splinters or cement residue
3. Dry the clamping cone of the implanted hip cup and check the position and function of the joint by inserting a trial implant; make sure that any screws that may have been used have been screwed in sufficiently deep
4. Remove the trial implant and rinse and dry it again
5. Hold one part of the rim of the ELEC®plus inlay using your thumb and index finger and carefully allow it to slide into the clamping cone of the hip cup; do not apply any additional pressure in this state
6. Check the correct flush position of the ELEC®plus inlay by running a finger along the rim of the cup. In doing so, check if the edge of the face of the cup and of the face of

the ELEC®plus inlay are flush; if necessary, correct any uneven areas by carefully applying targeted pressure

7. Once this is done, firmly press the ELEC®plus inlay centrally into the dome using your thumb and wedge it into the clamping cone. Again check that it is positioned flush by running your finger along the rim and - if required - reposition it using the instrument authorised by the manufacturer of the endoprosthesis system.

Additional instructions for use to ensure the safe functioning of BIOLOX®delta inlays

As a vital prerequisite for safe functioning, maximum care must be taken when inserting the BIOLOX®delta inlay into the cup housing. In order to ensure a central mounting of the BIOLOX®delta inlay in the metal cup, the following procedure should be followed during installation:

The metal cup should be implanted into the acetabulum in accordance with the surgical technique. Trial implants should then be used to verify the correct fit of the implant and the movements of the joint. It must be ensured that any screws used are completely recessed in the metal cup. Afterwards, the trial implant must be removed and the metal cup rinsed and dried thoroughly. There must be no bone or tissue residues remaining in the clamping area of the metal cup. The BIOLOX®delta inlay is inserted by hand or with the insertion instrument. Whenever the insertion instrument is used, the respective instructions for use (see 1.3 Accessories) must be complied with. When inserted by hand, the inlay is held with two fingers and inserted into the cup housing. The positioning of the BIOLOX®delta inlay is checked by feeling the edge of the cup and correcting it if necessary. The connection between the edge of the metal and edge of the ceramic must be flush! Once the BIOLOX®delta inlay is positioned correctly, press down with the thumb on the base of the insert. To fix the implant in its final position, tightly clamp the BIOLOX®delta inlay into the cup housing with one or more moderate hammer blows on the impactor in the direction of the axis of symmetry using an impactor that is suitable for BIOLOX®delta inlays and recommended by the endoprosthesis manufacturer.

Caution: Never use a metal hammer to hit the BIOLOX®delta inlay or the ceramic insert!

3. Packaging and Sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner

bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative Planning and Postoperative Care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant. Recognised procedures must be used for postoperative care.

5. Indications

Significant damage to the hip joint as a result of:

- Primary or secondary osteoarthritis (coxarthrosis) or
- Rheumatoid arthritis
- Dislocated femoral neck fracture
- Avascular necrosis of the femoral head (ARCO stage > II)
- Previous unsuccessful procedures, such as non-union of fractures of the hip cup or femoral neck after osteosynthesis. (Complications as a result of orthopaedic endoprostheses, implants)

Total hip replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total hip replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

6. Contraindications

- Active infection of the hip joint to be operated on
- Active systemic infection
- Active infection of the skin in the area to be operated on
- Allergy to the implant materials
- Paralysis or lasting or non-reversible weakening of the musculature, which makes walking impossible (in the absence of pain)
- Rapidly progressing neurological disease (does not apply for treatment of dislocated femoral neck fracture)

Additional contraindications when using the ceramic insert:

- Revision involving leaving the stem and/or hip cup in situ (risk of component fracture due to damaged fixation cones)
- Inclination of more than 40°–45°, anteversion of more than 10°–20° (risk of dislocation, subluxation, impingement)
- Retroversion

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Comorbidities associated with an increased mortality risk
- BMI $\geq 30\text{kg/m}^2$ or $< 18.5\text{kg/m}^2$
- Smoking
- Diabetes mellitus
- Mental illnesses
- Anaemia
- Intra-articular corticosteroid injections in the affected joint three months before the planned joint replacement
- Generally increased perioperative risk
- Younger age is associated with higher revision rates

Risk factors and conditions specific to hip arthroplasty:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Muscle spasms or other spastic conditions
- Expected extreme stresses, e.g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities that complicate fixation of the implant
- Weakening of the supporting structures due to a tumour
- Concomitant pain in the region of the lower back

8. Undesirable Effects

The undesirable effects and risks listed below can occur in the course of a total hip replacement:

- Noise generated due to a ceramic on ceramic articulation
- Foreign body reactions, osteolysis, loosening
- Toxic reactions
- Sensitisation
- Limited ROM
- Instability
- Luxation, dislocation, dissociation
- Implant failure
- Bone fracture
- Loosening
- Drifting/ tilting
- Heterotopic ossification
- Difference in leg length

- Tissue damage
- Iliopsoas syndrome / irritation
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack / stroke

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient Information and Documentation

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. Based on the current state of technology, the expected service life is approx. 15 years.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant.

Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

The patient can download the patient information from the website <https://www.ohst.de/patient-information/>.

A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

10. Implant passport

Following surgery, the patient must be given an implant passport containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant passport must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website.

The patient's details (patient name or patient ID) must be documented in the implant passport, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the passport. The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

11. Explanation of the Label Symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 185).

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