OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Deutschland / Germany / Allemagne /

Germania / Alemania / Duitsland /

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IMPLANT

SPC Femoral Stem

Before using the product, the user is under obligation to carefully study and comply with the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance.

These implants may be used only by physicians with appropriate experience, practice and skills in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

1. Product description and implant materials

The SPC femoral stem is an endoprosthesis hip replacement for cementless anchorage in the femur. The principle is based on the findings of Prof. Lorenzo Spotorno and has been tried and tested for more than 20 years.

Produced from titanium wrought alloy (ISO 5832-11) and fitted with a 12/14 cone, the stem can be combined with both metal and ceramic femoral heads. For the various femoral anatomies of the patients the SPC femoral stem is available in 13 sizes, each of which has three different CCD angles 125°, 135° und 145°.

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Reference number
SPC Femoral Stem size 5 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-205
SPC Femoral Stem size 6 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-206
SPC Femoral Stem size 7 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-207
SPC Femoral Stem size 8 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-208
SPC Femoral Stem size 9 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-209
SPC Femoral Stem size 10 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-210
SPC Femoral Stem size 11.25 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-211
SPC Femoral Stem size 12.5 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-212
SPC Femoral Stem size 13.75 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-213
SPC Femoral Stem size 15 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-215
SPC Femoral Stem size 16.25 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-216
SPC Femoral Stem size 17.5 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-217
SPC Femoral Stem size 20 CCD 125° cementless	ISO 5832-11 Ti6Al7Nb	271-220

Name	Material	Reference number
SPC Femoral Stem size 5 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-105
SPC Femoral Stem size 6 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-106
SPC Femoral Stem size 7 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-107
SPC Femoral Stem size 8 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-108
SPC Femoral Stem size 9 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-109
SPC Femoral Stem size 10 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-110
SPC Femoral Stem size 11.25 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-111
SPC Femoral Stem size 12.5 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-112
SPC Femoral Stem size 13.75 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-113
SPC Femoral Stem size 15 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-115
SPC Femoral Stem size 16.25 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-116
SPC Femoral Stem size 17.5 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-117
SPC Femoral Stem size 20 CCD 135° cementless	ISO 5832-11 Ti6Al7Nb	271-120
SPC Femoral Stem size 5 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-005
SPC Femoral Stem size 6 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-006
SPC Femoral Stem size 7 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-007
SPC Femoral Stem size 8 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-008
SPC Femoral Stem size 9 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-009
SPC Femoral Stem size 10 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-010
SPC Femoral Stem size 11.25 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-011
SPC Femoral Stem size 12.5 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-012
SPC Femoral Stem size 13.75 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-013
SPC Femoral Stem size 15 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-015
SPC Femoral Stem size 16.25 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-016
SPC Femoral Stem size 17.5 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-017
SPC Femoral Stem size 20 CCD 145° cementless	ISO 5832-11 Ti6Al7Nb	271-020

12 Instrument overview

The instruments of OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Reference number
Basic Instrument Set, modular	367-116
Rasp Instrument Set SPC Femoral Stem	271-324
Femoral Stem Extraction Instrument Set	206-010

13 Accessories

Name	Reference number
Surgical Technique SPC	50000566
X-ray Template SPC CCD 125°	271-125

Name	Reference number
X-ray Template SPC CCD 135°	271-135
X-ray Template SPC CCD 145°	271-145
Instructions for use stem reducer	50000427
Instructions for Use Extraction Instrument Set	50000428
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance with the appropriate hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must therefore not come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on the functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. After explantation, previous damages which exclude any reuse must therefore be assumed.

2.2 Authorised component combinations

The following femoral heads must be used for assembly with the implant:

Nominal Ø	Name	Neck length	Reference number
Ø22 mm	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L	-4 / 0 / +4 mm	030-2200 to 030-2202
Ø28 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL / XXL	-3,5 / 0 / +3,5 / +7 / +10,5 mm	020-2800 to 020-2804
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL / XXL	-3,5 / 0 / +3,5 / +7 / +10,5 mm	030-2800 to 030-2804

Nominal Ø	Name	Neck length	Reference number
	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-907 to 367-909
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-1140 to 367-1142
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	384-001 to 384-003
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	110230, 110240, 110250
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	013-001 to 013-003
Ø32 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL / XXL	-4 / 0 / +4 / +8 / +12 mm	020-3200 to 020-3204
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL / XXL	-4 / 0 / +4 / +8 / +12 mm	030-3200 to 030-3204
	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	367-920 to 367-912
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	367-1143 to 367-1145, 367-1149
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	384-004 to 384-006
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	110260, 110270, 110280, 110291
	ELEC®plus Femoral Head ISO 6474-, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	013-004 to 013-007
Ø36 mm	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	367-930 to 367-932
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	367-1146 to 367-1148, 367-1150

Nominal Ø	Name	Neck length	Reference number
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	384-007 to 384-009
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	110300, 110310, 110320, 110330
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	013-008 to 013-011
Ø40 mm to Ø60 mm	Fracture head size S implant steel ISO 5832-9, 12/14	-4 mm	155-140 to 155-160
	Fracture head size M implant steel ISO 5832-9, 12/14	0 mm	155-040 to 155-060
	Fracture head size L implant steel ISO 5832-9, 12/14	+4 mm	155-240 to 155-260

We only guarantee compatibility of our products in combination with our own CE-marked products and with the products we have approved for combined use, for which the competent authority has issued its approval. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

2.3 Information for use

The application of the implant is done without cement. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads the corresponding instruction for use is also to be considered. After reduction, the definitive stability, mobility and muscle tension must be checked.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available. Before inserting the implant, the implant bed must be irrigated sufficiently. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

The porous coating (TPS, Bonit®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

After opening the joint capsule and dislocating the femoral head from the acetabulum, the latter is resected as determined in the preoperative planning and the femoral head is removed in full.

The height of the osteotomy is determined by the anteversion of the femoral neck: the stronger the anteversion, the lower the resection plane. It has proven advantageous to retain 1 to 1.5 cm of the femoral neck to ensure that the proximal ribbed section can adjust better.

Then open the medullary cavity using the cavity chisel. The chisel has to be applied to the far lateral and dorsal side in order to facilitate the subsequent driving in of the rasps in the direction of the femoral axis. The cavity chisel is placed relative to the desired antetorsion. Fracturing of the greater trochanter must be avoided. The medullary cavity opening can be extended with the femoral rasp

Once the medullary cavity has been opened, the femoral stem is rasped. The first rasp determines the orientation of the subsequent sizes. Care must therefore be taken to ensure correct antetorsion already during the first rasping. The antetorsion can be reviewed by means of the guide rod and rasp handle and is usually 10°-15°. The procedure is started with the smallest rasp size, which is connected to the rasp handle. After that, the subsequent rasps are used in ascending order until the preoperatively determined size is reached. The sizes of the rasps conform to the sizes of the implants. The correct position of the rasp inside the femur can be verified with an image converter. In order to guarantee a press-fit, the contours of the rasps and the implants are identical, with the exception of the ribs. The ribs are not prepared in advance. If, during the rasping procedure, the stem size does not correspond to the preoperatively determined size and if there is a difference of two or more sizes, the alignment of the axis may be incorrect or there may be an osseous obstruction present. In this case, the selected stem might be too small and would therefore not provide the required stability.

The design of the rasps is such that they can also be used for trial reductions. The leg length, the range of motion and the tension of the ligaments can be checked by attaching the various trial cones and trial heads. Remove the handle from the rasp, which remains inside the femur, and attach the trial cone to the rasp. Trial cones with three different angles are offered for trial reductions of the various stem variants. The trial cones snap into place as soon as they are in their correct position. Then attach the trial head to the trial cone by hand. Trial heads with various diameters are offered in the neck lengths S to XXL for trial reductions. After reduction, the definitive stability, mobility and muscle tension should be checked.

The stem size to be implanted is selected on the basis of the last rasp size used. The implant is positioned centrally by hand and inserted as deep into the femur as possible.

Caution: During the insertion process, the implant has to be in a central position.

After that, insert the prosthetic device into the implant bed up to the resection plane using the stem impactor. After that, thoroughly clean and dry the prosthetic cone and attach the femoral head with the previously determined neck length (S to XXL).

Caution: No femoral heads with a neck length greater than XXL may be used!

Then reduce the stem with the femoral head into the cup and check the range of motion, leg length and tension of the ligaments. The operation is routinely completed with the layer-by-layer wound closure.

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar data is indispensable and provides important information about suitable implants, placement, possible component combinations and enables the size of the implant to be used to be preselected. The operation may be performed only if the patient's material compatibility has been determined. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available in case another size is required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

- Advanced wear of the hip joint due to degenerative, post-traumatic or rheumatoid arthritis (coxarthrosis)
- Fracture or avascular necrosis of the femoral head
- Hip dysplasia
- Sequelae of earlier surgical procedures, e.g. osteosynthesis, articular reconstruction, arthrodesis
- Young and active patients with good bone quality
- Trumpet-shaped femoral morphology

- This system may also be indicated for the revision of earlier surgical attempts which have been unsuccessful

Nevertheless, the surgical intervention is more likely to succeed if certain criteria are met. In 1985 an indication protocol was drawn up which awards points for various parameters on the basis of clinical and radiological examinations. For more information about the indication protocol see the surgical technique for the SPC Femoral Stem (1.3. Accessories).

6. Contraindications

The following direct contraindications result from the indication protocol described above:

- Acute or chronic, local or systemic infection
- Severe muscle, nerve or vascular diseases endangering the extremity concerned;
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Severe osteoporosis, severe metabolic disorders with effects on the bone metabolism
- Severe hip dysplasia
- Loss of the ligamenture
- Bone morphologies which are not suitable for proximal anchorage (morpho-cortical index MCI 0-4 points)
- Any underlying condition that might compromise the function of the implant
- Hypersensitivity to the materials used

If femoral heads with a neck (XL and XXL) are used, the range of motion is reduced by about 30° and achieves values between 80° and 100° in flexion and extension.

7. Risks and conditions that can impair the success of the operation

Potential risks associated with the operation are:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fracture of the femoral bone
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Growth in children and adolescents
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

8. Undesirable effects

The adverse events listed below are among the most typical and commonly occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Muscle spasm
- Stiffness,
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Inflammation
- Oedem/fluids
- Metallosis
- Elevated metal ions in blood
- Coxa vara
- Osteolysis
- Heterotopic ossification
- Pseudotumours

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, even though the contraindications described above partially apply to the patient, the patients should be informed about the possible effects of these circumstances and the anticipated risks. It should be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. It is also recommended that the patient is informed of activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 192).

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