# OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 3385 5420-0

Fax: +49 3385 5420-99

www.ohst.de

CE 0482

---

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## IMPLANT

### Expersus Femoral Stem

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for hip arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website <https://www.ohst.de/ifu-instructions/>.

#### 1. Product Description and Implant Materials

The Expersus hip stem is an endoprosthesis (femoral component) for cementless or cemented anchorage in total hip arthroplasty (THA) or hip hemiarthroplasty.

The cementless variant is manufactured from titanium alloy Ti6Al4V (ISO 5832-3) and the cemented variant from implant steel (ISO 5832-9). The Expersus femoral stem comes with a 12/14 cone. The stem can be combined with both metal and ceramic femoral heads. Apart from the polished neck area, the cementless Expersus femoral stem has an all-round titanium plasma spray and hydroxyapatite coating. The body of the cemented femoral stem has a satin-finished, lightly blasted surface.

To account for the various femoral anatomies of the patients, the femoral stem comes in 9 sizes with a standard 125° and a standard 135° variant for cementless anchorage as well as in 9 sizes with a standard 125° and in 10 sizes with a standard 135° variant for cemented anchorage.

The Expersus femoral stem can generally be used in patients in whom total hip replacement surgery is indicated. The Expersus femoral stem may only be used to treat skeletally mature patients.

All implant variants of the Expersus femoral stem can be inserted with the same instrument set.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

#### 1.1 Overview of Implants

Designation	Material	Reference number
Expersus femoral stem cemented		
Expersus femoral stem size 1 STD, 135°, cemented (*)	ISO 5832-9 Implant steel	367-1556
Expersus femoral stem size 2 STD, 135°, cemented (*)	ISO 5832-9 Implant steel	367-1557
Expersus femoral stem size 3 STD, 135°, cemented (**)	ISO 5832-9 Implant steel	367-1558
Expersus femoral stem size 4 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1559
Expersus femoral stem size 5 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1560

Designation	Material	Reference number
Expersus femoral stem size 6 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1561
Expersus femoral stem size 7 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1562
Expersus femoral stem size 8 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1563
Expersus femoral stem size 9 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1564
Expersus femoral stem size 10 STD, 135°, cemented	ISO 5832-9 Implant steel	367-1565
Expersus femoral stem size 2 STD, 125°, cemented (*)	ISO 5832-9 Implant steel	367-1547
Expersus femoral stem size 3 STD, 125°, cemented (**)	ISO 5832-9 Implant steel	367-1548
Expersus femoral stem size 4 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1549
Expersus femoral stem size 5 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1550
Expersus femoral stem size 6 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1551
Expersus femoral stem size 7 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1552
Expersus femoral stem size 8 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1553
Expersus femoral stem size 9 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1554
Expersus femoral stem size 10 STD, 125°, cemented	ISO 5832-9 Implant steel	367-1555

(\*) Weight limited to 65 kg (see Contraindications)

(\*\*) Weight limited to 90 kg (see Contraindications)

Expersus femoral stem cementless		
Expersus femoral stem size 2 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1400
Expersus femoral stem size 3 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1401
Expersus femoral stem size 4 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1402
Expersus femoral stem size 5 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1403
Expersus femoral stem size 6 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1404
Expersus femoral stem size 7 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1405
Expersus femoral stem size 8 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1406
Expersus femoral stem size 9 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1407
Expersus femoral stem size 10 STD, 135°, cementless	ISO 5832-3 Ti6Al4V	367-1408
Expersus femoral stem size 2 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1409
Expersus femoral stem size 3 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1410
Expersus femoral stem size 4 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1411
Expersus femoral stem size 5 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1412
Expersus femoral stem size 6 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1413
Expersus femoral stem size 7 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1414
Expersus femoral stem size 8 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1415

Expersus femoral stem size 9 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1416
Expersus femoral stem size 10 STD, 125°, cementless	ISO 5832-3 Ti6Al4V	367-1417

## 1.2 Overview of accessories

### 1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
Expersus femoral stem instrument set	367-1576
Expersus femoral stem rasp instrument set	367-1492
Expersus femoral stem optional instrument set	367-1577

### 1.2.2 Accessories

Designation	Reference number
X-ray template Expersus femoral stem cementless STD 135° KD28	367-246
X-ray template Expersus femoral stem cementless STD 125° KD28	367-247
X-ray template Expersus femoral stem cementless STD 135° KD32	367-074
X-ray template Expersus femoral stem cementless STD 125° KD32	367-075
X-ray template Expersus femoral stem cemented STD 135° KD28	367-244
X-ray template Expersus femoral stem cemented STD 125° KD28	367-245
X-ray template Expersus femoral stem cemented STD 135° KD32	367-072
X-ray template Expersus femoral stem cemented STD 125° KD32	367-073
Cement restrictor size 1, external Ø24mm	506-101
Cement restrictor size 2, external Ø27mm	506-102

## 1.3 Accompanying reference documents

Designation	Reference number
Surgical technique Expersus femoral stem cementless	50000351
Surgical technique Expersus femoral stem cemented	50000352
Instructions for use Extraction Instrument Set	50000428
Instructions for use Stem Repositioner	50000427
Implant Card	50000572
Patient information Hip joint replacement	50000841
MRI safety information	50000851
List of symbols	50000859

## 2. Handling

### 2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for

use (50000354), which is available for download on the website <https://www.ohst.de/fu-instructions/>, must be observed.

**Caution:** Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

**Caution:** the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint. The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

## **2.2** Authorised Component Combinations

The following femoral heads must be used for assembly with the implant:

Nominal Ø	Name	Neck length	Reference number
Ø22 mm	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L	-4 / 0 / +4 mm	030-2200 to 030-2202
Ø28 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL	-3,5 / 0 / +3,5 / +7 mm	020-2800 to 020-2803

Nominal Ø	Name	Neck length	Reference number
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL	-3,5 / 0 / +3,5 / +7 mm	030-2800 to 030-2803
	BIOLOX®delta Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-1140 to 367-1142
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	013-001 to 013-003
Ø32 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	020-3200 to 020-3203
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	030-3200 to 030-3203
	BIOLOX®delta Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	367-1143 to 367-1145, 367-1149
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	013-004 to 013-007
Ø36 mm	BIOLOX®delta Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	367-1146 to 367-1148, 367-1150
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	013-008 to 013-011
Ø40 mm to Ø60 mm	Fracture head size S implant steel ISO 5832-9, 12/14	-4 mm	155-140 to 155-160
	Fracture head size M implant steel ISO 5832-9, 12/14	0 mm	155-040 to 155-060
	Fracture head size L implant steel ISO 5832-9, 12/14	+4 mm	155-240 to 155-260

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

## 2.3 Instructions for Use

The application of the implant is done with or without cement. A cement restrictor must be used for cemented applications. Before using the cement restrictor, the most up-to-date version of the respective instructions for use, which is available for download on the website <https://www.ohst.de/fu-instructions/>, must be observed. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads, the most up-to-date version of the respective instructions for use, which is available for download on the website <https://www.ohst.de/fu-instructions/>, must also be observed. After reduction, the definitive stability, mobility and muscle tension must be checked.

**Caution:** Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

**Caution:** If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

**Caution:** The patient weight limitation stated in the contraindications must be observed for the Expersus femoral stem sizes listed under 1.1.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

The porous coating (TPS, BONIT®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

**Caution:** Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free gloves.

For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

**Caution:** when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

If femoral heads with a neck (XL and XXL) are used, the range of motion is reduced by about 30° and achieves values between 80° and 100° in flexion and extension.

**Caution:** Component malposition or the use of a constrained hip or a total hip with a modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement, premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

## **2.4** Surgical Technique

After opening the joint capsule and dislocating the femoral head from the cup, the femoral head should be resected as determined during the preoperative planning and completely removed.

A femoral neck resection guide is optionally available for marking the resection plane. Then open the medullary cavity using the cavity chisel. Apply the chisel to the far lateral and dorsal side in order to facilitate the subsequent driving in of the rasps in the direction of the femoral axis. The cavity chisel should be positioned relative to the desired antetorsion. Avoid fracturing the greater trochanter. The medullary cavity opening can be enlarged with the surgical awl.

Once the medullary cavity has been opened, rasp the femoral stem. The first rasp determines the alignment of the subsequent sizes. For this reason, take care to ensure that the antetorsion is already correct during the first rasping step. The antetorsion can be checked using the guide rod and rasp handle and is usually 10°-15°. Start the procedure with the smallest rasp size, which is connected to the rasp handle. After that, use the subsequent rasps in ascending order until the preoperatively determined size is reached. The sizes of the rasps match the sizes of the implants. Verify the correct position of the rasp inside the femur with an image converter. If during the rasping procedure the stem size does not correspond to the preoperatively determined size and there is a difference of two or more sizes, the alignment of the axis may be incorrect or there may be an osseous obstruction present. In this case, the selected stem might be too small and would therefore not provide the required stability. On request, special rasp handles are available for the various access routes. The rasps are designed in such a way that they can also be used for trial reductions. The leg length, the range of motion and the tension of the ligaments can be checked by attaching the various trial cones and trial heads. Remove the handle from the rasp, which remains inside the femur, and attach the trial cone to the rasp. Trial cones with two offsets are offered for trial reductions of the various stem versions. The trial cones snap into place as soon as they are in their correct position. Then attach the trial head to the trial cone by hand. Trial heads with various diameters are available for trial reductions in the neck lengths S to XL. After reduction, the definitive stability, mobility and muscle tension should be checked.

With the cemented Expersus femoral stem, the cement restrictor is implanted first and the cement inserted afterwards; this surgical step is not applicable to the cementless variant of the Expersus femoral stem.

The cement restrictor is used depending on the preoperatively determined diaphyseal diameter of the medullary cavity. A special insertion instrument is available for controlled insertion to the optimal depth. The size 1 cement restrictor is recommended for medullary cavity diameters between 14 mm and 19 mm, or size 2 for diameters between 18 mm and 22 mm. Screw the cement restrictor onto the threaded rod of the insertion instrument and insert it into the medullary cavity. The insertion depth can be read off the scale on the threaded rod and should exceed the stem length of the prosthetic femoral stem to be implanted by at least 5 mm, but never more than 20 mm. It is also possible to gauge the position of the cement restrictor inside the femoral stem by means of X-ray control based on the position of the radiopaque ring. Thereafter, remove the insertion instrument by rotating it anti-clockwise. After inserting the cement restrictor, rinse and dry the implant bed. The cement restrictor should be implanted before the bone cement is mixed. The cement is to be inserted

using a modern cementing technique; in doing so, the instructions for use of the cement manufacturer must be observed.

Insert the prosthetic device into the implant bed up to the resection plane using the stem impactor. When using the cemented variant, complete this step immediately after injecting the bone cement and gently apply pressure while waiting for the bone cement to harden. Then thoroughly clean and dry the prosthetic cone and attach the femoral head with the previously determined neck length (S to XL).

**Caution:** Do not use femoral heads with a neck length greater than LXI

Then reduce the stem with the femoral head into the cup and check the range of motion, leg length and tension of the ligaments. The operation is routinely completed with layer-by-layer wound closure.

### **3. Packaging and Sterility**

Depending on the sterilisation method used, implants are packaged in a double or triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

**Caution:** the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

### **4. Preoperative Planning and Postoperative Care**

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size. Surgery may only be performed once it has been established that the patient is able to tolerate the implant material. In addition, the patient weight limitation stated in the contraindications must be observed for the Expersus femoral stem sizes named under 1.1. The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format. To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant. Recognised procedures must be used for postoperative care.

### **5. Indications**

- Advanced degeneration of the hip joint due to degenerative, post-traumatic rheumatoid arthritis or congenital hip dysplasia
- Fracture (e.g. femoral neck fracture) or avascular necrosis of the femoral head

- Sequelae of earlier surgical procedures, e.g. internal fixation, articular reconstruction, arthrodesis
- Certain cases of ankylosis

Total hip replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total hip replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

## 6. Contraindications

- Acute or chronic local or systemic infection
- Severe muscle, nerve or vascular diseases that endanger the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Any underlying condition that might compromise the function of the implant
- Hypersensitivity to the materials used
- Patients weighing more than 65 kg for the cemented variant of the Expersus femoral stem size 1 STD 135°, size 2 STD 135° and size 2 STD 125° (REF 367-1556, 367-1557 and 367-1547)
- Patients weighing more than 90 kg for the cemented variant of the Expersus femoral stem size 3 STD 135° and size 3 STD 125° (REF 367-1558 and 367-1548)

## 7. Risk factors and conditions that may affect the success of the surgery

**Caution:** Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Use of high-dose cortisone or cytotoxic drugs
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Muscle spasms or other spastic conditions
- Growth in children and adolescents

- Expected extreme stresses, e.g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities that complicate fixation of the implant
- Weakening of the supporting structures due to a tumour

## **8. Undesirable Effects**

The undesirable effects and risks listed below can occur in the course of a total hip arthroplasty (THA) or hemiarthroplasty:

- Foreign body reactions, osteolysis, loosening
- ARMD / trunnionosis / metallosis
- Toxic reactions
- Luxation / dislocation / dissociation
- Limited ROM
- Loosening
- Implant failure
- Difference in leg length
- Stem subsidence
- Bone fracture
- Instability
- Tissue damage
- Heterotopic ossification
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack / stroke
- Bone cement reaction / bone cement implantation syndrome (BCIS) (e.g. cardiac arrhythmia, increased pulmonary vascular resistance)

Caution: Revision surgery may be required if specific undesirable effects occur.

## **9. Patient Information and Documentation**

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. Based on the current state of technology, the expected service life is approx. 10 to 15 years.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant. Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks. The user can download MRI safety information from the website <https://www.ohst.de/professionals/>. The patient can download the patient information from the website <https://www.ohst.de/patient-information/>. A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

## **10.** Implant Card

Following surgery, the patient must be given an implant card containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant card must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website. The patient's details (patient name or patient ID) must be documented in the implant card, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the card. The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

## **11.** Explanation of the Label Symbols

The user can download the list of symbols from the website <https://www.ohst.de/professionals/>.

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