OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 3385 5420-0

Fax: +49 3385 5420-99

www.ohst.de

CE 0482

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IMPLANT

Flat-head screw

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for hip arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website <https://www.ohst.de/ifu-instructions/>.

1. Product Description and Implant Materials

The flat-head screws are used for optional fixation of acetabular prostheses in the context of total hip arthroplasty.

The flat-head screws are made of a Ti6Al4V alloy (ISO 5832-3). They have a head diameter of 8 mm and a cancellous bone screw thread with \varnothing 6.5 mm. The flat-head screws come in lengths of 15 mm to 60 mm with 5-mm increments.

The flat-head screws may be optionally used only in patients who require hip arthroplasty and who fulfil the relevant indications, taking into account the contraindications of the entire prosthetic system selected. In addition, the optional use of flat head screws is possible only to treat skeletally mature patients.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

1.1 Overview of Implants

Designation	Material	Reference number
Flat-head Screw \varnothing 6.5 x 15 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-15
Flat-head Screw \varnothing 6.5 x 20 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-20
Flat-head Screw \varnothing 6.5 x 25 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-25
Flat-head Screw \varnothing 6.5 x 30 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-30
Flat-head Screw \varnothing 6.5 x 35 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-35
Flat-head Screw \varnothing 6.5 x 40 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-40
Flat-head Screw \varnothing 6.5 x 45 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-45
Flat-head Screw \varnothing 6.5 x 50 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-50
Flat-head Screw \varnothing 6.5 x 55 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-55
Flat-head Screw \varnothing 6.5 x 60 mm, self-tapping	ISO 5832-3 Ti6Al4V	000-290-60

1.2 Overview of accessories

1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
Hexagonal Screwdriver SW 3.5mm grey with silicone handle, L=250 mm	250161-5
Cardan Screwdriver, SW 3.5mm with silicone handle grey, L= 273mm	00-092-10

1.2.2 Accessories

Designation	Reference number
Screw gauge curved L=232mm, angle: 135°	367-115
Screw Holding Forcep	367-1021

1.3 Accompanying reference documents

Designation	Reference number
Primaro Cup surgical technique	50000349
Implant Card	50000572
Patient information Hip joint replacement	50000841
MRI safety information	50000851
List of symbols	50000859

2. Handling

2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>, must be observed.

Caution: Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

Caution: the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants

is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explanation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint.

The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

2.2 Authorised Component Combinations

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

2.3 Instructions for Use

Caution: The application instructions of the implant system being used with the flat-head screws have to be observed.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

Caution: when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can

cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

Caution: Component malposition or the use of a constrained hip or a total hip with a modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement, premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

2.4 Surgical Technique

When using self-tapping flat-head screws, the bone has to be pre-drilled with a drill of diameter 3.2 mm according to the screw length being used. The depth of the drill-hole can be verified with the screw gauge. After that, the flat-head screws are inserted according to the surgical technique of the implant system being used.

Caution: The position of the screw holes of the implant/screw combination to be used, as well as the corresponding screw lengths, must be selected in a way that precludes damage to vessels and nerves or unintentional perforation of structures.

Caution: When tightening the screws, ensure that the alignment of the cup component does not change.

Caution: The screws must only be tightened with the appropriate screwdriver.

3. Packaging and Sterility

Depending on the sterilisation method used, implants are packaged in a double or triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative Planning and Postoperative Care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant. Recognised procedures must be used for postoperative care.

5. Indications

- Optional fixation of acetabular components in the context of primary hip arthroplasty or revision.

In addition, the indications and recommendations for the prosthetic system being used, of which flat-head screws are integral components, are to be taken into account.

6. Contraindications

- Allergy to the materials used
- Bone substance unsuitable for ensuring that the screw holds

In addition, the contraindications and recommendations for the prosthetic system being used, of which flat-head screws are integral components, are to be taken into account.

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

- Overweight
- Smoking
- Diabetes mellitus
- Mental illness
- Anaemia
- Intra-articular corticosteroid injections in the affected joint, less than three months prior to the operation
- Generally increased risk factors for an operation

8. Undesirable Effects

The undesirable effects listed below can occur when using the flat-head screws:

- Neurovascular damage
- Tissue damage
- Deep infections

- Screw breakage, screw bending

In addition, the following undesirable effects of the implant system using flat-head screws can occur in total hip arthroplasty (THA):

- Noise generated due to a ceramic on ceramic articulation
- Foreign body reactions, osteolysis, loosening
- Toxic reactions
- Sensitisation
- Limited ROM
- Instability
- Luxation, dislocation, dissociation
- Implant failure
- Bone fracture
- Loosening
- Drifting/ tilting
- Heterotopic ossification
- Difference in leg length
- Tissue damage
- Iliopsoas syndrome / irritation
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack / stroke

Caution: Several of these undesirable effects are accompanied by pain. Revision surgery may be required if these undesirable effects occur.

9. Patient Information and Documentation

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. All information provided to the patient must be documented in writing by the surgeon performing the surgery. Based on the current state of technology, the expected service life is approx. 15 years.

Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant.

Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

The user can download MRI safety information from the website <https://www.ohst.de/professionals/>. The patient can download the patient information from the website <https://www.ohst.de/patient-information/>. A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

10. Implant Card

Following surgery, the patient must be given an implant card containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant card must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website.

The patient's details (patient name or patient ID) must be documented in the implant card, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the card.

The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

11. Explanation of the Label Symbols

The user can download the list of symbols from the website <https://www.ohst.de/professionals/>.

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