OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 3385 5420-0

Fax: +49 3385 5420-99

www.ohst.de

CE 0482

Müller Geradschaft		Müller право стебло	
Deutsch (de)	3	Български (bg)	82
Müller Straight Stem		Пряма ніжка ендопротеза Müller	
English (en)	13	Українська (uk)	92
Tige droite de Müller		„Müller“ tiesus stiebas	
Français (fr)	22	Lietuvių k. (lt)	102
Stelo dritto Müller		Прямая ножка ендопротеза Müller	
Italiano (it)	32	Русский (ru)	111
Vástago recto Müller		الجدع المستقيم Müller	
Español (es)	42	العربية (ar)	123
Haste reta Müller		Müller endoprotēzes kājiņa	
Português (pt)	52	Latviešu (lv)	131
Ευθύγραμμο στέλεχος Müller		Прав стем Müller	
Ελληνικά (el)	62	Македонски (mk)	140
Müller rovny dřík			
Česky (cs)	73		

IMPLANTAT

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IMPLANT

Müller Straight Stem

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for hip arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website <https://www.ohst.de/ifu-instructions/>.

1. Product Description and Implant Materials

The Müller Straight Stem is a femoral stem prosthesis and can be used with suitable femoral heads for total hip arthroplasty or hemiarthroplasty.

The Müller Straight Stem is produced from CoCrMo cast alloy (ISO 5832-4) and fitted with a 12/14 cone, and can be combined with both metal and ceramic femoral heads. The standard and lateral variants are each available in 11 sizes. The stems differ from one another in terms of length and offset. All Müller Straight Stems have the same CCD angle of 135°.

The Müller Straight Stem is inserted with dual fixation (dorsoventral cementing, mediolateral clamping). The Müller Straight Stem may only be used in patients who require total hip arthroplasty or hip hemiarthroplasty and who fulfil the relevant indications, taking into account the contraindications. In addition, the Müller Straight Stem may only be used to treat skeletally mature patients.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

1.1 Overview of Implants

Designation	Material	Reference number
Müller Straight Stem size 6.25 Standard	ISO 5832-4 CoCrMo	110-065
Müller Straight Stem size 7.5 Standard	ISO 5832-4 CoCrMo	110-075
Müller Straight Stem size 8.75 Standard	ISO 5832-4 CoCrMo	110-085
Müller Straight Stem size 10 Standard	ISO 5832-4 CoCrMo	110-100
Müller Straight Stem size 11.25 Standard	ISO 5832-4 CoCrMo	110-115
Müller Straight Stem size 12.5 Standard	ISO 5832-4 CoCrMo	110-125
Müller Straight Stem size 13.75 Standard	ISO 5832-4 CoCrMo	110-135
Müller Straight Stem size 15 Standard	ISO 5832-4 CoCrMo	110-150
Müller Straight Stem size 16.25 Standard	ISO 5832-4 CoCrMo	110-165
Müller Straight Stem size 17.5 Standard	ISO 5832-4 CoCrMo	110-175

Designation	Material	Reference number
Müller Straight Stem size 20 Standard	ISO 5832-4 CoCrMo	110-200
Müller Straight Stem size 6.25 Lateral (*)	ISO 5832-4 CoCrMo	111-065
Müller Straight Stem size 7.5 Lateral	ISO 5832-4 CoCrMo	111-075
Müller Straight Stem size 8.75 Lateral	ISO 5832-4 CoCrMo	111-085
Müller Straight Stem size 10 Lateral	ISO 5832-4 CoCrMo	111-100
Müller Straight Stem size 11.25 Lateral	ISO 5832-4 CoCrMo	111-115
Müller Straight Stem size 12.5 Lateral	ISO 5832-4 CoCrMo	111-125
Müller Straight Stem size 13.75 Lateral	ISO 5832-4 CoCrMo	111-135
Müller Straight Stem size 15 Lateral	ISO 5832-4 CoCrMo	111-150
Müller Straight Stem size 16.25 Lateral	ISO 5832-4 CoCrMo	111-165
Müller Straight Stem size 17.5 Lateral	ISO 5832-4 CoCrMo	111-175
Müller Straight Stem size 20 Lateral	ISO 5832-4 CoCrMo	111-200

(*) Weight limit see Contraindications

1.2 Overview of accessories

1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
Instrument set Müller Straight Stem basic tray	506-2109
Instrument set Müller Straight Stem rasps	506-2110
Instrument set Müller Straight Stem starter set	506-2111
Instrument set Müller Straight Stem optional	506-2112
Femoral Stem Extraction Instrument Set	206-010

1.2.2 Accessories

Designation	Reference number
X-ray Templates Müller Straight Stem Standard CoCrMo, Dual Fixation	5500-10021
X-ray Templates Müller Straight Stem Lateral CoCrMo, Dual Fixation	5500-20021
Cement restrictor size 1, external Ø24mm	506-101
Cement restrictor size 2, external Ø27mm	506-102

1.3 Accompanying reference documents

Designation	Reference number
Surgical Technique for Müller Straight Stem CoCrMo Hook Coupling Connector	50000017
Instructions for use Extraction Instrument Set	50000428
Instructions for use Stem Repositioner	50000427
Implant Card	50000572

Designation	Reference number
Patient information Hip joint replacement	50000841
MRI safety information	50000851
List of symbols	50000859

2. Handling

2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>, must be observed.

Caution: Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

Caution: the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint. The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

2.2 Authorised Component Combinations

The following femoral heads must be used for assembly with the implant:

Nominal Ø	Name	Neck length	Reference number
Ø22 mm	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L	-4 / 0 / +4 mm	030-2200 to 030-2202
Ø28 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL / XXL	-3,5 / 0 / +3,5 / +7 / +10,5 mm	020-2800 to 020-2804
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL / XXL	-3,5 / 0 / +3,5 / +7 / +10,5 mm	030-2800 to 030-2804
	BILOX®delta Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-1140 to 367-1142
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	013-001 to 013-003
Ø32 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL / XXL	-4 / 0 / +4 / +8 / +12 mm	020-3200 to 020-3204
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL / XXL	-4 / 0 / +4 / +8 / +12 mm	030-3200 to 030-3204
	BILOX®delta Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	367-1143 to 367-1145, 367-1149
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	013-004 to 013-007
Ø36 mm	BILOX®delta Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	367-1146 to 367-1148, 367-1150
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	013-008 to 013-011
Ø40 mm to Ø60 mm	Fracture head size S implant steel ISO 5832-9, 12/14	-4 mm	155-140 to 155-160
	Fracture head size M implant steel ISO 5832-9, 12/14	0 mm	155-040 to 155-060
	Fracture head size L implant steel ISO 5832-9, 12/14	+4 mm	155-240 to 155-260

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

2.3 Instructions for Use

The application of the implant is done with cement. A cement restrictor must be used for cemented applications. Before using the cement restrictor, the most up-to-date version of the respective instructions for use, which is available for download on the website <https://www.ohst.de/fu-instructions/>, must be observed. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads, the most up-to-date version of the respective instructions for use, which is available for download on the website <https://www.ohst.de/fu-instructions/>, must also be observed. After reduction, the definitive stability, mobility and muscle tension must be checked.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

Caution: when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

If femoral heads with a neck (XL and XXL) are used, the range of motion is reduced by about 30° and achieves values between 80° and 100° in flexion and extension.

Caution: Component malposition or the use of a constrained hip or a total hip with a modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement, premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

2.4 Surgical Technique

The Müller Straight Stem is used in the usual positionings and approaches for cemented femoral stem endoprostheses. Its implantation consists of the following steps.

Resection of the femoral neck: Osteotomy with the oscillating saw is done according to the pre-operative planning.

Preparation of the acetabulum and implantation of the acetabular cup.

Opening of the proximal femur: A box chisel is provided for opening the proximal femur. The chisel should be applied to the far lateral and dorsal side in order to facilitate the subsequent driving in of the rasps.

Preparation of the implant bed with the rasps: Gradual rasping of the femoral prosthesis stem bed is carried out, starting with the smallest rasp. Special attention must be paid to the correct, central positioning, antetorsion and the relationship to the planned centre of rotation. The stem size to be implanted corresponds to the last rasp size used. If the planned rasp size or depth is not achieved, the femur rasp can be used to create space proximally.

Trial reduction: Using rasps with hook coupling, the trial reduction is performed by fitting standard or lateralised trial heads directly onto the hook coupling of the rasp. Trial heads with various head diameters are available in the neck lengths S to XXL. **Preparation of cementation:** Based on the measured diameter of the medullary cavity usually a blocker made of PE (cement restrictor) is inserted. A special inserter allows the controlled insertion of the blocker at optimal depth. The cement restrictor should be implanted before the bone cement is mixed. Cement restrictor size 1 is recommended for medullary cavity diameters between 14 mm and 19 mm, or size 2 for diameters between 18 mm and 22 mm. The cement restrictor should be positioned 0.5 to 2 cm distal to the tip of the prosthesis. During an X-ray examination, the position of the cement restrictor can be assessed from the position of the radiopaque ring. After insertion of the restrictor, cleaning lavage and subsequent drying of the prosthesis stem bed are performed.

Introduction of the definitive shaft: The cement must be inserted according to the modern cementing technique and the instructions for use of the cement manufacturer must be observed. It must be ensured that the viscosity of the cement allows a deep penetration into the cancellous bone. This can be supported additionally by sustained impression (pressurisation). The Müller Straight Stem is inserted into the prepared prosthesis bed with the stem impactor. Check for the correct axial alignment as well as the insertion depth and the antetorsion while doing this. The final blocking can be supported with the use of uniform light hammer blows.

Reduction: After careful removal of unnecessary cement residue, cleaning lavage and further, successful trial reduction, the cone of the implanted Müller Straight Stem is to be thoroughly cleaned and dried. Before the final reduction, an original 12/14 femoral head will be attached to the stem taper and will be fixed. It should be noted that no heads with a neck length greater than XXL (+ 12 mm) may be used! For further information please see Operation Techniques - Müller Straight Stem (cf. 1.3 Accompanying reference documents).

3. Packaging and Sterility

Depending on the sterilisation method used, implants are packaged in a double or triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/fifu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative Planning and Postoperative Care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognised procedures must be used for postoperative care.

5. Indications

- Advanced wear of the hip joint due to degenerative, post-traumatic or rheumatoid arthritis
- Fracture or avascular necrosis of the femoral head
- State following previous operations, e. g. internal fixation, joint reconstruction, arthrodesis

Total hip replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total hip replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

6. Contraindications

- Acute or chronic, local or systemic infection
- Severe muscle, nerve or vascular diseases, which endanger the affected limb

- Lack of bone substance or defective bone quality, which endangers stable sit of the prosthesis
- Any concomitant disease that can endanger the function of the implant
- Hypersensitivity to the materials used
- Patients weighing more than 75 kg for the Müller Straight Stem size 6.25 LAT (111-065)

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

- Overweight
- Smoking
- Diabetes mellitus
- Mental illness
- Anaemia
- Intra-articular corticosteroid injections in the affected joint, less than three months prior to the operation
- Generally increased risk factors for an operation

8. Undesirable Effects

The undesirable effects and risks listed below can occur in the course of a total hip arthroplasty (THA) or hemiarthroplasty:

- Foreign body reactions, osteolysis, loosening
- Toxic reactions
- Luxation / dislocation / dissociation
- Implant failure
- Limited ROM
- Difference in leg length
- Instability
- ARMD / trunnionosis / metallosis
- Tissue damage
- Heterotopic ossification
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack / stroke
- Bone cement reaction / bone cement implantation syndrome (BCIS) (e.g. cardiac arrhythmia, increased pulmonary vascular resistance)

Caution: Several of these undesirable effects are accompanied by pain. Revision surgery may be required if these undesirable effects occur.

9. Patient Information and Documentation

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. Based on the current state of technology, the expected service life is approx. 15 years.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, healing up of the implant, induction of electrical currents, loosening of the implant.

Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks. The user can download MRI safety information from the website <https://www.ohst.de/professionals/>. The patient can download the patient information from the website <https://www.ohst.de/patient-information/>. A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

10. Implant Card

Following surgery, the patient must be given an implant card containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant card must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website. The patient's details (patient name or patient ID) must be documented in the implant card, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the card.

The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

11. Explanation of the Label Symbols

The user can download the list of symbols from the website <https://www.ohst.de/professionals/>.

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