**OHST**   
medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 3385 5420-0

Fax: +49 3385 5420-99

www.ohst.de

**CE** 0482

---

|                                      |     |
|--------------------------------------|-----|
| ZEN Kniesystem                       |     |
| Deutsch (de)                         | 3   |
| ZEN Knee System                      |     |
| English (en)                         | 23  |
| ZEN Système de genou                 |     |
| Français (fr)                        | 43  |
| ZEN Sistema per ginocchio            |     |
| Italiano (it)                        | 63  |
| Sistema de rodilla ZEN               |     |
| Español (es)                         | 83  |
| <b>Σύστημα γονάτου ZEN</b>           |     |
| Ελληνικά (el)                        | 104 |
| <b>Система колінного суглоба ZEN</b> |     |
| Українська (uk)                      | 125 |
| نظام الركبة ZEN                      |     |
| العربية (ar)                         | 147 |

## IMPLANTAT

## ZEN Kniesystem

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## IMPLANT

### ZEN Knee System

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for knee arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website <https://www.ohst.de/ifu-instructions/>.

#### 1. Product Description and Implant Materials

The ZEN Knee System is a tricompartmental uncoupled system for total knee replacement, for cemented, cementless or hybrid use, to alleviate pain and improve the functioning of the knee joint in case of severe joint damage.

If secure fixation of the standard components is no longer possible due to insufficient bone substance of the proximal tibia, the ZEN Tibia Modular should be used.

To reduce the allergy potential, there are components with TiNbN coating.

The ZEN Knee System is available as a highly congruent mobile-bearing variant and as a fixed-bearing variant. All components are available in various sizes in order to adapt to the individual conditions of each patient.

The ZEN knee prostheses are implants that remain in the body permanently and either partially or fully replace the damaged knee. The knee joint in its entirety is replaced by three components: the femoral component, the tibial component and the patellar component.

Irrespective of the shape and fixation manner of the individual joint components, they always perform the same functions (ISO 7207-1, ISO 21536):

- |                     |                                                                                                                                                                                                                                      |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Femoral component:  | the component of a total knee joint replacement intended to be secured to the femur to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.    |
| Tibial component:   | the component of a total knee joint replacement intended to be secured to the tibia to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.    |
| Patellar component: | the component of a total or partial knee joint replacement which is used to replace the articulating surface of the patella. These implants can be manufactured as one component or a set of components to be assembled by the user. |

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

## 1.1 Overview of Implants

### 1.1.1 Femoral components

Femur STD:

| Size | Version | Reference number  |                            |                                  |
|------|---------|-------------------|----------------------------|----------------------------------|
|      |         | cemented          |                            | cementless                       |
|      |         | ISO 5832-4 CoCrMo | ISO 5832-4 CoCrMo<br>TiNbN | ISO 5832-4 CoCrMo<br>TPS, BONIT® |
| 1    | right   | 482-002           | 482-122                    | 482-014                          |
| 2    | right   | 482-004           | 482-124                    | 482-016                          |
| 3    | right   | 482-006           | 482-126                    | 482-018                          |
| 4    | right   | 482-008           | 482-128                    | 482-020                          |
| 5    | right   | 482-010           | 482-130                    | 482-022                          |
| 6    | right   | 482-012           | 482-132                    | 482-024                          |
| 1    | left    | 482-001           | 482-121                    | 482-013                          |
| 2    | left    | 482-003           | 482-123                    | 482-015                          |
| 3    | left    | 482-005           | 482-125                    | 482-017                          |
| 4    | left    | 482-007           | 482-127                    | 482-019                          |
| 5    | left    | 482-009           | 482-129                    | 482-021                          |
| 6    | left    | 482-011           | 482-131                    | 482-023                          |
| 2,5  | right   | 482-212           | 482-220                    | 482-228                          |
| 3,5  | right   | 482-214           | 482-222                    | 482-230                          |
| 4,5  | right   | 482-216           | 482-224                    | 482-232                          |
| 5,5  | right   | 482-218           | 482-226                    | 482-234                          |
| 2,5  | left    | 482-211           | 482-219                    | 482-227                          |
| 3,5  | left    | 482-213           | 482-221                    | 482-229                          |
| 4,5  | left    | 482-215           | 482-223                    | 482-231                          |
| 5,5  | left    | 482-217           | 482-225                    | 482-233                          |

### 1.1.2 Tibial components

Tibia STD, mobile bearing:

| Size | Version | Reference number  |                            |                                  |
|------|---------|-------------------|----------------------------|----------------------------------|
|      |         | cemented          |                            | cementless                       |
|      |         | ISO 5832-4 CoCrMo | ISO 5832-4 CoCrMo<br>TiNbN | ISO 5832-4 CoCrMo<br>TPS, BONIT® |
| 1    | right   | 482-098           | 482-134                    | 482-026                          |
| 2    | right   | 482-100           | 482-136                    | 482-028                          |
| 3    | right   | 482-102           | 482-138                    | 482-030                          |

| Size | Version | Reference number  |                            |                                  |
|------|---------|-------------------|----------------------------|----------------------------------|
|      |         | cemented          |                            | cementless                       |
|      |         | ISO 5832-4 CoCrMo | ISO 5832-4 CoCrMo<br>TiNbN | ISO 5832-4 CoCrMo<br>TPS, BONIT® |
| 4    | right   | 482-104           | 482-140                    | 482-032                          |
| 5    | right   | 482-106           | 482-142                    | 482-034                          |
| 6    | right   | 482-108           | 482-144                    | 482-036                          |
| 1    | left    | 482-097           | 482-133                    | 482-025                          |
| 2    | left    | 482-099           | 482-135                    | 482-027                          |
| 3    | left    | 482-101           | 482-137                    | 482-029                          |
| 4    | left    | 482-103           | 482-139                    | 482-031                          |
| 5    | left    | 482-105           | 482-141                    | 482-033                          |
| 6    | left    | 482-107           | 482-143                    | 482-035                          |
| 2,5  | right   | 482-236           | 482-244                    | 482-252                          |
| 3,5  | right   | 482-238           | 482-246                    | 482-254                          |
| 4,5  | right   | 482-240           | 482-248                    | 482-256                          |
| 5,5  | right   | 482-242           | 482-250                    | 482-258                          |
| 2,5  | left    | 482-235           | 482-243                    | 482-251                          |
| 3,5  | left    | 482-237           | 482-245                    | 482-253                          |
| 4,5  | left    | 482-239           | 482-247                    | 482-255                          |
| 5,5  | left    | 482-241           | 482-249                    | 482-257                          |

Tibia STD, fixed bearing:

| Size | Version | Reference number  |                            |                                  |
|------|---------|-------------------|----------------------------|----------------------------------|
|      |         | cemented          |                            | cementless                       |
|      |         | ISO 5832-4 CoCrMo | ISO 5832-4 CoCrMo<br>TiNbN | ISO 5832-4 CoCrMo<br>TPS, BONIT® |
| 1    | right   | 522-902           | 522-926                    | 522-914                          |
| 2    | right   | 522-904           | 522-928                    | 522-916                          |
| 3    | right   | 522-906           | 522-930                    | 522-918                          |
| 4    | right   | 522-908           | 522-932                    | 522-920                          |
| 5    | right   | 522-910           | 522-934                    | 522-922                          |
| 6    | right   | 522-912           | 522-936                    | 522-924                          |
| 1    | left    | 522-901           | 522-925                    | 522-913                          |
| 2    | left    | 522-903           | 522-927                    | 522-915                          |
| 3    | left    | 522-905           | 522-929                    | 522-917                          |
| 4    | left    | 522-907           | 522-931                    | 522-919                          |
| 5    | left    | 522-909           | 522-933                    | 522-921                          |
| 6    | left    | 522-911           | 522-935                    | 522-923                          |
| 2,5  | right   | 522-702           | 522-710                    | 522-718                          |
| 3,5  | right   | 522-704           | 522-712                    | 522-720                          |
| 4,5  | right   | 522-706           | 522-714                    | 522-722                          |

| Size | Version | Reference number  |                            |                                  |
|------|---------|-------------------|----------------------------|----------------------------------|
|      |         | cemented          |                            | cementless                       |
|      |         | ISO 5832-4 CoCrMo | ISO 5832-4 CoCrMo<br>TiNbN | ISO 5832-4 CoCrMo<br>TPS, BONIT® |
| 5,5  | right   | 522-708           | 522-716                    | 522-724                          |
| 2,5  | left    | 522-701           | 522-709                    | 522-717                          |
| 3,5  | left    | 522-703           | 522-711                    | 522-719                          |
| 4,5  | left    | 522-705           | 522-713                    | 522-721                          |
| 5,5  | left    | 522-707           | 522-715                    | 522-723                          |

Tibia Modular, mobile bearing:

| Size | Reference number                        |                                                  |
|------|-----------------------------------------|--------------------------------------------------|
|      | cemented                                |                                                  |
|      | ISO 5832-4 CoCrMo<br>ISO 5832-3 Ti6Al4V | ISO 5832-4 CoCrMo<br>ISO 5832-3 Ti6Al4V<br>TiNbN |
| 2    | 482-170                                 | 482-202                                          |
| 3    | 482-171                                 | 482-203                                          |
| 4    | 482-172                                 | 482-204                                          |
| 5    | 482-173                                 | 482-205                                          |
| 6    | 482-174                                 | 482-206                                          |

Tibia Modular, fixed bearing:

| Size | Reference number                        |                                                  |
|------|-----------------------------------------|--------------------------------------------------|
|      | cemented                                |                                                  |
|      | ISO 5832-4 CoCrMo<br>ISO 5832-3 Ti6Al4V | ISO 5832-4 CoCrMo<br>ISO 5832-3 Ti6Al4V<br>TiNbN |
| 2    | 522-937                                 | 522-942                                          |
| 3    | 522-938                                 | 522-943                                          |
| 4    | 522-939                                 | 522-944                                          |
| 5    | 522-940                                 | 522-945                                          |
| 6    | 522-941                                 | 522-946                                          |

Augment:

| Size | Height<br>[mm] | Version                     | Reference number   |
|------|----------------|-----------------------------|--------------------|
|      |                |                             | cemented           |
|      |                |                             | ISO 5832-3 Ti6Al4V |
| 2    | 5              | left lateral / right medial | 482-182            |
| 3    | 5              | left lateral / right medial | 482-183            |
| 4    | 5              | left lateral / right medial | 482-184            |
| 5    | 5              | left lateral / right medial | 482-185            |
| 6    | 5              | left lateral / right medial | 482-186            |
| 2    | 10             | left lateral / right medial | 482-187            |
| 3    | 10             | left lateral / right medial | 482-188            |
| 4    | 10             | left lateral / right medial | 482-189            |
| 5    | 10             | left lateral / right medial | 482-190            |
| 6    | 10             | left lateral / right medial | 482-191            |
| 2    | 5              | left medial / right lateral | 482-192            |
| 3    | 5              | left medial / right lateral | 482-193            |
| 4    | 5              | left medial / right lateral | 482-194            |
| 5    | 5              | left medial / right lateral | 482-195            |
| 6    | 5              | left medial / right lateral | 482-196            |
| 2    | 10             | left medial / right lateral | 482-197            |
| 3    | 10             | left medial / right lateral | 482-198            |
| 4    | 10             | left medial / right lateral | 482-199            |
| 5    | 10             | left medial / right lateral | 482-200            |
| 6    | 10             | left medial / right lateral | 482-201            |

Stem extension:

| Length | Version     | Reference number   |
|--------|-------------|--------------------|
|        |             | cemented           |
|        |             | ISO 5832-3 Ti6Al4V |
| 40     | straight    | 482-175            |
| 75     | straight    | 482-176            |
| 110    | straight    | 482-177            |
| 110    | with offset | 482-178            |

Centraliser:

| Size | Reference number  |
|------|-------------------|
|      | cemented          |
|      | ISO 5834-2 UHMWPE |
| Ø10  | 482-179           |
| Ø12  | 482-180           |
| Ø14  | 482-181           |

Insert STD, mobile bearing:

| Size | Height<br>[mm] | Reference number                             |                                            | Size | Height<br>[mm] | Reference number                             |                                            |
|------|----------------|----------------------------------------------|--------------------------------------------|------|----------------|----------------------------------------------|--------------------------------------------|
|      |                | ISO 5834-2<br>UHMWPE<br>ISO 5832-1<br>1.4441 | Vitamin-E<br>XL-PE<br>ISO 5832-1<br>1.4441 |      |                | ISO 5834-2<br>UHMWPE<br>ISO 5832-1<br>1.4441 | Vitamin-E<br>XL-PE<br>ISO 5832-1<br>1.4441 |
| 1    | 10             | 482-054                                      | 482-301                                    | 4    | 10             | 482-078                                      | 482-310                                    |
| 1    | 12             | 482-055                                      | 482-302                                    | 4    | 12             | 482-079                                      | 482-311                                    |
| 1    | 14             | 482-056                                      | 482-303                                    | 4    | 14             | 482-080                                      | 482-312                                    |
| 2    | 10             | 482-062                                      | 482-304                                    | 5    | 10             | 482-086                                      | 482-313                                    |
| 2    | 12             | 482-063                                      | 482-305                                    | 5    | 12             | 482-087                                      | 482-314                                    |
| 2    | 14             | 482-064                                      | 482-306                                    | 5    | 14             | 482-088                                      | 482-315                                    |
| 3    | 10             | 482-070                                      | 482-307                                    | 6    | 10             | 482-094                                      | 482-316                                    |
| 3    | 12             | 482-071                                      | 482-308                                    | 6    | 12             | 482-095                                      | 482-317                                    |
| 3    | 14             | 482-072                                      | 482-309                                    | 6    | 14             | 482-096                                      | 482-318                                    |

Insert UC, mobile bearing:

| Size | Height<br>[mm] | Reference number                             |                                            | Size | Height<br>[mm] | Reference number                             |                                            |
|------|----------------|----------------------------------------------|--------------------------------------------|------|----------------|----------------------------------------------|--------------------------------------------|
|      |                | ISO 5834-2<br>UHMWPE<br>ISO 5832-1<br>1.4441 | Vitamin-E<br>XL-PE<br>ISO 5832-1<br>1.4441 |      |                | ISO 5834-2<br>UHMWPE<br>ISO 5832-1<br>1.4441 | Vitamin-E<br>XL-PE<br>ISO 5832-1<br>1.4441 |
| 1    | 10             | 482-049                                      | 482-319                                    | 4    | 10             | 482-073                                      | 482-334                                    |
| 1    | 12             | 482-050                                      | 482-320                                    | 4    | 12             | 482-074                                      | 482-335                                    |
| 1    | 14             | 482-051                                      | 482-321                                    | 4    | 14             | 482-075                                      | 482-336                                    |
| 1    | 17             | 482-052                                      | 482-322                                    | 4    | 17             | 482-076                                      | 482-337                                    |
| 1    | 20             | 482-053                                      | 482-323                                    | 4    | 20             | 482-077                                      | 482-338                                    |
| 2    | 10             | 482-057                                      | 482-324                                    | 5    | 10             | 482-081                                      | 482-339                                    |
| 2    | 12             | 482-058                                      | 482-325                                    | 5    | 12             | 482-082                                      | 482-340                                    |
| 2    | 14             | 482-059                                      | 482-326                                    | 5    | 14             | 482-083                                      | 482-341                                    |
| 2    | 17             | 482-060                                      | 482-327                                    | 5    | 17             | 482-084                                      | 482-342                                    |
| 2    | 20             | 482-061                                      | 482-328                                    | 5    | 20             | 482-085                                      | 482-343                                    |
| 3    | 10             | 482-065                                      | 482-329                                    | 6    | 10             | 482-089                                      | 482-344                                    |
| 3    | 12             | 482-066                                      | 482-330                                    | 6    | 12             | 482-090                                      | 482-345                                    |
| 3    | 14             | 482-067                                      | 482-331                                    | 6    | 14             | 482-091                                      | 482-346                                    |
| 3    | 17             | 482-068                                      | 482-332                                    | 6    | 17             | 482-092                                      | 482-347                                    |
| 3    | 20             | 482-069                                      | 482-333                                    | 6    | 20             | 482-093                                      | 482-348                                    |

Insert STD, fixed bearing:

| Size | Height [mm] | Reference number  |                 | Size | Height [mm] | Reference number  |                 |
|------|-------------|-------------------|-----------------|------|-------------|-------------------|-----------------|
|      |             | ISO 5834-2 UHMWPE | Vitamin-E XL-PE |      |             | ISO 5834-2 UHMWPE | Vitamin-E XL-PE |
| 1    | 10          | 522-947           | 522-801         | 4    | 10          | 522-962           | 522-816         |
| 1    | 12          | 522-948           | 522-802         | 4    | 12          | 522-963           | 522-817         |
| 1    | 14          | 522-949           | 522-803         | 4    | 14          | 522-964           | 522-818         |
| 1    | 17          | 522-950           | 522-804         | 4    | 17          | 522-965           | 522-819         |
| 1    | 20          | 522-951           | 522-805         | 4    | 20          | 522-966           | 522-820         |
| 2    | 10          | 522-952           | 522-806         | 5    | 10          | 522-967           | 522-821         |
| 2    | 12          | 522-953           | 522-807         | 5    | 12          | 522-968           | 522-822         |
| 2    | 14          | 522-954           | 522-808         | 5    | 14          | 522-969           | 522-823         |
| 2    | 17          | 522-955           | 522-809         | 5    | 17          | 522-970           | 522-824         |
| 2    | 20          | 522-956           | 522-810         | 5    | 20          | 522-971           | 522-825         |
| 3    | 10          | 522-957           | 522-811         | 6    | 10          | 522-972           | 522-826         |
| 3    | 12          | 522-958           | 522-812         | 6    | 12          | 522-973           | 522-827         |
| 3    | 14          | 522-959           | 522-813         | 6    | 14          | 522-974           | 522-828         |
| 3    | 17          | 522-960           | 522-814         | 6    | 17          | 522-975           | 522-829         |
| 3    | 20          | 522-961           | 522-815         | 6    | 20          | 522-976           | 522-830         |

Insert UC, fixed bearing:

| Size | Height [mm] | Reference number | Size | Height [mm] | Reference number |
|------|-------------|------------------|------|-------------|------------------|
|      |             | Vitamin-E XL-PE  |      |             | Vitamin-E XL-PE  |
| 1    | 10          | 522-831          | 4    | 10          | 522-846          |
| 1    | 12          | 522-832          | 4    | 12          | 522-847          |
| 1    | 14          | 522-833          | 4    | 14          | 522-848          |
| 1    | 17          | 522-834          | 4    | 17          | 522-849          |
| 1    | 20          | 522-835          | 4    | 20          | 522-850          |
| 2    | 10          | 522-836          | 5    | 10          | 522-851          |
| 2    | 12          | 522-837          | 5    | 12          | 522-852          |
| 2    | 14          | 522-838          | 5    | 14          | 522-853          |
| 2    | 17          | 522-839          | 5    | 17          | 522-854          |
| 2    | 20          | 522-840          | 5    | 20          | 522-855          |
| 3    | 10          | 522-841          | 6    | 10          | 522-856          |
| 3    | 12          | 522-842          | 6    | 12          | 522-857          |
| 3    | 14          | 522-843          | 6    | 14          | 522-858          |
| 3    | 17          | 522-844          | 6    | 17          | 522-859          |
| 3    | 20          | 522-845          | 6    | 20          | 522-860          |

### 1.1.3 Patellar components

Patella STD, mobile bearing:

| Size | Reference number                                             |                                                                             |
|------|--------------------------------------------------------------|-----------------------------------------------------------------------------|
|      | cemented                                                     | cementless                                                                  |
|      | ISO 5832-12 CoCrMo<br>ISO 5834-2 UHMWPE<br>ISO 5832-1 1.4441 | ISO 5832-12 CoCrMo<br>TPS, BONIT®<br>ISO 5834-2 UHMWPE<br>ISO 5832-1 1.4441 |
| 1    | 482-109                                                      | 482-115                                                                     |
| 2    | 482-110                                                      | 482-116                                                                     |
| 3    | 482-111                                                      | 482-117                                                                     |
| 4    | 482-112                                                      | 482-118                                                                     |
| 5    | 482-113                                                      | 482-119                                                                     |
| 6    | 482-114                                                      | 482-120                                                                     |

Patella FIX, fixed bearing:

| Size | Reference number                       |
|------|----------------------------------------|
|      | cemented                               |
|      | ISO 5834-2 UHMWPE<br>ISO 5832-1 1.4441 |
| 1    | 482-163                                |
| 2    | 482-164                                |
| 3    | 482-165                                |
| 4    | 482-166                                |
| 5    | 482-167                                |
| 6    | 482-168                                |

## 12 Overview of Accessories

### 1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

| Designation                                       | Reference number |
|---------------------------------------------------|------------------|
| ZEN size 1 & 6 instrument set                     | 482-541          |
| ZEN Fixed Bearing size 1 & 6 instrument set       | 522-120          |
| ZEN Femur instrument set                          | 482-542          |
| ZEN General instrument set                        | 482-543          |
| ZEN Fixed Bearing General instrument set          | 522-121          |
| ZEN Trial Components instrument set               | 482-544          |
| ZEN Fixed Bearing Trial Components instrument set | 522-122          |
| ZEN Patella instrument set                        | 482-545          |
| ZEN Tibia instrument set                          | 482-546          |
| ZEN Tibia Modular instrument set I                | 521-103          |

| Designation                              | Reference number |
|------------------------------------------|------------------|
| ZEN Tibia Modular instrument set II      | 521-104          |
| ZEN Anterior Cut First instrument set    | 482-742          |
| ZEN Optional instrument set              | 482-775          |
| ZEN instrument set intermediate sizes I  | 482-781          |
| ZEN instrument set intermediate sizes II | 482-782          |

### 1.2.2 Accessories

| Designation                                                                                | Reference number |
|--------------------------------------------------------------------------------------------|------------------|
| ZEN Femur + Insert X-ray templates                                                         | 482-507          |
| ZEN Femur X-ray templates                                                                  | 482-559          |
| ZEN Tibia X-ray templates                                                                  | 482-508          |
| ZEN Tibia Fixed Bearing with Insert X-ray templates                                        | 522-501          |
| ZEN Tibia Modular X-ray templates                                                          | 482-738          |
| ZEN Tibia Modular Fixed Bearing with Insert and Stem Extension Straight X-ray templates    | 522-502          |
| ZEN Tibia Modular Offset Stem X-ray templates                                              | 482-739          |
| ZEN Tibia Modular Fixed Bearing with Insert and Stem Extension with offset X-ray templates | 522-503          |

## 1.3 Accompanying reference documents

| Designation                                         | Reference number |
|-----------------------------------------------------|------------------|
| ZEN Knee System Mobile Bearing surgical technique   | 50000434         |
| ZEN Tibia Modular Mobile Bearing surgical technique | 50000435         |
| ZEN Knee System Fixed Bearing surgical technique    | 50000656         |
| ZEN Tibia Modular Fixed Bearing surgical technique  | 50000657         |
| Implant Card                                        | 50000572         |
| Patient Information Knee Joint Replacement          | 50000836         |
| MRI safety information                              | 50000851         |
| List of symbols                                     | 50000859         |

## 2. Handling

### 2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>, must be observed.

**Caution:** Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

**Caution:** the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding side is marked on the implants with "L" for the left knee and "R" for the right knee. The orientation of the implants must always correspond to the side of the body of the knee to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right knee.

The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

## 2.2 Authorised Component Combinations

Together with its instruments, the ZEN Knee System is designed as a complete system and does not allow for the substitution of components from other systems of manufacturers. Only the above-mentioned system components may be used for implantation.

In the compatibility instructions, the combination options for the sizes of system components are presented according to the following scheme.

= Sizes can be combined

= Sizes **must not** be combined

# Compatibility instructions ZEN Mobile Bearing:

The compatibility of ZEN Mobile Bearing components should be inferred from the following overview.

|                       |     | Femur |   |     |   |     |   |     |   |     |   |
|-----------------------|-----|-------|---|-----|---|-----|---|-----|---|-----|---|
|                       |     | 1     | 2 | 2,5 | 3 | 3,5 | 4 | 4,5 | 5 | 5,5 | 6 |
| Tibia mobile bearing  | 1   | ✓     | ✓ |     | x |     | x |     | x |     | x |
|                       | 2   | ✓     | ✓ |     | ✓ |     | x |     | x |     | x |
|                       | 2,5 | ✓     | ✓ |     | ✓ |     | x |     | x |     | x |
|                       | 3   | ✓     | ✓ |     | ✓ |     | ✓ |     | x |     | x |
|                       | 3,5 | ✓     | ✓ |     | ✓ |     | ✓ |     | x |     | x |
|                       | 4   | ✓     | ✓ |     | ✓ |     | ✓ |     | ✓ |     | x |
|                       | 4,5 | ✓     | ✓ |     | ✓ |     | ✓ |     | ✓ |     | x |
|                       | 5   | ✓     | ✓ |     | ✓ |     | ✓ |     | ✓ |     | ✓ |
|                       | 5,5 | ✓     | ✓ |     | ✓ |     | ✓ |     | ✓ |     | ✓ |
|                       | 6   | ✓     | ✓ |     | ✓ |     | ✓ |     | ✓ |     | ✓ |
|                       |     | 1     | 2 |     | 3 |     | 4 |     | 5 |     | 6 |
| Insert mobile bearing |     |       |   |     |   |     |   |     |   |     |   |

# Compatibility instructions ZEN Fixed Bearing:

The compatibility of ZEN Fixed Bearing components should be inferred from the following overview.

|                      |     | Femur |   |     |   |     |   |     |   |     |   |
|----------------------|-----|-------|---|-----|---|-----|---|-----|---|-----|---|
|                      |     | 1     | 2 | 2,5 | 3 | 3,5 | 4 | 4,5 | 5 | 5,5 | 6 |
| Tibia fixed bearing  | 1   | ✓     | ✓ |     | x |     | x |     | x |     | x |
|                      | 2   | ✓     | ✓ |     | ✓ |     | x |     | x |     | x |
|                      | 2,5 | ✓     | ✓ |     | ✓ |     | x |     | x |     | x |
|                      | 3   | x     | ✓ |     | ✓ |     | ✓ |     | x |     | x |
|                      | 3,5 | x     | ✓ |     | ✓ |     | ✓ |     | x |     | x |
|                      | 4   | x     | x |     | ✓ |     | ✓ |     | ✓ |     | x |
|                      | 4,5 | x     | x |     | ✓ |     | ✓ |     | ✓ |     | x |
|                      | 5   | x     | x |     | x |     | ✓ |     | ✓ |     | ✓ |
|                      | 5,5 | x     | x |     | x |     | ✓ |     | ✓ |     | ✓ |
|                      | 6   | x     | x |     | x |     | x |     | ✓ |     | ✓ |
|                      |     | 1     | 2 |     | 3 |     | 4 |     | 5 |     | 6 |
| Insert fixed bearing |     |       |   |     |   |     |   |     |   |     |   |

# General compatibility instructions:

The femoral components, stem extensions and augments can be used for both the Fixed Bearing variant and for the Mobile Bearing variant.

The compatibility of patellar components should be inferred from the following overview.

|         |   | Femur |   |     |   |     |   |     |   |     |   |
|---------|---|-------|---|-----|---|-----|---|-----|---|-----|---|
|         |   | 1     | 2 | 2,5 | 3 | 3,5 | 4 | 4,5 | 5 | 5,5 | 6 |
| Patella | 1 | ✓     | x |     | x |     | x |     | x |     | x |
|         | 2 | x     | ✓ |     | x |     | x |     | x |     | x |
|         | 3 | x     | x |     | ✓ |     | x |     | x |     | x |
|         | 4 | x     | x |     | x |     | ✓ |     | x |     | x |
|         | 5 | x     | x |     | x |     | x |     | ✓ |     | x |
|         | 6 | x     | x |     | x |     | x |     | x |     | ✓ |

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

### 2.3 Instructions for Use

Before the instruments are used, the enclosed product information and instructions for use must be observed.

The surfaces with a porous coating (TPS, BONIT®) and the roughened implant surfaces must not come in contact with clothing or other fibre-shedding materials.

**Caution:** Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free gloves.

Follow the manufacturer's instructions on using the bone cement to ensure correct performance of the cementing step. The use of bone cement mixed under vacuum is recommended.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed. In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

**Caution:** when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

## **2.4 Surgical Technique**

### **2.4.1 ZEN Knee**

#### Approach

- Parapatellar incision with medial arthrotomy or another approach, e.g. for a revision procedure or moderate valgus malposition

#### Flexion gap - tibia processing step

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Mounting and positioning of the tibial goniometer
- Adjustment of rotation in the frontal plane
- Adjustment of the slope
- Adjustment of the cutting height
- Attachment of the tibial cutting block
- Performance of the tibial cut
- Verification of the tibial cut

#### Flexion gap - femur processing step

- Determination of the femur size
- Intramedullary opening
- Align and fix femoral cutting block in a.-p. position
- Anterior and posterior femoral resection

Caution: This step is a temporary step which has to be finalised in a later work step.

- Review of the flexion gap

#### Extension gap - distal femur processing step

- Distal positioning of the cutting block
- Review of the angle and distal cut
- Review of the extension gap

#### Finalisation of the femur

- Attachment of the femoral cutting block finisher
- Execution of the cuts

Caution: It is important to carry out all cuts precisely, as the previous A/P incision was only temporary.

- Verification of the cuts
- Preparation of the trochlear groove

#### Tibia processing step

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Selection of the plateau
- Tibia finalisation

### Patella processing step

For reduced bone thickness, the thinner full-poly patella should be employed, as it guarantees greater preservation of the bone. Prosthetic patellar replacement is contraindicated where the patella is too thin or in the case of a flat patella.

- Removal of the osteophytes
- Align and fix patellar resection guide perpendicular to the patellar tendon
- Measurement of the patellar thickness (at least 13 mm of bone must remain after deducting the thickness of the patellar component)

| Components | Height of patella STD | Height of patella FIX |
|------------|-----------------------|-----------------------|
| Size 1     | 11.2 mm               | 9.7 mm                |
| Size 2     | 11.2 mm               | 9.7 mm                |
| Size 3     | 11.2 mm               | 9.7 mm                |
| Size 4     | 11.2 mm               | 9.7 mm                |
| Size 5     | 11.7 mm               | 10.2 mm               |
| Size 6     | 11.7 mm               | 10.2 mm               |

- Adjustment of the height of the drilling gauge and insertion of the pins
- Removal of the patellar resection guide
- Positioning of the patellar cutting block
- Performance of the cut
- Alignment of the base plate of the trial patella
- Clamping of the patellar drilling guide into the patella forceps
- Positioning of the patella base plate and compression of the patellar forceps
- Drilling of the pin holes
- Removal of the patellar forceps and drilling guide
- Insertion of the trial patella insert

### Attachment of the trial components

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Components are inserted in the following order: Trial tibia, trial insert, trial femur and, if applicable, the trial patella

**Caution:** The implant components must not be used together in any random combination. Compatibility instructions are provided in section 2.2.

**Caution:** Should the height of the insert have to be changed or the trial components removed after successfully checking the overall mobility of the joint, the components are removed in the reverse order of the insertion. In other words: first remove the trial femur, followed by the trial insert.

### Implantation

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Components are inserted in the following order: Tibial components, insert, femoral component and, if applicable, the patellar component

**Caution:** The UC (ultra-congruent) inserts are to be used only after the resection of both cruciate ligaments. They do not provide any room for a posterior cruciate ligament that is still intact.

## 2.4.2 ZEN Tibia Modular

### Access

- Parapatellar incision with medial arthrotomy or another form of access, e.g. for a revision procedure or moderate valgus malposition

### Preparing the tibia - tibial incision

#### Intramedullary alignment:

- Medial-lateral centred opening of the tibia, aligned 1/3 anterior and 2/3 posterior
- Insertion of the IM rod
- If insufficient grip of the IM rod is achieved: Use of the medullary drills in ascending order, starting with the smallest diameter
- Mounting of the intramedullary (IM) alignment guide
- Mounting of the IM alignment guide on the in-situ IM rod/medullary drill
- Positioning of the IM alignment guide with the help of the axial alignment rod
- Adjusting the rotation in the frontal plane
- Performing varus-valgus adjustment
- Fixing the rotation and varus-valgus adjustment
- Adjustment of cutting height with tibial palp
- Checking the incision with a cortical palp
- Fixation of the tibial cutting block with pins and removal of the IM alignment guide and the IM rod/medullary drill
- Checking the alignment
- Performing the tibial incision
- Checking the tibial incision
- Optional subsequent femur preparation

### Preparing the tibia - final preparation with straight stem

- Selection of the plateau
- If a stem extension with offset is to be used, the further procedure according to the section "Preparing the tibia - final preparation with offsetstem" must be taken into account.
- Positioning and fixing of the trial tibia modular
- Drilling the tibial shaft
- Drilling the stem extension
- Preparing the tibial wings

### Preparing the tibia – final preparation with offset stem

- Selection of the plateau
- Positioning and fixing of the trial tibia modular
- Drilling the tibial shaft
- Drilling the stem extension
- Preparing the tibial wings

### Attachment of the trial components

- Sequence: Trial insert, trial femur (the trial tibia and possibly the trial patella had already been inserted in the previous step)

- Selection of the trial insert and trial femur according to the previously determined femur size.
- Caution: The implant components must not be used together in any random combination. Compatibility instructions are provided in section 2.2.
- Combine trial inserts for size and height and insert them into trial tibia
  - Place the trial femur with the femoral impactor on the femur, which is bent by at least 90°
  - Axis check
  - Insert trial patella if necessary
  - Check the overall mobility
- Caution: Should the height of the insert have to be changed or the trial components removed after successfully checking the overall mobility of the joint, remove the components in the reverse order of the insertion. In other words: first remove the trial femur and then the trial insert.

## Implantation

### When using a straight stem:

- The stem extension is placed onto the tibia and fixed in place with the screw provided and then tightened with the torque screwdriver
- With stem extension L110, the corresponding centraliser is also mounted

### When using an offset stem:

- The alignment of the stem extension with offset on the tibia modular is performed using the mounting device - then tightening of the screw provided
- Optional augments are attached using the torque screwdriver
- The components are inserted in the following order: Tibial component with stem extension and optional augments and centralisers, the insert, the femoral component and - where applicable - the patellar component

Caution: The UC (ultra-congruent) inserts are to be used only after the resection of both cruciate ligaments. They do not provide any room for a posterior cruciate ligament that is still intact.

## 2.4.3 General information

It is important to ensure that the surfaces are not damaged. Appropriate impactors are available for impacting the implants. Position the impactor at the centre and align it parallel to the mechanical axis. The patella is pushed into the bone using the patellar forceps and the pressing aid. In order to prevent excessive heating of the bone, all resections must be performed while cooling with physiological saline solution. The thickness of the saw blades must be adjusted to the instruments used! The smaller height of the fixed-bearing patellar component may provide greater bone preservation. For further information see the surgical technique of the knee system (see 1.3.). The cement is inserted using a modern cementing technique; in doing so, the instructions for use of the cement manufacturer must be observed.

## **3. Packaging and Sterility**

Depending on the sterilisation method used, implants are packaged in a double or triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made

of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/fu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

**Caution:** the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

#### **4. Preoperative Planning and Postoperative care**

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognised procedures must be used for postoperative care.

#### **5. Indications**

Total prosthesis is indicated in patients who have lost their mobility for any of the following reasons in a very advanced phase of the disease:

- Osteoarthritis of the knee (gonarthrosis)
- Osteonecrosis of the knee

Degree of severity:

- Knee pain
- Evidence of structural damage in the knee
- Failure of conservative treatment measures
- Limitation of quality of life related to the knee condition
- Subjective suffering related to the knee condition

This system may also be indicated for the correction of earlier surgical attempts which have been unsuccessful.

Total knee replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total knee replacement should be the course of action only where all other treatment approaches which could preserve

the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

## **6. Contraindications**

- Active bacterial infection in the knee (absolute)
- Circumstances that typically also apply as contraindications to other elective operations (e.g. infection, acute cardiovascular event, ASA classification > II) (absolute)
- Known allergy to implant materials (absolute)
- Insufficiency of one or both lateral ligaments in the knee joint (absolute)
- Severe muscle, nerve or vascular diseases that endanger the extremity concerned (e.g. primary myopathies, dystonia) (relative)
- Obesity: Severity BMI  $\geq 40$  (relative)
- Comorbidities that lead to shortened life expectancy (relative)
- Severe varus/valgus malposition of the knee (relative)

Additional ZEN Tibia Modular:

- Severely limited bone quality (relative)
- Moderate or severe bone losses of the proximal tibia (AORI classification > type 2) (relative)

## **7. Risk factors and conditions that may affect the success of the surgery**

**Caution:** Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of knee arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Use of high-dose cortisone or cytotoxic drugs
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to knee arthroplasty:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Growth in children and adolescents
- Expected extreme stresses, e.g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture
- Joint deformities that complicate fixation of the implant
- Weakening of the supporting structures due to a tumour

## **8. Undesirable Effects**

The residual risks listed below are among the most typical and commonly occurring consequences of knee arthroplasty:

- Soft tissue impingement / patella clunk
- Iliotibial ligament syndrome
- Foreign body reactions, osteolysis, loosening
- ARMD, trunnionosis, metallosis
- Toxic reactions
- Luxation, dislocation, dissociation
- Instability
- Implant failure
- Bone fracture
- Patellofemoral complications
- Limited ROM
- Ligamentous apparatus injury/rupture
- Instability of the ligamentous apparatus
- Bone cement reaction (bone cement implantation syndrome (BCIS))
- Loosening
- Tissue irritation
- Tissue damage
- Neurovascular damage
- Arthrofibrosis
- Heterotopic ossification
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack / stroke

Caution: Revision surgery may be required if specific undesirable effects occur.

## **9. Patient Information and Documentation**

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a knee replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. Based on the current state of technology, the expected service life is approx. 5 years.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant. Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

The user can download MRI safety information from the website <https://www.ohst.de/professionals/>.

The patient can download the patient information from the website <https://www.ohst.de/patient-information/>.

A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

## **10. Implant Card**

Following surgery, the patient must be given an implant card containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure; therefore, the implant card must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website. The patient's details (patient name or patient ID) must be documented in the implant card, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the card.

The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

## **11. Explanation of the Label Symbols**

The user can download the list of symbols from the website <https://www.ohst.de/professionals/>.

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