OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Germany

Tel.: +49 3385 5420-0

Fax: +49 3385 5420-99

www.ohst.de

CE 0482

Müller II Pfanne		Acetabul Müller II	
Deutsch (de)	3	Română (ro)	103
Müller II Cup		Müller II ацетабуларна чашка	
English (en)	15	Български (bg)	116
Cotyle Müller II		Чашка-западина Müller II	
Français (fr)	27	Українська (uk)	130
Cotile Müller II		Gūžduobė „Müller II“	
Italiano (it)	39	Lietuvių k. (lt)	146
Cótilo Müller II		Вертлужная впадина Müller II	
Español (es)	51	Русский (ru)	160
Müller II Copo		Müller II تجویف	
Português (pt)	63	العربية (ar)	179
Müller II Κυπέλιο		Müller II kauss	
Ελληνικά (el)	75	Latviešu (lv)	190
Kyčelní jamka Müller II		Капа Müller II	
Česky (cs)	88	Македонски (mk)	201

IMPLANTAT

Müller II Pfanne

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IMPLANT

Müller II Cup

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product. The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for hip arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

The most up-to-date version of the instructions for use must always be observed and is available for download on the website <https://www.ohst.de/ifu-instructions/>.

1. Product Description and Implant Materials

The Müller II Cup is a polyethylene acetabulum cup designed for cemented fixation in the acetabulum. It is manufactured using UHMWPE according to ISO 5834-2 and is available in the variants "Standard", "Flat", "Snap" and "Dysplasia". The Standard, Snap and Flat variants are available for spherical head diameters of 22 mm, 28 mm and 32 mm, while the Dysplasia variant is designed to retain femoral heads with a diameter of 28 mm and 32 mm. To ensure radiographic visibility, all implants are equipped with a ring made of implant steel (ISO 5832-1).

The Müller II Cup can generally be used in patients in whom hip arthroplasty is indicated and who fulfil the relevant indications and contraindications. The Müller II Cup is only approved for skeletally mature patients.

The product, the package contents and the materials used are specified on the labels. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

1.1 Overview of Implants

Designation	Material	Reference number
Müller II Cup Standard		
Mueller II Cup (Standard) Ø36, KD22	ISO 5834-2 UHMWPE	1300-22-36
Mueller II Cup (Standard) Ø38, KD22	ISO 5834-2 UHMWPE	1300-22-38
Mueller II Cup (Standard) Ø40, KD22	ISO 5834-2 UHMWPE	1300-22-40
Mueller II Cup (Standard) Ø42, KD22	ISO 5834-2 UHMWPE	1300-22-42
Mueller II Cup (Standard) Ø44, KD22	ISO 5834-2 UHMWPE	1300-22-44
Mueller II Cup (Standard) Ø42, KD28	ISO 5834-2 UHMWPE	1300-28-42
Mueller II Cup (Standard) Ø44, KD28	ISO 5834-2 UHMWPE	1300-28-44
Mueller II Cup (Standard) Ø46, KD28	ISO 5834-2 UHMWPE	1300-28-46
Mueller II Cup (Standard) Ø48, KD28	ISO 5834-2 UHMWPE	1300-28-48
Mueller II Cup (Standard) Ø50, KD28	ISO 5834-2 UHMWPE	1300-28-50

Designation	Material	Reference number
Mueller II Cup (Standard) Ø52, KD28	ISO 5834-2 UHMWPE	1300-28-52
Mueller II Cup (Standard) Ø54, KD28	ISO 5834-2 UHMWPE	1300-28-54
Mueller II Cup (Standard) Ø56, KD28	ISO 5834-2 UHMWPE	1300-28-56
Mueller II Cup (Standard) Ø58, KD28	ISO 5834-2 UHMWPE	1300-28-58
Mueller II Cup (Standard) Ø60, KD28	ISO 5834-2 UHMWPE	1300-28-60
Mueller II Cup (Standard) Ø62, KD28	ISO 5834-2 UHMWPE	1300-28-62
Mueller II Cup (Standard) Ø64, KD28	ISO 5834-2 UHMWPE	1300-28-64
Mueller II Cup (Standard) Ø66, KD28	ISO 5834-2 UHMWPE	1300-28-66
Mueller II Cup (Standard) Ø68, KD28	ISO 5834-2 UHMWPE	1300-28-68
Mueller II Cup (Standard) Ø70, KD28	ISO 5834-2 UHMWPE	1300-28-70
Mueller II Cup (Standard) Ø44, KD32	ISO 5834-2 UHMWPE	1300-32-44
Mueller II Cup (Standard) Ø46, KD32	ISO 5834-2 UHMWPE	1300-32-46
Mueller II Cup (Standard) Ø48, KD32	ISO 5834-2 UHMWPE	1300-32-48
Mueller II Cup (Standard) Ø50, KD32	ISO 5834-2 UHMWPE	1300-32-50
Mueller II Cup (Standard) Ø52, KD32	ISO 5834-2 UHMWPE	1300-32-52
Mueller II Cup (Standard) Ø54, KD32	ISO 5834-2 UHMWPE	1300-32-54
Mueller II Cup (Standard) Ø56, KD32	ISO 5834-2 UHMWPE	1300-32-56
Mueller II Cup (Standard) Ø58, KD32	ISO 5834-2 UHMWPE	1300-32-58
Mueller II Cup (Standard) Ø60, KD32	ISO 5834-2 UHMWPE	1300-32-60
Mueller II Cup (Standard) Ø62, KD32	ISO 5834-2 UHMWPE	1300-32-62
Mueller II Cup (Standard) Ø64, KD32	ISO 5834-2 UHMWPE	1300-32-64
Mueller II Cup (Standard) Ø66, KD32	ISO 5834-2 UHMWPE	1300-32-66
Mueller II Cup (Standard) Ø68, KD32	ISO 5834-2 UHMWPE	1300-32-68
Mueller II Cup (Standard) Ø70, KD32	ISO 5834-2 UHMWPE	1300-32-70
Müller II Cup Flat		
Mueller II Cup (Flat) Ø36, KD22	ISO 5834-2 UHMWPE	1301-22-36
Mueller II Cup (Flat) Ø38, KD22	ISO 5834-2 UHMWPE	1301-22-38
Mueller II Cup (Flat) Ø40, KD22	ISO 5834-2 UHMWPE	1301-22-40
Mueller II Cup (Flat) Ø42, KD22	ISO 5834-2 UHMWPE	1301-22-42
Mueller II Cup (Flat) Ø44, KD22	ISO 5834-2 UHMWPE	1301-22-44
Mueller II Cup (Flat) Ø42, KD28	ISO 5834-2 UHMWPE	1301-28-42
Mueller II Cup (Flat) Ø44, KD28	ISO 5834-2 UHMWPE	1301-28-44
Mueller II Cup (Flat) Ø46, KD28	ISO 5834-2 UHMWPE	1301-28-46
Mueller II Cup (Flat) Ø48, KD28	ISO 5834-2 UHMWPE	1301-28-48
Mueller II Cup (Flat) Ø50, KD28	ISO 5834-2 UHMWPE	1301-28-50
Mueller II Cup (Flat) Ø52, KD28	ISO 5834-2 UHMWPE	1301-28-52
Mueller II Cup (Flat) Ø54, KD28	ISO 5834-2 UHMWPE	1301-28-54
Mueller II Cup (Flat) Ø56, KD28	ISO 5834-2 UHMWPE	1301-28-56

Designation	Material	Reference number
Mueller II Cup (Flat) Ø58, KD28	ISO 5834-2 UHMWPE	1301-28-58
Mueller II Cup (Flat) Ø60, KD28	ISO 5834-2 UHMWPE	1301-28-60
Mueller II Cup (Flat) Ø62, KD28	ISO 5834-2 UHMWPE	1301-28-62
Mueller II Cup (Flat) Ø64, KD28	ISO 5834-2 UHMWPE	1301-28-64
Mueller II Cup (Flat) Ø66, KD28	ISO 5834-2 UHMWPE	1301-28-66
Mueller II Cup (Flat) Ø68, KD28	ISO 5834-2 UHMWPE	1301-28-68
Mueller II Cup (Flat) Ø70, KD28	ISO 5834-2 UHMWPE	1301-28-70
Mueller II Cup (Flat) Ø44, KD32	ISO 5834-2 UHMWPE	1301-32-44
Mueller II Cup (Flat) Ø46, KD32	ISO 5834-2 UHMWPE	1301-32-46
Mueller II Cup (Flat) Ø48, KD32	ISO 5834-2 UHMWPE	1301-32-48
Mueller II Cup (Flat) Ø50, KD32	ISO 5834-2 UHMWPE	1301-32-50
Mueller II Cup (Flat) Ø52, KD32	ISO 5834-2 UHMWPE	1301-32-52
Mueller II Cup (Flat) Ø54, KD32	ISO 5834-2 UHMWPE	1301-32-54
Mueller II Cup (Flat) Ø56, KD32	ISO 5834-2 UHMWPE	1301-32-56
Mueller II Cup (Flat) Ø58, KD32	ISO 5834-2 UHMWPE	1301-32-58
Mueller II Cup (Flat) Ø60, KD32	ISO 5834-2 UHMWPE	1301-32-60
Mueller II Cup (Flat) Ø62, KD32	ISO 5834-2 UHMWPE	1301-32-62
Mueller II Cup (Flat) Ø64, KD32	ISO 5834-2 UHMWPE	1301-32-64
Mueller II Cup (Flat) Ø66, KD32	ISO 5834-2 UHMWPE	1301-32-66
Mueller II Cup (Flat) Ø68, KD32	ISO 5834-2 UHMWPE	1301-32-68
Mueller II Cup (Flat) Ø70, KD32	ISO 5834-2 UHMWPE	1301-32-70
Müller II Cup Snap		
Mueller II Cup (Snap) Ø36, KD22	ISO 5834-2 UHMWPE	1310-22-36
Mueller II Cup (Snap) Ø38, KD22	ISO 5834-2 UHMWPE	1310-22-38
Mueller II Cup (Snap) Ø40, KD22	ISO 5834-2 UHMWPE	1310-22-40
Mueller II Cup (Snap) Ø42, KD22	ISO 5834-2 UHMWPE	1310-22-42
Mueller II Cup (Snap) Ø44, KD22	ISO 5834-2 UHMWPE	1310-22-44
Mueller II Cup (Snap) Ø42, KD28	ISO 5834-2 UHMWPE	1310-28-42
Mueller II Cup (Snap) Ø44, KD28	ISO 5834-2 UHMWPE	1310-28-44
Mueller II Cup (Snap) Ø46, KD28	ISO 5834-2 UHMWPE	1310-28-46
Mueller II Cup (Snap) Ø48, KD28	ISO 5834-2 UHMWPE	1310-28-48
Mueller II Cup (Snap) Ø50, KD28	ISO 5834-2 UHMWPE	1310-28-50
Mueller II Cup (Snap) Ø52, KD28	ISO 5834-2 UHMWPE	1310-28-52
Mueller II Cup (Snap) Ø54, KD28	ISO 5834-2 UHMWPE	1310-28-54
Mueller II Cup (Snap) Ø56, KD28	ISO 5834-2 UHMWPE	1310-28-56
Mueller II Cup (Snap) Ø58, KD28	ISO 5834-2 UHMWPE	1310-28-58
Mueller II Cup (Snap) Ø60, KD28	ISO 5834-2 UHMWPE	1310-28-60
Mueller II Cup (Snap) Ø62, KD28	ISO 5834-2 UHMWPE	1310-28-62

Designation	Material	Reference number
Mueller II Cup (Snap) Ø64, KD28	ISO 5834-2 UHMWPE	1310-28-64
Mueller II Cup (Snap) Ø66, KD28	ISO 5834-2 UHMWPE	1310-28-66
Mueller II Cup (Snap) Ø68, KD28	ISO 5834-2 UHMWPE	1310-28-68
Mueller II Cup (Snap) Ø70, KD28	ISO 5834-2 UHMWPE	1310-28-70
Mueller II Cup (Snap) Ø44, KD32	ISO 5834-2 UHMWPE	1310-32-44
Mueller II Cup (Snap) Ø46, KD32	ISO 5834-2 UHMWPE	1310-32-46
Mueller II Cup (Snap) Ø48, KD32	ISO 5834-2 UHMWPE	1310-32-48
Mueller II Cup (Snap) Ø50, KD32	ISO 5834-2 UHMWPE	1310-32-50
Mueller II Cup (Snap) Ø52, KD32	ISO 5834-2 UHMWPE	1310-32-52
Mueller II Cup (Snap) Ø54, KD32	ISO 5834-2 UHMWPE	1310-32-54
Mueller II Cup (Snap) Ø56, KD32	ISO 5834-2 UHMWPE	1310-32-56
Mueller II Cup (Snap) Ø58, KD32	ISO 5834-2 UHMWPE	1310-32-58
Mueller II Cup (Snap) Ø60, KD32	ISO 5834-2 UHMWPE	1310-32-60
Mueller II Cup (Snap) Ø62, KD32	ISO 5834-2 UHMWPE	1310-32-62
Mueller II Cup (Snap) Ø64, KD32	ISO 5834-2 UHMWPE	1310-32-64
Mueller II Cup (Snap) Ø66, KD32	ISO 5834-2 UHMWPE	1310-32-66
Mueller II Cup (Snap) Ø68, KD32	ISO 5834-2 UHMWPE	1310-32-68
Mueller II Cup (Snap) Ø70, KD32	ISO 5834-2 UHMWPE	1310-32-70
Müller II Cup Dysplasia		
Mueller II Cup (Dysplasia) Ø42, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-42
Mueller II Cup (Dysplasia) Ø44, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-44
Mueller II Cup (Dysplasia) Ø46, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-46
Mueller II Cup (Dysplasia) Ø48, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-48
Mueller II Cup (Dysplasia) Ø50, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-50
Mueller II Cup (Dysplasia) Ø52, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-52
Mueller II Cup (Dysplasia) Ø54, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-54
Mueller II Cup (Dysplasia) Ø56, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-56
Mueller II Cup (Dysplasia) Ø58, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-58
Mueller II Cup (Dysplasia) Ø60, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-60
Mueller II Cup (Dysplasia) Ø62, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-62
Mueller II Cup (Dysplasia) Ø64, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-64
Mueller II Cup (Dysplasia) Ø66, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-66
Mueller II Cup (Dysplasia) Ø68, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-68
Mueller II Cup (Dysplasia) Ø70, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-70
Mueller II Cup (Dysplasia) Ø44, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-44
Mueller II Cup (Dysplasia) Ø46, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-46
Mueller II Cup (Dysplasia) Ø48, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-48
Mueller II Cup (Dysplasia) Ø50, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-50

Designation	Material	Reference number
Mueller II Cup (Dysplasia) Ø52, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-52
Mueller II Cup (Dysplasia) Ø54, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-54
Mueller II Cup (Dysplasia) Ø56, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-56
Mueller II Cup (Dysplasia) Ø58, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-58
Mueller II Cup (Dysplasia) Ø60, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-60
Mueller II Cup (Dysplasia) Ø62, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-62
Mueller II Cup (Dysplasia) Ø64, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-64
Mueller II Cup (Dysplasia) Ø66, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-66
Mueller II Cup (Dysplasia) Ø68, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-68
Mueller II Cup (Dysplasia) Ø70, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-70

1.2 Overview of accessories

1.2.1 Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
Acetabulum Reamer Instrument Set Müller II	367-1578
Aceatabular Size Testers Instrument Set	367-1573
Mueller II Cup Insertion Instrument Set	367-400
Optional:	
Optional Instruments Müller II	367-1580

1.2.2 Accessories

Designation	Reference number
X-ray Templates Mueller II Cup Dysplasia Ø28/Ø32	367-2001
X-ray Templates Mueller II Cup Standard Ø22	367-2002
X-ray Templates Mueller II Cup Standard Ø28/Ø32	367-2003
X-ray Templates Mueller II Cup Flat Profile Ø28/Ø32	367-2004
X-ray Templates Mueller II Cup Snap Fit Ø28/Ø32	367-2005
X-ray Templates Mueller II Cup Snap Fit Ø22	367-2006
X-ray Templates Mueller II Cup Flat Profile Ø22	367-2010

1.3 Accompanying reference documents

Designation	Reference number
Surgical Technique for Müller II Cup	50000402
Implant Card	50000572
Patient information Hip joint replacement	50000841
MRI safety information	50000851
List of symbols	50000859

2. Handling

2.1 General information

This implant is part of a system and must only be used with the corresponding original system components. Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before using the instruments, the most up-to-date version of the respective instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>, must be observed.

Caution: Implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

Caution: the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding orientation is marked on the implants with "L" for the left side of the body and "R" for the right side of the body. The orientation of the implants must always correspond to the side of the body of the joint to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right joint.

The packaging components and implants are to be passed for waste recycling in accordance with their materials and the statutory provisions.

In agreement with the manufacturer, these implants may also be returned to the manufacturer for proper disposal at no charge. The return consignment to the manufacturer must be marked with the words "Return for disposal"; implants must be cleaned and sterilised and proof of decontamination or a hygienic clearance certificate must be provided.

Any serious incidents that occur in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient resides.

2.2 Authorised Component Combinations

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In

this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

2.3 Instructions for Use

The application of the implant is done with cement. For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

Caution: The ball diameter of the femoral head prosthesis must agree perfectly with the spherical nominal diameter of the part of the acetabular cup implant which articulates as the sliding partner.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed.

Caution: when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

When using Snap or Dysplasia cups or cup inserts, the movement range in flexion and extension is reduced by approximately 40° (Snap) / 8° (Dysplasia) and in abduction and adduction by approximately 22° (Snap) / 4° (Dysplasia) as compared to the standard application.

Caution: Component malposition or the use of a constrained hip or a total hip with a modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement, premature dislocation or revision. In these

cases, the surgeon should inform the patient that high range of motion activities should be avoided.

2.4 Surgical Technique

Any approach to the hip joint deemed suitable by the surgeon can be used. The surgeon should have a good view of the anatomical structures so that correct working with the instruments is not impeded. After opening the joint capsule and dislocating the femoral head from the acetabulum, resect the femoral head as determined in the preoperative planning of the hip stem.

Starting with the smallest, ream the acetabulum using spherical acetabulum reamers in ascending order; in doing so, make sure to always let the reamers move freely and never apply force. Always allow the reamer to move freely and do not distort the reamer sizes by applying pressure laterally! The base of the cup should be completely free from cartilage and the subchondral bone should bleed evenly. The bony mass collected from the last reaming process can be used to fill in any potential gaps between the implant and the acetabulum.

Caution: The diameter (nominal size) of the last acetabulum reamer used has to be at least 4 mm larger than the smallest Müller II Cup available to ensure a cement coating of 2 mm all around.

Next, insert the respective equatorial size tester based on the size of the reamer and use it to check the alignment (45° inclination and 10° anteversion) of the hip cup. The equatorial size tester should remain in place through suction and be surrounded by a sufficient quantity of bone. After these surgical steps are complete, select the appropriate implant size.

Follow the cement manufacturer's instructions for use to prepare the bone cement using a state-of-the-art cementing technique and insert it into the prepared implant bed.

Caution: Make sure that a cement coating of around 2 mm is achieved. Use utmost care when removing excess cement to prevent damaging the cup or to prevent cement particles from penetrating the inner articulation surface.

The size of the implant to be used has to be 4 mm smaller than the size of the equatorial size tester selected to ensure an adequate cement coating. Implant the Müller II Cup with the corresponding insertion head that is screwed on to the curved handle. To this purpose, press the implant into the bone cement and align it in accordance with the pre-operative planning process. Then apply light pressure and wait until the bone cement has set. This completes the insertion of the cup. In order to protect the cup against damage during the remainder of the operation it is advisable to cover it up with a sterile compress. After all of the other components have been implanted, close the wound using the routine procedure.

When using the snap variant of the Müller II cup, greater resistance has to be overcome to reduce the femoral head into the cup than with the other Müller II cup variants. When using the snap variant of the Müller II cup, the femoral head can only be reduced after flexion/anteversion of the thigh until the chamfer in the neck area of the femoral head protrudes into the internal contour of the cup, allowing air to escape from the joint body. Otherwise an air pocket would form, thus preventing the femoral head from snapping into place. After reducing the femoral head, verify the correct fit of the femoral head inside the cup. The connection between the two implants is correct when the femoral head noticeably snaps into the cup.

3. Packaging and Sterility

Depending on the sterilisation method used, implants are packaged in a double or triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. The instruments are supplied unsterile in protective packaging and must, prior to use, be cleaned and sterilised as per the most up-to-date version of the instructions for use (50000354), which is available for download on the website <https://www.ohst.de/ifu-instructions/>. The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative Planning and Postoperative Care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

The X-ray templates must be used to plan the surgery. They are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a scale of 1:1 are available in digital format.

To check the correct fit (where applicable), trial prostheses and additional implant components should be kept to hand in the event that other sizes are required, or where it is not possible to use the planned implant. Recognised procedures must be used for postoperative care.

5. Indications

- Advanced degeneration of the hip joint due to degenerative, post-traumatic or rheumatoid arthritis
- Avascular necrosis of the femoral head
- Displaced intracapsular hip fracture
- Sequelae of earlier surgical procedures (osteosynthesis, articular reconstruction, arthrodesis, etc.)

Total hip replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total hip replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

6. Contraindications

- Acute local (in relation to the joint undergoing surgery) or active systemic infection
- Insufficient bone substance or quality that does not allow for stable fixation of the implant
- Hypersensitivity to the materials used
- Severe osteoporosis
- Severe malformations, congenital hip dislocations
- Local bone tumours that preclude stable fixation of the implant
- Growth in children and adolescents

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Use of high-dose cortisone or cytotoxic drugs
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour (e.g. bone cyst, non-ossifying fibroma)
- Loss of the ligamenture

8. Undesirable Effects

The undesirable effects and risks listed below can occur in the course of a total hip arthroplasty:

- Foreign body reactions, osteolysis, loosening
- Toxic reactions
- Sensitisation
- Limited ROM
- Instability
- Luxation, dislocation, dissociation

- Implant failure
- Loosening
- Bone cement reaction/bone cement implantation syndrome (BCIS) (e.g. cardiac arrhythmia, increased pulmonary vascular resistance)
- Drifting/tilting
- Heterotopic ossification
- Difference in leg length
- Tissue damage
- Deep vein thrombosis
- Blood loss
- Infections
- Pulmonary embolism
- Cardiac arrest
- Heart attack/stroke

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient Information and Documentation

The identification details of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks.

Patients who receive a hip replacement must be informed that the service life of the implant depends on a variety of factors; therefore, it is not possible to specify an expected service life. The service life depends on the patient's weight and degree of activity, the existing bone quality, accompanying diseases, the tribological pairing chosen, the quality of the implantation as well as unexpected complications, for example, due to falls or accidents. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances. Based on the current state of technology, the expected service life is approx. 15 years.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant.

Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

The user can download MRI safety information from the website <https://www.ohst.de/professionals/>. The patient can download the patient information from the website <https://www.ohst.de/patient-information/>. A short report about the safety and clinical performance is available from the Eudamed database (<https://ec.europa.eu/tools/eudamed>). Until the database is launched, the short report can be provided on request.

10. Implant Card

Following surgery, the patient must be given an implant card containing all the necessary information about the implant. Several components of a system are used during the initial surgical procedure;

therefore, the implant card must be obtained directly from OHST Medizintechnik AG. Adhesive labels to document the implant used are enclosed with the products. These labels contain the product designation, the article numbers (REF), the serial numbers (SN), the UDI code and the manufacturer's details including the website.

The patient's details (patient name or patient ID) must be documented in the implant card, along with the date of implantation as well as the name and the address of the implanting health care facility; for each implanted component, a label must be affixed to the respective section of the card.

The user must inform the patient that any additional or updated information intended to ensure the safe use of the product by the patient is available on the indicated website.

11. Explanation of the Label Symbols

The user can download the list of symbols from the website <https://www.ohst.de/professionals/>.

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