

PE and resorbable cement restrictors

⇒ a cement restrictor, also known as cement plug is used in orthopaedic surgeries, prior to the introduction of bone cement into the medullary cavity and implantation of a prosthesis

Cement restrictor UHMWPE

| Article No. | Sizes |
|-------------|-------|
| 506-101 | 1 |
| 506-102 | 2 |



- ⇒ Material: UHMWPE (ISO 5834-2)
- ⇒ available in two sizes
- ⇒ one single setting instrument with measuring scale for controlled insertion in optimal depth



Only one setting instrument is necessary for implantation of PE cement restrictor.

Setting Instrument

| Article No. | Designation |
|-------------|--|
| 506-100 | Insertion Instrument for Cement Restrictor |



Neptun® Plugs

| Article No. | Sizes |
|-------------|-------|
| NP-770008 | 8 |
| NP-770010 | 10 |
| NP-770012 | 12 |
| NP-770014 | 14 |
| NP-770016 | 16 |
| NP-770018 | 18 |



- ⇒ made from elastic gelatine-based biomaterials
- ⇒ fully resorbed within days
- ⇒ simplifies cement removal during revisions
- ⇒ available in six sizes
- ⇒ precise setting instrument including trials for size testing available



Our instruments for Neptun® Plugs, will be supplied in an easily handle PP tray, which can be added to

Setting Instruments

| Article No. | Designation |
|-------------|---|
| 506-100 | Insertion Instrument for Cement Restrictor |
| 018-001 | Slider for Insertion Instrument |
| 206-021 | Holder for Trial Cement Restrictors |
| 367-080 | Trial medullary plug with attachment sz. 8 |
| 367-081 | Trial medullary plug with attachment sz. 10 |
| 367-082 | Trial medullary plug with attachment sz. 12 |
| 367-083 | Trial medullary plug with attachment sz. 14 |
| 367-084 | Trial medullary plug with attachment sz. 16 |
| 367-085 | Trial medullary plug with attachment sz. 18 |