

FlexiBit system

The FlexiBit is a modular-flexible drill with AO connection. The quick coupling system allows for easy tool change.

The FlexiBit fulfils all requirements for cleaning and sterilisation within the scope of reprocessing as a medical device.

In order to avoid incorrect usage the Flexi Bit shaft should always be used in combination with 367-165 Drill Guide.

FlexiBit-Shaft

Art. no. Description

367-1457 FlexiBit-shaft, L 110 mm



FlexiBit-Drill

Art. no.	Ø drill	L	LS
367-1447	3,2	25	15
367-1448	3,2	30	20
367-1449	3,2	40	30
367-1450	3,2	50	40
367-1451	3,2	60	50
367-1452	3,2	70	60
367-1453	3,2	80	70



Drill Guide for screws Ø 3,2 mm

Art.-No. 367-165

- total length: 185 mm
- guide length: 25 mm



Hexagon screw driver SW 3,5

Art.-No. 250161-5

- total length: 250 mm
- length of handle: 104 mm



Hexagon screw driver SW 3,5
conical 5°

Art.-No. 206-007

- length: 390 mm
- length of handle: 104 mm



Screw measuring instrument, curved

Art.-No. 367-115

- length: 160 mm
- outer diameter: 8 mm
- measuring range: 60 mm



Screw holding forcep

Art.-No. 367-1021

- total length: 180 mm
- usage: for Flat-head screw
Ø 6,5 mm



Cardan screw driver,
wrench size 3,5 mm

Art.-No. 00-092-10

- total length: 273 mm
- length of handle: 105 mm
- wrench size: 3,5 mm

