

The insertion instrument for BIOLOX® ceramic inserts

1/1

BIOLOX® inserts from CeramTec has been represented the state of the art since many years. We have developed a special insertion instrument that provides a safe placement of the ceramic inserts of the BIOLOX® series XLW 18 and XLW 18 zero.

The complete insertion instrument set consists of the three components mounting device, connector and liner handle.

CeramTec
 THE CERAMIC EXPERTS
**Liner handle**

- length: 384 mm, height: 55 mm
- with silicone handle
(standard colour grey, other colours available for orders exceeding 25 pieces)

**Connector**

- for head sizes: Ø 28, 32, 36, 40

**Mounting device**

- for insert sizes: Ø 35, 36, 37, 38, 39, 40, 41, 42, 44, 46, 48, 50, 52



Liner handle		Art. no. 510-000			
Connector		Ø28 (Art. no. 510-100)	Ø32 (Art. no. 510-200)	Ø36 (Art. no. 510-300)	Ø40 (Art. no. 510-400)
Mounting device XLW 18					
510-001	Ø35 XLW 18	XLW 18 28/35G			
510-002	Ø36 XLW 18	XLW 18 28/36G			
510-004	Ø37 XLW 18	XLW 18 28/37G			
510-006	Ø39 XLW 18	XLW 18 28/39G	XLW 18 32/39G		
510-008	Ø41 XLW 18	XLW 18 28/41G	XLW 18 32/41G		
510-010	Ø44 XLW 18	XLW 18 28/44G	XLW 18 32/44G	XLW 18 36/44G	
510-012	Ø48 XLW 18	XLW 18 28/48G	XLW 18 32/48G	XLW 18 36/48G	XLW 18 40/48G
510-014	Ø52 XLW 18	XLW 18 28/52G	XLW 18 32/52G	XLW 18 36/52G	XLW 18 40/52G
Mounting device XLW 18 zero					
510-003	Ø36 XLW 18 zero	XLW 18 28/36 zero			
510-005	Ø38 XLW 18 zero	XLW 18 28/38 zero			
510-007	Ø40 XLW 18 zero	XLW 18 28/40 zero	XLW 18 32/40 zero		
510-009	Ø42 XLW 18 zero		XLW 18 32/42 zero		
510-010	Ø44 XLW 18 zero		XLW 18 32/44 zero	XLW 18 36/44 zero	
510-011	Ø46 XLW 18 zero		XLW 18 32/46 zero	XLW 18 36/46 zero	
510-012	Ø48 XLW 18 zero		XLW 18 32/48 zero	XLW 18 36/48 zero	XLW 18 40/48 zero
510-013	Ø50 XLW 18 zero		XLW 18 32/50 zero	XLW 18 36/50 zero	XLW 18 40/50 zero
510-014	Ø52 XLW 18 zero		XLW 18 32/52 zero	XLW 18 36/52 zero	XLW 18 40/52 zero