

PE and resorbable cement restrictors

- ⇒ a cement restrictor, also known as cement plug is used in orthopaedic surgeries, prior to the introduction of bone cement into the medullary cavity and implantation of a prosthesis

Cement restrictor UHMWPE	
Article No.	Sizes
506-101	1
506-102	2



- ⇒ Material: UHMWPE (ISO 5834-2)
- ⇒ available in two sizes
- ⇒ one single setting instrument with measuring scale for controlled insertion in optimal depth



Only one setting instrument is required to implant the PE cement restrictor.


Setting Instrument	
Article No.	Designation
506-100	Insertion Instrument for Cement Restrictor



Neptun-Plug®	
Article No.	Sizes
NP-770008	8
NP-770010	10
NP-770012	12
NP-770014	14
NP-770016	16
NP-770018	18



- ⇒ made from elastic gelatine-based biomaterials
- ⇒ fully resorbed within days
- ⇒ simplifies cement removal during revisions
- ⇒ available in six sizes
- ⇒ precise setting instrument including trials for size testing available



Our Neptun-Plug® instruments are supplied in an easy-to-handle PP tray.

Setting Instruments	
Article No.	Designation
506-100	Insertion Instrument for Cement Restrictor
018-001	Slider for Insertion Instrument
206-021	Holder for Trial Cement Restrictors
367-080	Trial medullary plug with attachment sz. 8
367-081	Trial medullary plug with attachment sz. 10
367-082	Trial medullary plug with attachment sz. 12
367-083	Trial medullary plug with attachment sz. 14
367-084	Trial medullary plug with attachment sz. 16
367-085	Trial medullary plug with attachment sz. 18