

Pins & Adapters

For precise cutting block fixation



Pins & Adapters

Areas of application

Stable fixation and high precision

Our pins enable firm and precise fixation of cutting blocks to the bone

Minimal movement

By anchoring directly to the bone, pins offer high resistance to displacement or loosening, even during prolonged surgery or under mechanical stress

Reusable

The pins are designed for multiple sterilisations and reuse. They are easy to insert and remove for an efficient surgical procedure

Biocompatible steel materials

Easy handling

Quick insertion and removal by simply pulling out

Areas of application for pins

- Orthopaedics
- Traumatology



Pins & adapters

Quick pin fixation without prior drilling

Impact pins can be inserted directly into the bone without prior drilling.

Standard in Ø3.2 and Ø2.4

Pins Ø3.2



514-016 Pin Impactor 3.2



482-556 Pin Adaptor 3.2 Zimmer-Connection



514-006 Pin Ø3.2x50



482-558 Pin Ø3.2x70



514-018 Head Pin Ø3.2x45

Pins Ø2.4



482-557 Pin Adapter 2.4 Zimmer connection



482-553 Pin Ø2.4x80

Available in different versions.
Enquire and test now!

Pins & adapters

Drill pin and adapter

Drill Pins are the preferred choice for attaching cutting blocks to sclerotic bone.

The integrated magnetic function ensures a secure hold and allows quick, tool-free pin changes via the convenient quick-change system.

Advantages at a glance:

Self-centring design supports precise positioning

Optimised geometry for smooth drilling

High concentricity for stable guidance and precise alignment



482-833-Z Drill Pin Adaptor Zimmer-Connection



482-834-5 Drill Pin Ø3,2x90



Pins & Adapters

Pin Extraction

Pin-Extractor for reliable pin removal during surgical procedures

- **Practical one-hand grip mechanism**
- **Reusable**

Compatible with common pin systems and adapters

For pin diameters 3.2 and 2.4



514-005 Pin Extractor

Pins & Adapters

Article list

Item no.	Designation	Pieces per Packaging Unit (PU)
514-016	Pin Impactor 3.2 (for 514-006 and 482-558)	1
514-006	Pin Ø3.2x50	1
482-558	Pin Ø3,2x70	1
482-556	Pin Adaptor 3.2 Zimmer-Connection (for 482-558)	1
514-018	Head Pin Ø3.2x45	1
482-553	Pin Ø2.4x80	1
482-557	Pin Adaptor 2.4 Zimmer connection	1

Drill pin and adapter

482-834-5	Drill Pin Ø3,2x90	5
482-833-Z	Drill Pin Adaptor Zimmer-Connection	1

Pin Extractor

514-005	Pin Extractor	1
---------	---------------	---

Contact & Order



+49 3385 5420 0



sales@ohst.de



OHST Medizintechnik AG
Grünauer Fenn 3
14712 Rathenow
Germany

