

FlexiBit System

The FlexiBit System is a **modular and flexible drilling system** featuring an AO coupling. Its quick-connect mechanism enables fast and easy instrument changes, enhancing efficiency during surgical procedures.

Designed for **repeated clinical use**, the FlexiBit System meets all applicable requirements for cleaning, disinfection, and sterilisation as part of the reprocessing of reusable medical devices.

To ensure proper use and optimal performance, the FlexiBit System should always be used together with the Drill Guide (Ref. No. 367-165 resp. 367-164).



FlexiBit-Shaft


Art. No. Description

367-1457 FlexiBit-Shaft, L 110 mm

367-1458 FlexiBit-Shaft, L 180 mm

FlexiBit-Drill

Art.-No.	Ø Drill	L	LS
367-1447	3,2	25	15
367-1448	3,2	30	20
367-1449	3,2	40	30
367-1450	3,2	50	40
367-1451	3,2	60	50
367-1452	3,2	70	60
367-1453	3,2	80	70
367-1454	4,5	25	15
367-1455	4,5	40	30
367-1456	4,5	60	50



Drill Guide

Art.-No. 367-165 for Drill Ø 3,2

Art.-No. 367-165 for Drill Ø 4,5

- total length: 185 mm
- guide length: 25 mm



Hexagon Screw Driver SW 3,5

Art.-No. 250161-5

- total length: 250 mm
- length of handle: 104 mm



Hexagon Screw Driver A/F 3,5 conical 5°

Art.-No. 206-007

- length: 390 mm
- handle length: 104 mm



Screw Measuring Instrument Curved

Art.-No. 367-115

- length: 160 mm
- outer diameter: 8 mm
- measuring range: 60 mm



Screw Holding Forceps

Art.-No. 367-1021

- total length: 182 mm
- for Flat-head screw Ø 6,5 mm



Cardan Screw Driver, Wrench 3,5mm

Art.-No. 00-092-10

- total length: 273 mm
- handle length: 105 mm
- wrench size: 3,5 mm

